



OTS DOUBLE ADJUSTABLE ANKLE JOINT



- All stainless steel construction for stronger, longer-lasting control
- Lightweight, total weight only 1.5 oz per hinge
- Taller, reversible stainless steel upright helps to corrugate plastic strut, reducing the need for thicker plastics
- Allows AFOs to be fabricated with plastic as thin as 3/32 in.
- Adjustable stops allow movement from 0° to 30° in both direction
- Ultra-small joint head to reduce overall bulk
- Super low-profile at just over 1/4 in. thick (lower profile than most polyurethane hinges)
- Contourable, distal stirrup attachment increases strength of connection to molded footplate

HINGE ADJUSTMENT INSTRUCTIONS

1. Determine proper position or ROM for hinge
2. Loosen two faceplate screws (1, 3)
3. Adjust set screws (A, B) as to provide desired position or ROM
4. Tighten two faceplate screws (1, 3) to lock adjustable set screws (A, B)
5. Take care to avoid over tightening of center axis screw (2) during assembly



Apply (medium strength) thread lock compound (not included) to all screws (1-5) during final assembly.

ORDERING INFORMATION

- DAAJ1** OTS Double Adjustable Ankle Joint, Small
DAAJ2 OTS Double Adjustable Ankle Joint, Large

SPECIFICATIONS

Weight rating for DAAJ1: 50 lbs.
Weight rating for DAAJ2: 80 lbs.
Thickness: 0.25 in.



OTS Corp.

220 Merrimon Ave
Weaverville, NC 28787
800.221.4769
828.658.8330
www.ots-corp.com

OTS DOUBLE ADJUSTABLE ANKLE JOINT



**FABRICATION
MANUAL**



Prior to filling impression, determine what method will be used to contour hinges:

- Method 1** 4-40 Threaded rod through positive model
Method 2 Contouring / Alignment fixture

METHOD 1: USING A THREADED ROD

Mold Preparation

- Locate medial and lateral joint axis and position a 4-40 threaded rod, accordingly, prior to filling cast.
- Modify positive model as desired.
- Minimum of 3/16 in. build up on both medial and lateral malleoli is recommended.




Figure 1

Hinge Preparation

- Thread a 4-40 nut onto both medial and lateral sides and tighten against positive model.
- Nuts should act as spacers to keep hinges approximately 1/8 in. off positive model.
- Remove center axis screw from hinges.
- Most proximal assembly screw on faceplate.
- Slide hinge (at center axis point) over threaded rod.




Figure 2

Continue with instructions on next page

METHOD 2: USING CONTOURING/ALIGNMENT FIXTURE

Mold Preparation

- Modify positive model as desired.
- Minimum of 3/16 in. build up on both medial and lateral malleoli is recommended.

Hinge Preparation

- Attach medial and lateral hinges to fixture and align both with the desired joint axis.

Continue with instructions on next page

INSTRUCTIONS CONTINUED (USE FOR BOTH METHODS)

1. Contour proximal uprights and "AP Stabilizer Tab" to positive model. See Figure 3.
2. Contour distal "Y" stirrups to positive model. See Figure 4.
3. Remove hinges and attach 1/8 in. cork spacers (provided) to ensure proper spacing between hinge and positive model. See Figure 5.
4. Affix hinges to positive model using nails or small wood screws. See Figure 6.
5. Fill any gaps between positive model and hinges with modeling clay or similar material.
6. Vacuum form using standard procedures.
7. After cooling separate proximal and distal sections and remove AFO as normal.

It is recommended that a hot knife be used over joint head, upright or stirrup during removal.

DO NOT USE CAST CUTTER OR SAW OVER STAINLESS STEEL JOINT HEAD, AS IT MAY DAMAGE THE HINGE.

8. Finish all plastic edges as normal.




Figure 3




Figure 4




Figure 5

9. If desired for additional stirrup or upright retention, mark any additional holes for assembly on "Y" stirrup attachment and proximal upright. See Figure 7.

Holes should be marked and drilled (using a # 21 drill bit) and located as close to joint head as plastic trim lines will allow.

DO NOT SKIP THIS STEP.

10. Mark and drill remaining assembly holes in plastic, for attachment of hinges (# 11 drill). The included 10-32 countersunk head screws are for this application. The screws are then inserted through the plastic into the threaded holes using thread locking compound on all screws to avoid loosening.
11. Assemble AFO by installing hinges.
12. It is recommended that proximal strut be lined with moderately firm EVA foam to smooth the interior surface.




Figure 6




Figure 7