

A white line drawing of a mechanical elbow component is centered on a dark teal background. The drawing shows a cylindrical body with a hexagonal nut on the front face, a flange with several screws, and a textured top section. A hose or cable is attached to the side.

E-200 Elbow

Product Manual

Fillauer®

Contents

Intended Use.....	4
Warnings and Precautions	4
Qualified Provider	5
Specifications and Preparations Before Use (Risk Management for Installation and Calibration)	5
Elbow Assembly	7
Compatibility.....	10
Disposal / Waste Handling	10
Warranty	10
User Instructions	11
Serious Incidents	11



To see other language options, visit [fillauer.com](https://www.fillauer.com).

Intended Use

The E-200 Series elbows are used with any transhumeral or higher upper extremity prosthesis. These elbows are appropriately sized for smaller adults and adolescents. These elbow units feature multiple body-powered locking positions and come in a variety of configurations for color and strength. All models feature adjustable-friction axial rotation through the proximal turntable, and the standard model can be purchased with or without a [forearm section](#) (see prefabricated forearm table).

Performance Characteristics

	E-200 Elbow	E-200HD Elbow
Weight	12 oz. (340 g)	12.5 oz. (369 g)
Yoke	Aluminum	Aluminum
Locking Positions	11	9
Saddle	Single Thickness Stainless Steel	Double Thickness Stainless Steel
Primary Materials	Stainless Steel, Steel, Aluminum	Stainless Steel, Steel, Aluminum

These devices are intended for single patient use only.

Storage and Handling

It is recommended that prosthetic elbows be stored in a cool, clean, dry environment away from harsh chemicals (chlorine, acids, acetone, etc.).

Warnings and Precautions



NOTICE: An upper-limb prosthetic device user's ability to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations.



WARNING: Body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis. Failure to maintain tension while controlling the steering wheel could cause serious injury or death.



NOTICE: To prevent malfunction of the elbow, the instructions below on trimming the lock cable and housing must be followed. Failure to do so

will likely result in immediate or premature failure that will not be covered under warranty of the device.



CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and are not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids, abrasives, vibration, activities which would damage the biological limb, or prolonged extreme temperatures (< -5 °C or > 50 °C). Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the device with fresh water and dry immediately after exposure.

Qualified Provider

Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

Specifications and Preparations Before Use (Risk Management for Installation and Calibration)

Alignment

Prosthetic elbows should be aligned to match the natural hanging angle of the contralateral elbow and provide the best possible work envelope for a patient's specific goals. Standard alignment begins at 5 degrees of flexion and the adduction angle should match the hanging angle of the contralateral limb while also allowing midline tasks and reaching the mouth as necessary. The height of the elbow axis may match that of the other arm, but amputees often prefer a slightly shorter prosthetic side for control and weight.

Installation

- The elbow turntable may be laminated by removing the elbow and cork friction disk from the turntable. They are connected by the $\frac{3}{16}$ inch slotted nut on the threaded stud visible in the center of the turntable.

- Place the turntable on the distal end of the humeral section mold (beeswax, foam, plaster or similar). A PVA bag should be used to separate the turntable from the humeral mold if foam or plaster are used. Wax the interior surface and the distal face of the turntable.
- Pack any voids between the turntable and mold with silicone fitting gel or similar to prevent it from filling with laminate.
- Mask the turntable distal to the tie-in groove on the knurled ring surface to keep all laminate clear of the distal end.
- Laminate with the appropriate materials for durability and finish as desired by the patient, being sure to tie each structural layer into the tie-in groove in the lamination ring. Carbon fiber tape is a good choice for reinforcing the connection to the lamination ring and is commonly used for distal to proximal strips tied in with circumferential wraps.
- When reattaching the elbow unit, be sure to install the cork disc between the turntable and elbow, and place both washers under the slotted nut before delivery.
- Tighten the nut to set the proper amount of friction with a $\frac{3}{16}$ wrench.


Shortening the Cable and Housing

The elbow lock cable AND housing should, in most cases, be adjusted in length to allow smooth and consistent cycling of the lock. To shorten the cable housing, the cable must be pulled out distally and the proximal housing ferrule removed. The housing is then cut from the proximal end. NEVER remove the lock cable housing from the elbow to trim it. The housing has been permanently installed on all elbows. If removed, it must be reattached correctly (threadlocker and mechanically altering the thread), or it will thread into the elbow and cause malfunction of the lock.

The lock cable can be replaced. However, it is recommended that only Fillauer-fabricated lock cables (PN 50724) be used.


Elbow Assembly

- | | | | | | |
|-----------|-------|----------------------|------------|---------------------|-------------------------------------|
| A. | 50646 | Yoke Assembly | M. | 50722 | Bearing |
| B. | 54155 | Internal Assembly | N. | 50615 | Cork |
| C. | 50576 | Cable Lever, E-200A | O. | 50614 | Turntable, E-200 and HD |
| D. | 50724 | Cable Assembly | | 50617 | Turntable, E-200A and HD |
| E. | 50744 | Guide, E-200A | P. | 50668 | Belleville Washer |
| F. | 50455 | Cable Housing | Q. | 50667 | Lock Nut |
| G. | 50423 | Ferrule | R. | 50643 | Saddle Assembly, E-200 |
| H. | 50457 | Cable Housing Cover | | 50912 | Saddle Assembly, E-200HD |
| I. | 50637 | Gear Sector, E-200 | S. | 50729 | Washer |
| | 50638 | Gear Sector, E-200HD | T. | 50739 | Outside Washer |
| J. | 50642 | Long Sleeve | U. | 50740 | Screw |
| K. | 50640 | Shaft | V. | Elbow Cap—See Table | |
| L. | 50641 | Sleeve Short | W. | 50742 | Screw |
| | | | X. | 50427 | Anchor |
| | | | Y. | 50747 | Hanger |
| | | | Z. | 50745 | Swivel Lever, E-200A |
| | | | AA. | 50634 | Cotter Pin, E-200A |
| | | | AB. | 50725 | Cable Housing Assembly with Ferrule |




A. 50646 Yoke Assembly

- AC. 50675 Nut
- AD. 50674 Bumper
- AE. 50621 Yoke




B. 54155 Internal Assembly

- AF. 50647 Internal Cage Assembly
- AG. 50680 Bolt
- AH. 50636 Locking Bar Spring
- AI. 50624 Base
- AJ. 50679 Screw



AF. 50647 Internal Cage Assembly

- AK. 50625 Cage
- AL. 50635 Locking Bar
- AM. 50813 Pin
- AN. 50628 Cam
- AO. 50629 Spring
- AP. 50626 Cam Bearing
- AQ. 50627 Screw
- AR. 50632 Ratchet
- AS. 50648 Lever Assembly




AS. 50648 Lever Assembly

AT. 50630 Pawl

AU. 50705 Spring

AV. 50574 Lever

AW. 50706 Pin



Elbow Cap

Color	Outside Cable	Inside Cable	Friction Elbow
Tan	50645	54437	54436
Light Brown	54435	54405	54433
Medium Brown	54432	54431	54430
Dark Brown	54429	54428	54427
Jet Black	54434	54406	54407

Compatibility

An Elbow Lift Assist (PN 50752 for left and 50609 for right) may be attached medially (standard) or laterally. For lateral placement, be sure to order the opposite side part number. A Nudge Control (PN 52521 for right and 52522 for left) may be used to lock and unlock the elbow when it is not possible through the harness. Lock Cables may be replaced (PN 50724). **Prefabricated Forearms** (see table below) may be used with all E-200 series elbows. The Northwestern Lift Tab Jig (PN 51040) may be used with E-400 and E-200 Series elbows to properly locate the lift tab assemblies.

Prefabricated Forearms

	Tan	Light Brown	Dark Brown	Jet Black
Medium for 1 ¾ in. (45 mm) Wrist	57940	57545	57548	57561
Medium for 1 ¾ in. (45 mm) Wrist with E-200 Elbow	57541	57504	57507	57566

Disposal / Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

All metal components may be removed and recycled at the appropriate recycling facility.

Warranty

This product has a 12-month warranty against manufacturer defects

User Instructions

The providing health care professional must review the following information directly with the user.

Warnings and Precautions for the User



NOTICE: The user should monitor their prosthesis daily and contact their health care professional if they experience changes in device performance or if it begins to make noise.



CAUTION: All maintenance should be performed by the qualified health care professional.



NOTICE: An upper-limb prosthetic device user's ability to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations.



WARNING: Body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis. Failure to maintain tension while controlling the steering wheel could cause serious injury or death.



CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and are not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids, abrasives, vibration, activities which would damage the biological limb, or prolonged extreme temperatures (< -5 °C or > 50 °C). Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the device with fresh water and dry immediately after exposure.

Serious Incidents

In the unlikely event of a serious incident, seek immediate medical help and contact your prosthetist at your earliest possible convenience. Clinicians should contact their local Fillauer representative immediately in the event of any device failure.

Fillauer®

www.fillauer.com



Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946



Fillauer Europe

Kung Hans väg 2
192 68 Sollentuna, Sweden
+46 (0)8 505 332 00



© 2021 Fillauer LLC
M092/11-03-21/Rev.1