

# make way for **CHANGE**

Fillauer function & quality in a new modular design

External Shuttle Lock

External Clutch Lock

Razor Lock

## External LOCKS by Fillauer



- Interchangeable with any 36mm thread
- Titanium core pyramid rated up to 300 lbs
- Low-profile design
- High-strength aluminum body
- Lamination kits available
- Available in three different models



## Product Description

# External Locks

by Fillauer LLC

## External Shuttle Lock

125100



## External Clutch Lock

140100



## Razor Lock

125465



## Fillauer External Lock System

Get the Fillauer function and quality in three new external modular designs. Fillauer has integrated its well-known clutch and shuttle assemblies in three external configurations, suitable for use in any 36mm threaded adapter. Using high strength titanium cores, Fillauer offers a 300 lb weight rating for high performance for each of the devices. The external shuttle lock uses the same simple yet rugged tool steel lock insert found in other Fillauer shuttle designs. The clutch insert has greater density of engagement teeth and a left-handed thread to avoid loosening. The adaptable Razor lock is a ultra-thin nickel plated design with through holes to allow any number of distal attachments.

## Product Features

- Interchangeable with any 36mm thread
- Titanium core pyramid rated up to 300 lbs
- Low-profile design
- High-strength aluminum body
- Lamination kits available
- Available in three different models

## Indications

- External 36mm threaded attachments
- Patients requiring 300 lb weight rating
- For use with Fillauer Shuttle or Clutch Pin Clutch

## External Shuttle Lock

<b>125100</b>	External Shuttle Lock
<b>809725</b>	Plunger 1 1/2", with 1/4" Thread
<b>809727</b>	Plunger 2", with 1/4" Thread
<b>809725mm</b>	Plunger 1 1/2", with M 10 Metric Thread
<b>809727mm</b>	Plunger 2", with M 10 Metric Thread
<b>809720mm</b>	Plunger 2 1/2", with M 10 Metric Thread
<b>809721mm</b>	Plunger 3", with M 10 Metric Thread
<b>809880</b>	Fabrication Kit for External Locks

## External Clutch Lock

<b>140100</b>	External Clutch Lock
<b>140040</b>	2" Clutch Lock Plunger with 1/4-20 Thread
<b>140041</b>	1 1/2" Clutch Lock Plunger with 1/4-20 Thread
<b>140043</b>	2" Clutch Lock Plunger with M10 Metric Thread
<b>140044</b>	1 1/2" Clutch Lock Plunger with M10 Metric Thread
<b>809880</b>	Fabrication Kit for External Locks

## Razor Lock

<b>125465</b>	Razor Lock
<b>125293</b>	Lamination Ring with SHCS
<b>809725</b>	Plunger 1 1/2", with 1/4" Thread
<b>809727</b>	Plunger 2", with 1/4" Thread
<b>809725mm</b>	Plunger 1 1/2", with M 10 Metric Thread
<b>809727mm</b>	Plunger 2", with M 10 Metric Thread
<b>809720mm</b>	Plunger 2 1/2", with M 10 Metric Thread
<b>809721mm</b>	Plunger 3", with M 10 Metric Thread
<b>125205</b>	Fabrication Kit for Lamination

## L-Code Information

<b>L5671</b>	Below Knee/Above Knee Suspension Locking Mechanism
--------------	--

*The following is the suggested L-Code for the Fillauer External Locks. The manufacturer and distributor of the External Locks do not assume any liability for L-Code recommendations.*