



SPL2

from Fillauer LLC

The most popular
stance control
orthotic knee joint
just got better



- Enhanced position dependent pendulum design
- Multiple locking positions for tentative gait
- Added free motion feature
- Easy fabrication with minimal follow-up maintenance
- Works with any orthotic ankle design
- Only design available for use as a rehab knee orthoses
- Bilateral or unilateral uprights for easy donning

*Available in a kit or as a Fillauer
Custom Fabrication orthosis*

U.S. Patent
6979304



Fillauer LLC

fillauer.com/products/spl2

1.800.251.6398

AD278/07-09/04-20-10

Swing Phase Lock 2 Orthotic Knee Joint Ordering Guide

INDICATIONS

- Post-polio
- Spinal involvement
- Cerebral vascular accident
- Peripheral paresis/paralysis
- Nerve inflammations
- Neurological failures
- Myopathies
- Multiple sclerosis (or similar diseases)

CONTRAINDICATIONS

- Knee flexion contractures > 10°
- Central paresis/paralysis
- Hip flexion contracture
- Weight bearing orthotic applications

WARRANTY


Fillauer orthotic components are warranted for 12 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor. Items under warranty will be replaced or repaired (at Fillauer's discretion) at no charge. The warranty will be void if the item has been fabricated or installed outside Fillauer's recommendations, if the item has been exposed to a corrosive environment, or if the item has been used in extremely abusive activities that could result in injury.

SWING PHASE LOCK KITS

IQ Swing Phase Lock 2 Kit includes the SPL2 Joint, SPC Joint, and Satellite (remote control)

| | |
|---------------|------------------------------------|
| IQ200L | IQ Swing Phase Lock 2 Kit, Left |
| IQ250L | IQ SPL2 Joint, Left |
| IQ120L | IQ SPC Joint, Left |
| IQ140 | IQ SPL Satellite (remote control) |
| IQ050L | IQ Single Sided Swing Phase Lock 2 |
| IQ250L | IQ SPL2 Joint, Left |
| IQ140 | IQ SPL Satellite (remote control) |
| IQ200R | IQ Swing Phase Lock 2 Kit, Right |
| IQ250R | IQ SPL2 Joint, Right |
| IQ120R | IQ SPC Joint, Right |
| IQ140 | IQ SPL Satellite (remote control) |
| IQ050R | IQ Single Sided Swing Phase Lock 2 |
| IQ250R | IQ SPL2 Joint, Right |
| IQ140 | IQ SPL Satellite (remote control) |

NEW FEATURES OF THE SPL2



SWING PHASE LOCK KITS

Fabrication Kits

| | |
|---------------|---------------------------------------|
| IQ150 | Alignment Kit |
| 023516 | Swing Lock Upright Kit, AL 3/16 x 3/4 |
| 023616 | Swing Lock Upright Kit, SS 3/16 x 3/4 |

Polycar-C™

| | |
|---------------|-------------------------------------|
| 076630 | PolyCar-C™ Insert, .060 Small, pair |
| 076632 | PolyCar-C™ Insert, .090 Med, pair |
| 076634 | PolyCar-C™ Insert, .090 Large, pair |
| 076636 | PolyCar-C™ Insert, .125 XL, pair |
| 020803 | Pre-cut PolyCar-C™, Large |

Carbon also available

Ankle Joints

| | |
|---------------|--|
| 020886 | Double Action, SS, Large, 3/4 channel |
| 020888 | Double Action w/Positioning Pins, SS, Large, 3/4 channel |
| 020882 | Double Action, SS, Large, 5/8 channel |
| 020884 | Double Action w/Positioning Pins, SS, Large, 5/8 channel |

Insert Stirrups for Plastic Orthotics

| | |
|---------------|---|
| 020099 | Insert Stirrup, Limited Motion, 2-5/16" |
| 020115 | Insert Stirrup, Limited Motion, 2-3/4" |
| 020123 | Insert Stirrup, Toe Lift, 2-5/16" |
| 020149 | Insert Stirrup, Toe Lift, 2-3/4" |
| 020094 | Insert Stirrup, Double Action, 2-3/4" |
| 020016 | Offset Stirrup, Limited, 2-1/4" x 3/8" |
| 020024 | Offset Stirrup, Limited, 2-1/4" x 1/2" |
| 020032 | Offset Stirrup, Limited, 2-1/2" x 3/8" |
| 020040 | Offset Stirrup, Limited, 2-1/2" x 1/2" |
| 020057 | Offset Stirrup, Toe Lift, 2-1/4" x 3/8" |
| 020065 | Offset Stirrup, Toe Lift, 2-1/4" x 1/2" |
| 020073 | Offset Stirrup, Toe Lift, 2-1/2" x 3/8" |
| 020081 | Offset Stirrup, Toe Lift, 2-1/2" x 1/2" |