

# Tag... You're it.

*Fillauer*  
COMPANIES

## **MIGHTYMITE<sup>®</sup>** **PEDIATRIC** **KNEE** from Hosmer

- Rugged, compact, rated to 132 lbs
- Length is 1-5/8" from distal end of socket to distal knee cap
- Gliding, easy to flex, polycentric design with friction control
- 4 different proximal attachments
- 3 different extension assist springs



**FC100510**  
**FC100500**

## MIGHTYMITE® KNEE OPTIONS

	<p><b>FC100500</b> MightyMite® 4-Bar Knee</p>		<p><b>FC100510</b> MightyMite® 4-Bar Knee with Proximal Pyramid Attachment</p>
---	---	--	--

## MIGHTYMITE® KNEE ADAPTERS

	<p><b>FC100501</b> MightyMite® Pyramid Adapter</p>		<p><b>FC100504</b> MightyMite® Threaded Adapter</p>
	<p><b>FC100503</b> MightyMite® Tapered 4-Hole Adapter</p>		<p><b>FC100130</b> MightyMite® AK Socket Adapter</p>
	<p><b>FC100502</b> MightyMite® Slide Adapter</p>		<p><b>FC100134</b> MightyMite® Medium Pyramid Receiver</p>
	<p><b>FC100505</b> MightyMite® Male Dovetail 4-Hole Adapter</p>		<p><b>FC100132</b> MightyMite® 4-Hole Adapter with Rotation</p>

## L-CODES\*

**L5611** 4-Bar Linkage with Friction Control

**L5920** Alignable System

*\*Suggested L-Codes are provided as a reference only. It is the responsibility of the practitioner to confirm this information.*

PRODUCT DESCRIPTION	WEIGHT	BUILD HEIGHT
MightyMite® 4-Bar Knee	270 g / 9.55 oz	63 mm / 3-3/4 in
MightyMite® 4-Bar Knee with Proximal Pyramid	294 g / 10.35 oz	70 mm / 2-1/2 in
MightyMite® 4-Bar Knee with Slide Adapter	323 g / 11.40 oz	84 mm / 3-3/8 in
MightyMite® 4-Bar Knee with Tapered 4-Hole Adapter	310 g / 11.07 oz	75 mm / 2-15/16 in
MightyMite® 4-Bar Knee with Threaded Adapter	288 g / 10.15 oz	67 mm / 2-5/8 in