



**CYLINDRICAL
LOCK**
from Fillauer LLC

- Dependable Fillauer shuttle design
- Retrofits other cylindrical lock sockets
- Strong, lightweight, fiber reinforced body
- Standard 4-hole attachment pattern
- Adjustable push button
- Includes 2 Fillauer 1 in. or 1-1/2 in. plungers



CYLINDRICAL LOCK

from Fillauer LLC



The Cylindrical Lock incorporates the same dependable Fillauer locking mechanism in a versatile cylindrical body with a robust 4-hole attachment base. This value-added design features a fiber reinforced, injection molded body, and deep pin well for easy insertion. The Fillauer Cylindrical Lock can be retrofitted to other cylindrical lock sockets and is optimized for thermoforming or lamination. The Cylindrical Lock has been designed to not bind like other designs when it is torqued. With its focus on value, the Cylindrical Lock comes ready to fit with two 1" Fillauer shuttle pins and attachment hardware. Great value, trusted quality.

FEATURES & BENEFITS

- Dependable Fillauer shuttle design
- Retrofits other cylindrical lock sockets
- Strong, lightweight, fiber reinforced body
- Will not bind when torqued
- Cost effective pricing
- Majority gel suspension liners
- Serrated pin type attachment
- Standard distal end profile
- Includes two Fillauer 1 in. plungers

SPECIFICATIONS

- Standard 4-hole attachment pattern
- Weight 0.193 lbs. / 88 g
- Rated up to 300 lbs.

CONTRAINDICATIONS

- Narrow distal end
- Previous pin marring issues
- Long limb lengths
- Offset alignment

SHUTTLE LOCK ASSEMBLY

125240	Cylindrical Lock, with 1 in. Plungers 809722mm, 1 in. Long Metric Plunger, 2 Each
125241	Cylindrical Lock, with 1-1/2 in. Plungers 809725mm, 1-1/2 in. Long Metric Plunger, 2 Each
125243	Cylindrical Lock, with 2 in. Plungers 809727mm, 2 in. Long Metric Plunger, 2 Each
809685	Latch Pin Button for Lock Systems
809873	Cylindrical Lock Body
880285	M6 x 1 x 25 mm Flat Socket Head Cap Screw, 4 Each

SHUTTLE LOCK ASSEMBLY (SOLD SEPARATELY)

809867	Cylindrical Housing Dummy Kit
809864	Fabrication Dummy Body
880275	M6-1x16 mm Socket Head Set Screw
880494	5/16-18x3 in. Socket Head Set Screw
809791	Button Shield Dummy

PLUNGERS (SOLD SEPARATELY)

809722	Plunger 1 in., with 1/4 in. Thread
809725	Plunger 1-1/2 in., with 1/4 in. Thread
809727	Plunger 2 in., with 1/4 in. Thread
809722mm	Plunger 1 in., with M10 Metric Thread
809725mm	Plunger 1-1/2 in., with M10 Metric Thread
809727mm	Plunger 2 in., with M10 Metric Thread
809720mm	Plunger 2-1/2 in., with M10 Metric Thread
809721mm	Plunger 3 in., with M10 Metric Thread

May also use fabrication kits 125204 and 125202 with this item