



Feel good in your skin

*Fillauer*  
COMPANIES

**DREAM  
SKIN**  
from Centri

- Micro Coated Vinyl reduces clothing friction
- Water, dirt, and stain resistant
- Cosmetic cover with a more natural appearance
- Available in solid or split toe
- Offered in 5 sizes and 18 color shades



# DREAM SKIN

from Centri



Centri Dream Skin is a protective leg covering that combines superior durability and a stain-resistant cosmesis with the elastic properties necessary for easy donning and application. Using our new molding technology, Dream Skin can offer the life-like cosmetic quality of a custom cosmesis, down to the micro finish details of skin texture, coloring, and toenails, at the price of a standard production covering.

The secret is a new lightweight vinyl material which is made tougher than other designs without sacrificing cosmetic quality. Dream Skin also integrates a special micro-coating that prevents most permanent stains and helps clothing slide easily over the prosthesis when dressing.

## PRODUCT FEATURES

- Micro coated vinyl resists marks and reduces clothing friction
- Cosmetic cover with a more natural appearance
- Superior durability
- Water resistant
- Available in 18 different skin colors (use Centri color guide)
- Expanded for easy donning

## WARRANTY

- 6 month warranty

## ORDERING INFORMATION

Choose from these color options (xx). Due to display variances colors may not be exact. Please reference Swatch Color Guide (PN 611200002).

01	02	03	Split Toe	Without Split Toe	Side	Foot Length	Calf Size Min.	Ankle Size Min.	Length for Split Toe	Length for Non Split Toe
04	05	06	421223xx4	421221xx4	Left	8.27–8.66 in.	12.01 in.	7.48 in.	19.69 in.	19.69 in.
07	08	09	421224xx4	421222xx4	Right	21–22 cm	30.5 cm	19.0 cm	50.0 cm	50.0 cm
10	11	12	423243xx4	423241xx4	Left	9.06–9.45 in.	12.60 in.	7.87 in.	21.65 in.	21.65 in.
13	14	15	423244xx4	423242xx4	Right	23–24 cm	32.0 cm	20.0 cm	55.0 cm	55.0 cm
16	17	18	425263xx4	425261xx4	Left	9.84–10.24 in.	13.39 in.	8.27 in.	23.62 in.	23.62 in.
			425264xx4	425262xx4	Right	25–26 cm	34.0 cm	21.0 cm	60.0 cm	60.0 cm
			427283xx4	427281xx4	Left	10.63–11.02 in.	13.98 in.	8.46 in.	23.62 in.	21.65 in.
			427284xx4	427282xx4	Right	27–28 cm	35.5 cm	21.5 cm	60.0 cm	55.0 cm
			429303xx4	429301xx4	Left	11.42–11.81 in.	14.57 in.	8.66 in.	25.59 in.	25.59 in.
			429304xx4	429302xx4	Right	29–30 cm	37.0 cm	22.0 cm	65.0 cm	65.0 cm