



Fillauer
COMPANIES

ABSOLUTE SUSPENSION SLEEVE

from Centri

- Preflexed suspension sleeve
- Maximum comfort
- Slides against clothing
- Available in extra soft
- Available in black or beige fabric
- Made of thermoplastic elastomer and covered with a reinforced Lycra® type material. The inside is also reinforced with Lycra®.



Centri • +46 8 505 332 00 • www.centri.se Fillauer LLC • 1.800.251.6398 • www.fillauer.com Hosmer • 1.800.827.0070 • www.hosmer.com

US Pat. 6,406,499, European Pat. 0955964, Canadian Pat. 2,190,764. AD240-02/08-06/06/13

ABSOLUTE SUSPENSION SLEEVE

from Centri

Centri's Absolute Suspension Sleeve optimizes comfort, durability, and range of motion. The sleeve is made of thermoplastic elastomer and covered with a reinforced Lycra® type material. The preflexed design allows maximum range of motion and alleviates posterior bunching when the knee is flexed. The Absolute Suspension Sleeve is available in multiple sizes and your choice of beige or black.

PRODUCT FEATURES

- Preflexed suspension sleeve
- Maximum comfort
- Slides against clothing
- Available in extra soft
- Available in black or beige fabric
- Made of Thermoplastic Elastomer and covered with a reinforced Lycra® type material. The inside is also reinforced with Lycra®.



US Pat. 6,406,499, European Pat. 0955964, Canadian Pat. 2,190,764
Licensed from Ohio Willow Wood

ORDERING INFORMATION

| Absolute Preflexed Suspension Sleeve, 3mm - Extra Soft | | | | |
|--|-------------------------|-------------------------|-----------|-----------|
| Size | Distal Circumference | Proximal Circumference* | Beige | Black |
| Small | 21-25 cm (08.3-09.8 in) | 26-31 cm (10.2-12.2 in) | 803471314 | 803461314 |
| Medium | 26-32 cm (10.2-12.6 in) | 31-38 cm (12.2-15.0 in) | 803471318 | 803461318 |
| Medium+ | 30-36 cm (11.8-14.2 in) | 38-44 cm (15.0-17.3 in) | 803471322 | 803461322 |
| Large | 32-40 cm (12.6-15.8 in) | 44-50 cm (17.3-19.7 in) | 803471326 | 803461326 |
| Large+ | 34-42 cm (13.4-16.5 in) | 50-60 cm (19.7-23.6 in) | 803471332 | 803461332 |
| X-Large | 40-49 cm (15.7-19.3 in) | 60-70 cm (23.6-27.5 in) | 803471340 | 803461340 |

| Absolute Straight Suspension Sleeve 3mm - Extra Soft | | | | |
|--|-------------------------|-------------------------|-----------|-----------|
| Size | Distal Circumference | Proximal Circumference* | Beige | Black |
| X-Small | 10-17 cm (3.9-6.7 in) | 16-24 cm (6.3-9.4 in) | 803371301 | 803361301 |
| Small | 15-22 cm (5.9- 8.7 in) | 23-34 cm (9.1-13.4 in) | 803371302 | 803361302 |
| Medium | 21-30 cm (8.3-11.8 in) | 29-43 cm (11.4-16.9 in) | 803371303 | 803361303 |
| Medium+ | 24-35 cm (9.4-13.8 in) | 33-47 cm (13-18.5 in) | 803371304 | 803361304 |
| Large | 28-42 cm (11-16.5 in) | 36-51 cm (14.2-20.1 in) | 803371305 | 803361305 |
| Large+ | 31-46 cm (12.2-18.1 in) | 39-55 cm (15.4-21.7 in) | 803371306 | 803361306 |
| X-Large | 34-50 cm (13.4-19.7 in) | 42-60 cm (16.5-23.6 in) | 803371307 | 803361307 |
| XX-Large | 36-55 cm (14.2-21.7 in) | 44-70 cm (17.3-27.6 in) | 803371308 | 803361308 |



*NOTE: The proximal circumference is measured 15 cm above the mid-patella tendon.