


smooth lines,
smooth rollover...



element
FOOT SYSTEM

 reddot design award
winner 2010

- > Unique 3rd Element technology
- > Innovative Micro-Slice technology for multi axial performance
- > Smooth rollover and superior stability
- > Very low profile - only 3.75 in. clearance (9.5 cm)
- > Rated for patients up to 300 lbs. (136 kg)
- > Durable; meets the ISO-22675 test standard
- > 3 year warranty/30 day patient trial



Patent
Pending

*e*MOTIS
A Fillauer Company

Fillauer LLC • Worldwide Distributor • 800.251.6398 • 423.624.0946 • www.fillauer.com

AD267/09-08-11/06-04-13

ELEMENT

from Emotis



The Element Foot sets itself apart from two spring prosthetic foot designs by adding an innovative 3rd Element. Emotis has designed the Element to support each phase of foot function in the gait cycle—from heel strike, to mid-stance, to toe off. The result is a foot with unsurpassed smoothness, stability and performance throughout the gait cycle.

FEATURES & BENEFITS

- Unique three Element design
- Smooth rollover and superior stability
- Adjustable heel wedge for fine tuning stiffness
- Split heel for multi-axial inversion and eversion
- Innovative Micro-Slice technology for multi-axial performance

SPECIFICATIONS

- Patient weight: Up to 300 lbs. (136 kg)
- Foot weight: 12 oz. (340 g)
- Build height: 3.75 in. (9.5 cm)
- Functional level: K3–K4
- Durable; meets ISO-22675 standard

L-CODES*

L5981 All lower extremity prostheses, flex-walk system or equal.

L5986 Multi-Axial Rotation Unit (“MCP” or Equal)

***Suggested L-Codes are provided as a reference only. It is the responsibility of the practitioner and/or patient care facility submitting this claim to determine and submit the correct coding within all established rules and guidelines.**

WARRANTY

- 3 years from the date of shipment to the practitioner
- Foot Shell: 6 months from the date of shipment to the practitioner

SATISFACTION GUARANTEE

- 30 days from the date of shipment to the practitioner

ORDERING INFORMATION

100-10-XX-XXX Element Foot System

Foot includes Spectra Sock. MCV Foot Shell sold separately. Foot selection is determined by foot size and patient weight. The foot model recommendations on the selection chart are based on average patients, engaged in normal activities. If your patient is much more active than average or is lifting heavy loads as part of their daily routine, then going up one weight category is advised.*

		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100-119 lbs.	45-54 kg	22-A11	23-A11	24-A11	*	*	*	*	*	*
120-139 lbs.	54-63 kg	22-A11	23-A11	24-A11	25-A32	26-A32	*	*	*	*
140-159 lbs.	64-72 kg	22-A11	23-A11	24-A11	25-A32	26-A32	27-A43	28-A43	*	*
160-179 lbs.	73-81 kg	22-A22	23-A22	24-A22	25-A43	26-A43	27-A43	28-A43	29-C64	30-C64
180-199 lbs.	82-90 kg	22-A22	23-A22	24-A22	25-A43	26-A43	27-B54	28-B54	29-C64	30-C64
200-224 lbs.	91-102 kg	22-A32	23-A32	24-A32	25-B54	26-B54	27-B54	28-B54	29-D75	30-D75
225-249 lbs.	102-113 kg	NA	NA	NA	25-B54	26-B54	27-C65	28-C65	29-D75	30-D75
250-275 lbs.	113-125 kg	NA	NA	NA	25-C65	26-C65	27-C65	28-C65	29-D76	30-D76
276-300 lbs.	125-136 kg	NA	NA	NA	25-C65	26-C65	27-D65	28-D65	29-D76	30-D76

ACCESSORIES

XX XX 13 CC 3 Micro Coated Vinyl Foot Shell

Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24–30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30

Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.

Spectra Sock

192500 Spectra Sock* (sold separately)

*Recommended option when using foot shell.

