



Dynamic Walk

Fillauer®



Reimagining Possibilities

The **Dynamic Walk AFO** is designed to allow people with dorsiflexor paresis or paralysis to confidently regain their stride. Whether their goal is to walk, hike, jog, or play, the Dynamic Walk's flexible dorsi-assist joints allow them to move freely. Featuring a heat-moldable, trimmable calf band and foot plate, the Dynamic Walk ensures a comfortable fit every time.

Fillauer®



Thermoformable
Carbon Composite

Flexible
Dorsi-assist
Joint

Narrow Foot Plate
Fits in Most Shoes



Dynamic Walk, Standard

- 700111121** Dynamic Walk, Small, Left
- 700111122** Dynamic Walk, Small, Right
- 700111141** Dynamic Walk, Medium, Left
- 700111142** Dynamic Walk, Medium, Right
- 700111161** Dynamic Walk, Large, Left
- 700111162** Dynamic Walk, Large, Right
- 700100100** Custom Fabricated*

Dynamic Walk, Wide Ankle

- 700111321** Dynamic Walk Wide Ankle, Small, Left
- 700111322** Dynamic Walk Wide Ankle, Small, Right
- 700111341** Dynamic Walk Wide Ankle, Medium, Left
- 700111342** Dynamic Walk Wide Ankle, Medium, Right
- 700111361** Dynamic Walk Wide Ankle, Large, Left
- 700111362** Dynamic Walk Wide Ankle, Large, Right
- 700100100** Custom Fabricated*

Padding Kit, Standard & Wide Ankle

- 700100220** Small
- 700100240** Medium
- 700100260** Large



Padding Kit, Single Side

- 700100321** Small, Left Lateral or Right Medial
- 700100322** Small, Right Lateral or Left Medial
- 700100341** Medium, Left Lateral or Right Medial
- 700100342** Medium, Right Lateral or Left Medial
- 700100361** Large, Left Lateral or Right Medial
- 700100362** Large, Right Lateral or Left Medial



Warranty

1 year warranty on AFO
3 month warranty on padding kit



Dynamic Walk, Single Side, Medial

- 700113121** Dynamic Walk, Medial, Small, Left
- 700113122** Dynamic Walk, Medial, Small, Right
- 700113141** Dynamic Walk, Medial, Medium, Left
- 700113142** Dynamic Walk, Medial, Medium, Right
- 700113161** Dynamic Walk, Medial, Large, Left
- 700113162** Dynamic Walk, Medial, Large, Right
- 700100100** Custom Fabricated*

Dynamic Walk, Single Side, Lateral

- 700114121** Dynamic Walk, Lateral, Small, Left
- 700114122** Dynamic Walk, Lateral, Small, Right
- 700114141** Dynamic Walk, Lateral, Medium, Left
- 700114142** Dynamic Walk, Lateral, Medium, Right
- 700114161** Dynamic Walk, Lateral, Large, Left
- 700114162** Dynamic Walk, Lateral, Large, Right
- 700100100** Custom Fabricated*

Custom-Fit Selection Guide

Shoe Size		Small	Medium	Large
U.S. & Canada	M	4 - 6.5	7 - 10	10.5 - 14
	W	5.5 - 8	8.5 - 11.5	12 - 15.5
U.K.	M	3.5 - 6	6.5 - 9.5	10 - 13.5
	W	3 - 5.5	6 - 9	9.5 - 13
Europe		35.5 - 38.5	39 - 43.5	44 - 48.5
Mexico		4.5 - 5	5.5 - 8.5	9 - 12.5
Japan	M	22 - 24.5	25 - 28	28.5 - 31.5
	W	21.5 - 24	24.5 - 27.5	28 - 31
Australia	M	3.5 - 6	6.5 - 9.5	10 - 13.5
	W	4 - 6.5	7 - 10	10.5 - 14
Korea (mm)		231 - 250	251 - 272	273 - 292
Length		9.1 in. (23 cm)	9.8 in. (25 cm)	10.6 in. (27 cm)
Height		12.6 in. (32 cm)	13.4 in. (34 cm)	14.6 in. (37 cm)
Regular Ankle Width		2.7 in. (6.8 cm)	2.9 in. (7.3 cm)	3.2 in. (8.1 cm)
Wide Ankle Width		2.9 in. (7.3 cm)	3.2 in. (8.1 cm)	3.5 in. (8.9 cm)

*See Custom Fabrication Order Form on back.

Fillauer[®]

www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936

© 2016 Fillauer LLC
AD364/01-28-16