



Patent
Pending



eMOTIS
A Fillauer Company

ibex™

FOOT SYSTEM

- > Innovative Micro-slice pylon
- > Multi-axis heel/toe for terrain conformance
- > Energy optimized pylon/heel plate design
- > Designed for patients up to 330 lbs.
- > Ibex XD rated for patients up to 440 lbs.
- > Durable; meets the ISO-22675 test standard
- > 3 year warranty/30 day patient trial

IBEX

from Fillauer LLC

The Ibx Foot System is a leap forward in multi-axial foot design. The innovative micro-slices in the Ibx pylon and the split heel plate provide controlled inversion/eversion for terrain conformance. The multi-axial design is energy optimized with a long carbon pylon and full length heel plate that allows the Ibx to reach foot flat sooner and store more energy. The unique components of the Ibx work together from heel strike to toe off to provide stability without sacrificing energy return. The Ibx foot's balanced performance gives amputees the confidence they need to meet life's challenges head on.



PRODUCT FEATURES

- Innovative Micro-Slice technology for multi-axial performance
- Energy optimized pylon/heel plate design
- Carbon fiber foot/shank design for outstanding energy return
- Adjustable heel wedge for fine tuning heel stiffness

IBEX FOOT SYSTEM

- Patient weight: Up to 330 lbs. (150 kg)
- Foot weight: 16.2 oz. (459 g)
- Build height: 6.75 in. (17 cm)
- Functional level: K3-K4
- Durable; meets ISO-22675 standard
- 3 year warranty/30 day patient trial

IBEX XD FOOT SYSTEM

- Patient weight: Up to 440 lbs. (200 kg)
- Foot weight: 21.2 oz. (601 g)
- Build Height: 7.63 in. (19 cm)

L-CODES

- L5980** All Lower Extremity Prostheses, Flex Foot System
L5986 All Lower Extremity Prostheses, Multi-Axial Rotation Unit ("MCP" or Equal)

Suggested L-Codes are provided as a reference only. It is the responsibility of the practitioner and/or patient care facility submitting this claim to determine and submit the correct coding within all established rules and guidelines.

WARRANTY

- 36 Months from the date of shipment to the practitioner
- Foot Shell—6 months from the date of shipment to the practitioner

SATISFACTION GUARANTEE

- 30 days from the date of shipment to the practitioner

PURCHASE INFORMATION

120-10-XX-XX(XX*) Ibx Foot System

Kit includes foot, elastomer heel wedges, Spectra Sock, and instruction manual. Foot selection requires foot length and weight. *Add "XD" for patients 330-440 lbs. as described in chart below.

		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100–119 lbs.	45–54 kg	22-A2	23-A2	24-A2	25-B4	26-B4	*	*	*	*
120–139 lbs.	54–63 kg	22-A2	23-A2	24-A2	25-B4	26-B4	27-C5	28-C5	*	*
140–159 lbs.	64–72 kg	22-B3	23-B3	24-B3	25-B5	26-B5	27-C5	28-C5	29-D6	30-D6
160–179 lbs.	73–81 kg	22-B3	23-B3	24-B3	25-B5	26-B5	27-C6	28-C6	29-D6	30-D6
180–199 lbs.	82–90 kg	22-B4	23-B4	24-B4	25-C6	26-C6	27-C6	28-C6	29-D7	30-D7
200–224 lbs.	91–102 kg	22-B4	23-B4	24-B4	25-C6	26-C6	27-D7	28-D7	29-D7	30-D7
225–249 lbs.	102–113 kg	NA	NA	NA	25-D7	26-D7	27-D7	28-D7	29-E8	30-E8
250–275 lbs.	113–125 kg	NA	NA	NA	25-D7	26-D7	27-E8	28-E8	29-E8	30-E8
276–300 lbs.	125–136 kg	NA	NA	NA	25-E8	26-E8	27-E8	28-E8	29-E9	30-E9
301–330 lbs.	137–150 kg	NA	NA	NA	25-E8	26-E8	27-E9	28-E9	29-E9	30-E9
330–440 lbs.	150–200 kg	NA	NA	NA	25-E8XD	26-E8XD	27-E9XD	28-E9XD	29-E9XD	30-E9XD

ACCESSORIES

XX XX 13 CC 3 Micro Coated Vinyl Foot Shell
Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24–30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30



Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.