



Step up your game

Fillauer
COMPANIES

**WAVE
SPORT**



Join us on
Facebook

Meet 2012 London Paralympic Shot Put Gold Medalist, Jackie Christiansen who won wearing his Wave Sport foot from Fillauer. The Wave Sport's full heel plate design enables patients like Jackie to wear a foot shell with any standard athletic shoe—putting the Wave Sport in a league of its own. To learn more about Jackie and the Wave Sport, visit www.fillauer.com.

Fillauer LLC • 800.251.6398 • 423.624.0946 • www.fillauer.com



Patent
Pending

WAVE SPORT

from Emotis

From extreme sports to light jogging, the Wave Sport has the performance active amputees desire. The Wave foot design uses an innovative wave spring between the pylon and heel plate to provide shock absorption and stability. The full heel plate allows it to be used in a foot shell and worn with standard athletic shoes. These features make the Wave Sport a versatile alternative to other high performance feet.



PRODUCT FEATURES

- Unique Wave-Spring heel technology
- Designed for foot shell and most athletic shoes
- Excellent shock absorption
- Full length heel for greater stability and dynamics
- “C-shaped” composite pylon for maximum energy return
- Custom built for each patient
- Modular adapters for easy alignment
- Accommodates a wide range of activity levels
- 2 year warranty/30 day patient trial

SPECIFICATIONS

- Patient weight: Up to 275 lbs. (125 kg)
- Foot weight: 19.4 oz. (550 g)
- Functional level: High K3–K4
- Durable; meets ISO-22675 standard

L-CODES*

L5987 All lower extremity prosthesis, shank foot system with vertical loading pylon.

*Suggested L-Codes are provided as a reference only, it is the responsibility of the practitioner to confirm this information.

ORDERING INFORMATION

Specify part numbers, activity level, weight, and foot size as described below.

Pylon Adapter	Product Number	Pylon Height*	Size Range	Clearance**	Adapter Options (sold separately)	
 4 Hole Plate	127-10-XX-XXXL	Low Profile	25-30	8.5 in. (23 cm)	 60254 Ti. Pyramid	 60293 Ti. Receiver
	127-10-XX-XXXS	Standard Height	22-24	9 in. (23 cm)		
	127-10-XX-XXXS	Standard Height	25-30	10 in. (25 cm)		
 36 mm Threaded	129-10-XX-XXXL	Low Profile	25-30	8.5 in. (23 cm)	 129-40-4000 Ti. Pyramid with 36 mm Threads	 129-40-4010 Ti. Receiver with 36 mm Threads
	129-10-XX-XXXS	Standard Height	22-24	9 in. (23 cm)		
	129-10-XX-XXXS	Standard Height	25-30	10 in. (25 cm)		
 Direct Mount Pyramid	129-12-XX-XXXL	Low Profile	25-30	7.75 in. (20 cm)		
	129-12-XX-XXXS	Standard Height	22-24	8.25 in. (21 cm)		
	129-12-XX-XXXS	Standard Height	25-30	9.25 in. (23 cm)		

*Practitioner should choose the standard height foot system for maximum performance as clearance allows. Consult Fillauer customer service for category selection. **Clearance includes foot shell (27 cm foot with 10 mm height). Measurement is made from mid-dome to floor.

ACCESSORIES

XX XX 13 CC 3 Micro Coated Vinyl Foot Shell
 Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24–30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30



Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.

Spectra Sock

192500 Spectra Sock* (sold separately)

*Recommended option when using foot shell.