

# Engineered for simplicity



*Fillauer*  
COMPANIES

**RAIZE**  
from Hosmer

The Raize is a unique microprocessor ankle/foot system in a low profile, lightweight design. Its plantar/dorsi-flexion range reduces damaging forces on the residual limb and greatly enhances stability on slopes and slippery surfaces. By simply pressing a button the wearer can control rollover resistance, lock or unlock the ankle, and adjust the heel height as needed. Visit us online at [www.fillauer.com](http://www.fillauer.com) to see the Raize in action.



Hosmer • 1.800.827.0070 • [www.hosmer.com](http://www.hosmer.com) Fillauer LLC • 1.800.251.6398 • [www.fillauer.com](http://www.fillauer.com)

AD323/10-11-11/07-12-12

# RAIZE

from Hosmer

## FEATURES & BENEFITS

- Microprocessor controlled hydraulic ankle/composite foot system
- Adjustable heel height and lock via control pod or remote fob
- Adjustable plantar/dorsiflexion range and resistance via control pod or remote fob
- Terrain accommodation mode aids gait on inclines and declines
- Simulates natural ankle movement and reduces shock
- Remote ankle lock for driving and donning shoes and socks
- Intuitive computer interface for initial setup by prosthetist

## ORDERING INFORMATION

### Raize™ Ankle/Foot System

Kit includes: Ankle/Foot, Software, Control Unit, Remote Fob, Battery, Charger, USB Cable, and Spectra Sock.

		24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100-130 lbs.	45-59 kg.	62869	62872	62876	62880	62884	62917	*
131-160 lbs.	60-72 kg.	62870	62873	62877	62881	62885	62918	62921
161-190 lbs.	73-86 kg.	62871	62874	62878	62882	62886	62919	62922
191-220 lbs.	87-100 kg.	*	62875	62879	62883	62916	62920	62923

\* Contact your distributor for recommendations.

## ACCESSORIES

**XX XX 13 CC 3** Micro Coated Vinyl Foot Shell  
 Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24–30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30



Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.

### Spectra Sock

**192500** Spectra Sock\* (sold separately)

\*Recommended option when using foot shell.

## L-CODES\*

**L5973** Endoskeletal ankle foot system, microprocessor controlled feature, dorsiflexion and/or plantar flexion control, includes power source.

\*Suggested L-Codes are provided as a reference only, it is the responsibility of the practitioner to confirm this information.

## WARRANTY

- 2 year limited warranty
- 6 month foot shell warranty (sold separately)

## SPECIFICATIONS

- Weight rating: 220 lbs. (100 kg)
- Foot sizes: 24–30 cm
- Foot weight: 1.62 lbs. (735 g) to 1.758 lbs. (797 g)
- Build height (sizes 24–28): 3.5 in. (9 cm)
- Build height (sizes 29–30): 3.75 in. (9.5 cm)
- Functional level: K3

