



**Everywhere is  
walking distance**

*Fillauer*  
COMPANIES

**AERIS**  
**SOLUTION**  
from Fillauer LLC

- Compliant Composite Technology (CCT)
- Advanced Energy Release (AER) pylon
- Enhanced plantarflexion for stability at heel strike
- Balanced compliance and dynamic response
- Durable; meets ISO-22675 standard



# AERIS SOLUTION

from Fillauer LLC



The AERIS SOLUTION foot system exceeds expectations by offering the best combination of value and performance in a 100% carbon composite foot. The AERIS Performance adds Compliant Composite Technology (CCT) that is derived by using multiple springs to increase the range of flexibility. The result is a unique Advanced Energy Release (AER) pylon design composed of nested parallel springs that enable the foot to react in a compliant manner at heel strike and then work together to create a dynamic transition from push-off to toe release. The AER pylon takes advantage of independent spring action to provide smooth, stable progression from heel-strike to toe-off. The AERIS SOLUTION foot system will exceed expectations for the low to moderate activity user.

## PRODUCT FEATURES

- 100% carbon composite energy storing foot
- Compliant Composite Technology (CCT)
- Advanced Energy Release (AER) pylon
- Enhanced plantar flexion for stability at heel strike
- Balanced compliance and dynamic response
- Smooth transitions from foot-flat to toe-off
- Proven durability enhances patient confidence
- Simple, maintenance-free design ensures long term worry-free use
- Low profile design fits majority of amputee population

## SPECIFICATIONS

- Patient weight: Up to 275 lbs. (125 kg)
- Foot weight: 16.0 oz. (450 g)
- Build height: 5.12 in. (13 cm)
- Functional level: K3
- Durable; meets ISO-22675 standard

## L-CODES\*

**L5976** All Lower Extremity Protheses, Energy Storing Foot

\*Suggested L-Codes are provided as a reference only, it is the responsibility of the practitioner to confirm this information.

## WARRANTY

- 24 Months from the date of shipment to the practitioner
- Foot Shell - 6 months from the date of shipment to the practitioner

## SATISFACTION GUARANTEE

- 30 days from the date of shipment to the practitioner

## ORDERING INFORMATION

**130-10-XX-X** AERIS SOLUTION Foot System

Kit includes foot, elastomer heel wedges, Spectra Sock, and instruction manual. Foot selection requires foot length and weight.

		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
99-130 lbs.	45-59 kg	22-A	23-A	24-A	25-B	26-B	27-B	28-B	29-D	30-D
131-163 lbs.	60-74 kg	22-A	23-A	24-A	25-B	26-B	27-C	28-C	29-D	30-D
164-198 lbs.	75-90 kg	22-B	23-B	24-B	25-C	26-C	27-C	28-C	29-D	30-D
199-235 lbs.	91-107 kg	22-B	23-B	24-B	25-C	26-C	27-D	28-D	29-E	30-E
236-275 lbs.	108-125 kg	NA	NA	NA	25-D	26-D	27-D	28-D	29-E	30-E

## ACCESSORIES

**XX XX 13 CC 3** Micro Coated Vinyl Foot Shell  
Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24-30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30

Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.

## Spectra Sock

**192500** Spectra Sock\* (sold separately)  
\*Recommended option when using foot shell.

