

DuraShock



Fillauer®

DuraShock






A shock and torsion reducing pylon simulates the rotation and shock absorption of the natural limb during stance phase. Beneficial for everyday use, the DuraShock is designed for common activities, including gardening, hiking, tennis, golf, and traversing stairs.

Product Features

- Up to 60° of axial rotation (30° in each direction)
- Shock absorption for increased socket comfort
- Maintenance-free design
- Titanium core technology

Warranty

- 1 year from date of fitting

	<p>DS003L DuraShock 30 mm, Long Pylon, Red, 144 lbs. (65 kg) DS002L DuraShock 30 mm, Long Pylon, Blue, 179 lbs. (81 kg) DS001L DuraShock 30 mm, Long Pylon, Black, 220 lbs. (99 kg) DS004L DuraShock 34 mm, Long Pylon, Heavy Duty, 250 lbs. (113 kg)</p>	<ul style="list-style-type: none"> • 30 mm pylon, 15 oz. (428 g) • 34 mm pylon, 21.7 oz. (609 g) • Length as shipped 11.25 in. (28.6 cm) • 30 mm pylon minimum length 5.75 in. (14.6 cm) • 34 mm pylon minimum length 6.75 in. (17.145 cm) • Dampening ring included
	<p>DS003S DuraShock 30 mm, Short Pylon, Red, 144 lbs. (65 kg) DS002S DuraShock 30 mm, Short Pylon, Blue, 179 lbs. (81 kg) DS001S DuraShock 30 mm, Short Pylon, Black, 220 lbs. (99 kg) DS004S DuraShock 34 mm, Short Pylon, Heavy Duty, 250 lbs. (113 kg)</p>	<ul style="list-style-type: none"> • 30 mm pylon, 11 oz. (318 g) • 34 mm pylon, 15.2 oz. (425 g) • Length as shipped 7.5 in. (19 cm) • 30 mm pylon minimum length 4.75 in. (11.4 cm) • 34 mm pylon minimum length 5 in. (12.7 cm) • Dampening ring included
	<p>DS003SC DuraShock 30 mm, Tube Clamp, Red, 144 lbs. (65 kg) DS002SC DuraShock 30 mm, Tube Clamp, Blue, 179 lbs. (81 kg) DS001SC DuraShock 30 mm, Tube Clamp, Black, 220 lbs. (99 kg)</p>	<ul style="list-style-type: none"> • Weight 11.6 oz. (330 g) • Length as shipped 6 in. (15.2 cm) • Minimum build height 6.75 in. (17 cm) • Dampening ring included
	<p>DS003SS DuraShock, Short, Red, 144 lbs. (65 kg) DS002SS DuraShock, Short, Blue, 179 lbs. (81 kg) DS001SS DuraShock, Short, Black, 220 lbs. (99 kg) DS004SS DuraShock, Short, Heavy Duty, 300 lbs. (136 kg)</p>	<ul style="list-style-type: none"> • Standard build height 3.5 in. (8.9 cm) • Heavy duty build height 4.25 in. (10.8 cm) • Dampening ring included
	<p>MS001SS MiniShock, 132 lbs. (60 kg)</p>	<ul style="list-style-type: none"> • Standard build height 2.38 in. (6 cm) • Dampening ring included • Pyramid size: Medium

Patents 6468315 & 6682569

Fillauer®

www.fillauer.com

Fillauer LLC

2710 Amnicola Hwy
 Chattanooga, TN 37406
 423.624.0946
 Fax 423.629.7936

AD371/08-14-15/06-06-17