

## General Instructions for MightyMite® Modular Shuttle Locks

- Avoid humid or wet environments and always dry the device thoroughly should it get wet.
- ISO Weight Rating – 60kg/132lbs (based upon the strength of the distal attachment)
- Single Use Device – Do not re-use on another patient or re-fabricate.
- The MightyMite Shuttle Locks were designed to mate with the existing MightyMite components. FC100204 mates with the MightyMite four hole pattern and FC100202 mates with MightyMite medium pyramid receivers.

## Instructions to the Practitioner

- The practitioner should discuss indications, contraindications, care and maintenance, and application instructions with the patient and/or caregiver before use. In selecting the appropriate device, the practitioner should carefully consider the abilities of the patient. The patient or caregiver must be able to determine if the plunger is locked into place during daily application.

*Fillauer*®

2710 Amnicola Highway • P.O. Box 5189 • Chattanooga, TN 37406 USA • 800 251 6398 • 423 624 0946  
European Representative: Basko Healthcare • Haspelsstraat 6 • 1063 TE Amsterdam • The Netherlands

PM156/06-02

