

**CAUTION - Be sure to use the proper size pyramids and receivers together.**

**Pyramid sizes are defined as:**

- **SMALL** Mates with 20mm spherical radius (commonly used throughout many pediatric systems including on top of the knee).
- **MEDIUM** Mates with 22mm spherical radius (standard pyramid for MightyMite system, also used at the ankle in some pediatric systems).
- **LARGE** Mates with 30mm spherical radius (commonly used throughout adult systems).

**MightyMite pyramid and receiver components are compatible with:**

FC100110	SACH Foot Adapter w/Med Pyramid	Medium Pyramid
FC100114	Four Hole Adapter w/Med Pyramid	Medium Pyramid
FC100118	Integrated Pylon w/Aluminum Receiver	Medium Pyramid
FC100119	Integrated Steel Pylon w/Titanium Receiver	Medium Pyramid
FC100120	Integrated Steel Pylon w/Aluminum Receiver	Medium Pyramid
FC100121	Integrated Steel Pylon w/Aluminum Receiver	Large Pyramid
FC100122	Aluminum Tube Clamp f/Medium Pyramid	Medium Pyramid
FC100123	Titanium Tube Clamp f/Medium Pyramid	Medium Pyramid
FC100124	Aluminum Tube Clamp f/Large Pyramid	Large Pyramid
FC100134	Rotatable Medium Pyramid Receiver	Medium Pyramid
FC100144	Rotatable Small Pyramid Receiver	Small Pyramid
FC100152	BK Socket Adapter	Medium Pyramid



Fillauer LLC • 2710 Amnicola Highway • Chattanooga, TN 37406 • Tel: 800.251.6398 • [www.fillauer.com](http://www.fillauer.com)  
Hosmer • 561 Division Street • Campbell, CA 95008 • Tel: 800.827.0070 • [www.hosmer.com](http://www.hosmer.com)  
Centri® • Kung Hans Väg 2 • 192 68 Sollentuna, Sweden • Tel: +46 8 505 332 00 • [www.centri.se](http://www.centri.se)



PM195/04-09