

# Hydraulic Ankle Feet

Hydraulic Ankle Foot	Sizes	Weight Limit	Weight	Activity Level	L-code*	Build Height
<b>Raize Microprocessor Foot</b> 	24 – 30 cm	220 lbs. 100 kg	27.7 oz. 785 g	K3	L5973	3.5 – 3.75 in. 9 – 9.5 cm
<b>Motion Foot MX</b> 	22 – 30 cm	220 lbs. 100 kg	28 oz. 794 g	K3 – K4	L5981 L5968 L5990 L5999**	3.25 in. 8.3 cm

# Composite Feet

Foot	Sizes	Weight Limit	Weight	Activity Level	L-code*	Build Height
<b>K2</b> 	22 – 30 cm	365 lbs. 165 kg	14.1 oz. 395 g	K2	L5972 L5986	3.5 in. 9 cm
<b>Wave Comfort 2</b> 	22 – 30 cm	330 lbs. 150 kg	13.2 oz. 375 g	K3	L5981 L5986	4.25 in. 10.8 cm
<b>Element</b> 	22 – 30 cm	300 lbs. 136 kg	12 oz. 340 g	K3 – K4	L5981 L5986	3.75 in. 9.5 cm
<b>Element DS</b> 	22 – 30 cm	275 lbs. 125 kg	22.3 oz. 632 g	K3 – K4	L5987 L5984	6.5 in. 16.5 cm
<b>Ibex</b> 	22 – 30 cm	330 lbs. 150 kg	16.2 oz. 459 g	K3 – K4	L5980 L5986	6.75 in. 17 cm
<b>Ibex XD</b> 	25 – 30 cm	440 lbs. 200 kg	21.2 oz. 602 g	K3 – K4	L5980 L5986	7.63 in. 19 cm
<b>Aeris Performance</b> 	22 – 30 cm	330 lbs. 150 kg	18.4 oz. 520 g	K3 – K4	L5980 L5986	7.25 in. 18.4 cm
<b>Wave Sport</b> 	22 – 30 cm	275 lbs. 125 kg	19.4 oz. 550 g	K3 – K4	L5987	8 – 10 in. 20 – 25 cm
<b>Wave Sport V3</b> 	25 – 30 cm	275 lbs. 125 kg	19.4 oz. 550 g	K3 – K4	L5987	4 in. 10 cm
<b>AllPro</b> 	19 – 30 cm	330 lbs. 150 kg	15.9 oz. 450 g	K3 – K4	L5987 L5986	6 in. (15.2 cm) 8 in. (20.3 cm) 10 in. (25.4 cm) 12 in. (30.5 cm)

\*Suggested L-Codes are provided as a reference only. It is the responsibility of the practitioner and/or patient care facility submitting this claim to determine and submit the correct coding within all established rules and guidelines.

\*\*Contact Fillauer if you require more information.