

# Fillauer Reciprocating Gait Orthosis (RGO)

Practitioner \_\_\_\_\_ Branch Name \_\_\_\_\_  
 PO# \_\_\_\_\_ Date: \_\_\_\_\_ Ship Via: \_\_\_\_\_  
 Address \_\_\_\_\_  
 City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

**Fillauer** LLC

Tel: 800.251.6398  
 Fax: 423.698.6076  
 2710 Amnicola Highway  
 Chattanooga, TN 37406



### Patient Information

Patient's Name \_\_\_\_\_  
 Age \_\_\_\_\_ Sex \_\_\_\_\_ Weight \_\_\_\_\_ Height \_\_\_\_\_  
 Ship To Address \_\_\_\_\_  
 \_\_\_\_\_  
 Fabricate Complete RGO       Fabricate Pelvic Section Only

### RGO Options

**Cable Type**  
 Horizontal Cable     Dual Cable  
 Rocker Bar (Metal Pelvic Section Only)

**Pelvic Section**  
 Plastic Band       Metal Band

**Uprights**  
 Single       Double

**Hip Joints**  
 Regular Thrust Bearing Hip Joints  
 Abduction Hip Joints  
 Non-Abduction Hip Joints

**Hip Joint Sizes**  
 Small (3/16x5/8)  
 Medium (1/4x13/16)  
 Large (5/16x7/8)


**Knee Joint**  
 Fillauer Cam Lock  
 Drop Locks       Bail Lock  
 Other \_\_\_\_\_



**Knee Joint Bar Size**  
 3/16 x 5/8  
 3/16 x 3/4       1/4 x 3/4

**Heel Height**  
 0"       1/4"       1/2"

**Additional Instructions**  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### Measurements

- Provide all measurements, especially the M-L diameters taken with caliper gauge
- Supply a profile tracing from axilla to feet. Please check that tracing matches caliper diameters.
- Plaster molds of each extremity are required.
  - Maintain ankle joint in neutral position for low heel shoe.
  - Correct ankle inversion or eversion as much as possible.
  - Reduce knee valgus or varus.
  - Use a right angle casting board.
- When a plastic pelvic band is desired, refer to special casting suggestions in the RGO Fabrication Manual.

M-L at Metatarsal Head