

# Fillauer Shower Leg

Practitioner \_\_\_\_\_ Branch Name \_\_\_\_\_  
 PO# \_\_\_\_\_ Date: \_\_\_\_\_ Ship Via: \_\_\_\_\_  
 Address \_\_\_\_\_  
 City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

*Fillauer* LLC

Tel: 800.251.6398  
 Fax: 423.698.6076  
 2710 Amnicola Highway  
 Chattanooga, TN 37406




This unique Copolymer Prosthesis presents a solution to slippery shower surfaces for the Transtibial Amputee. It is waterproof, lightweight, and easy to don and doff with a suspension sleeve.

The I.D. of the Pylon Section is consistent to just distal of the Socket Flair so that it is possible to shorten the shank a considerable amount and still retain internal locking of the pod. Therefore, the Shower Leg is fabricated longer than the measurement given so the prosthetist can determine the proper length.

As a custom made item, a negative of a modified model (preferred) or an original negative wrap, patient information, fabrication information, measurements, and choice of transfer paper is required.

## Measurements



Remarks: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

## Patient Information

Patient's Name \_\_\_\_\_  
 Patient's Weight \_\_\_\_\_ Patient's Height \_\_\_\_\_  
 Amputation Level: \_\_\_\_\_  Right  Left

## Fabrication Information














The following has been provided for Fabrication:

- Original Negative Wrap  A Negative of a Modified Model (Preferred)

Is Alignment Normal?  No, If No:  Yes, Normal Alignment  
 Flexion Angle \_\_\_\_\_ °  
 AD&AB Angle \_\_\_\_\_ °

## Transfer Paper

Transfer Paper is also available in 30" width, sold per yard

- |                                                                                                                                                 |                                                                                                                                         |                                                                                                                                            |                                                                                                                                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <br><input type="checkbox"/> 999000<br>Animals & Blocks      | <br><input type="checkbox"/> 999001<br>Butterflies   | <br><input type="checkbox"/> 999002<br>Jungle Animals   | <br><input type="checkbox"/> 999003<br>Camouflage   |
| <br><input type="checkbox"/> 999005<br>Multi-Color           | <br><input type="checkbox"/> 999006<br>Chickies      | <br><input type="checkbox"/> 999008<br>Optical Illusion | <br><input type="checkbox"/> 999009<br>Martians      |
| <br><input type="checkbox"/> 999011<br>Caucasian            | <br><input type="checkbox"/> 999020<br>Leopard Print | <br><input type="checkbox"/> 999021<br>Solid Blue       | <br><input type="checkbox"/> 999022<br>American Flag |
| <br><input type="checkbox"/> 999024<br>Extreme Skateboarding | Other: _____                                                                                                                            |                                                                                                                                            |                                                                                                                                         |