

THE XD LINE

Rated to 440 pounds



IBEX XD

- Innovative Micro-Slice technology for multi-axial performance
- Energy optimized pylon/heel plate design
- Carbon fiber foot/shank design for outstanding energy return
- Adjustable heel wedge for fine tuning heel stiffness
- Build height: 7.63 in. (19 cm)

ORDERING INFORMATION

120-10-25-E8XD	Ibex XD, Size 25
120-10-26-E8XD	Ibex XD, Size 26
120-10-27-E9XD	Ibex XD, Size 27
120-10-28-E9XD	Ibex XD, Size 28
120-10-29-E9XD	Ibex XD, Size 29
120-10-30-E9XD	Ibex XD, Size 30

L-CODES

- L5980** All lower extremity prostheses, flex-foot system
- L5986** All lower extremity prostheses, multi-axial rotation unit ("MCP" or equal)

Suggested L-Codes are provided as a reference only. It is the responsibility of the practitioner and/or patient care facility submitting this claim to determine and submit the correct coding within all established rules and guidelines.

MCV FOOT SHELL

XX XX 13 CC 3 Micro Coated Vinyl Foot Shell (Sold Separately)
Example: 45 24 13 13 3 = Left, Size 24, Color 13

To order, select the side (left or right) and foot length (24–30 cm) from the chart below. Then, choose the color (CC 03, 09, or 13).

	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 25	46 26	46 27	46 28	46 29	46 30



Custom colors are available as a special order with three weeks lead time. Please contact Fillauer Customer Service for more details.

ACCESSORIES



126100 XD 34 mm Tube Clamp, Titanium



126120 XD 4-Prong Lamination Adapter



126133 XD 34 mm Pylon

Fillauer
COMPANIES

Fillauer LLC • 800.251.6398 • 423.624.0946 • www.fillauer.com
 Hosmer • 800.827.0070 • 408.379.5151 • www.hosmer.com