


THIS INSERT REPLACES PAGE 2 OF THE SPL MANUAL

Features and Benefits

- Automatic Locking Prior To Stance Phase
- Closely Mimics Normal Gait Patterns
- Free Flexion During Swing Phase
- Minimizes Energy Consumption
- Eliminates Circumduction Vaulting
- Toe Clearance During Swing Phase
- Patient Security
 - Increased Flexion Angle of 151 degrees
 - Mechanical Security Stops at 15 and 32 degrees of Flexion*
 - Extended wear Bumper*
 - Teflon Bushings for Longer Wear*
 - PEEK High Wear Insert*
- 3 Modes of Control
- No Cables or Special Heel Connections
- Easy Fabrication

New Features of the SPL 2



Warranty

Fillauer orthotic components are warranted for 12 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor. Items under warranty will be replaced or repaired (at Fillauer's discretion) at no charge. The warranty will be void if the item has been fabricated or installed outside Fillauer's recommendations, if the item has been exposed to a corrosive environment, or if the item has been used in extremely abusive activities that could result in injury.

Indications:

- Post-Polio
- Spinal Involvement
- Cerebral Vascular Accident
- Peripheral Paresis/Paralysis
- Nerve Inflammations
- Neurological Failures
- Myopathies
- Multiple Sclerosis (or similar diseases)

Contraindications:

- Knee Flexion Contracture greater than 10°
- Central Paralysis
- Hip Flexion Contracture
- Hip Musculature Involvement
- Poor Balance or Coordination
- Knee Hyperextension greater than 10°