



obsidian

*Fillauer*®

# Ready. Set. GO!

Get ready for race day with the **Obsidian** blade line.


The **Obsidian Running Blades** are available in a posterior mount option for kids and adults. The adult direct mount option offers a pyramid slide for a variety of alignment locations and easy adjustments.

The new **Obsidian Sprint Blade**, intended for elite athletes, provides a curve hugging split toe along with a snap at toe-off that will push the elite to a new level.




## Split-Toe Sole

The Obsidian has an exclusive split-toe design for improved multi-axial function and cornering performance. The running blades offer a tire-grade tread that can grip surfaces from the track to the trail.




## STF Technology

The unique STF carbon enhances durability which allows athletes to push their limits with a blade that can keep up.



## Dynamic Carbon Composite

Obsidian offers outstanding energy return through the curved shape and dynamic, waterproof carbon composite material.



# obsidian



## Specifications

	Running Blade Direct Mount	Pediatric Running Blade Posterior Mount	Running Blade Posterior Mount	Sprint Blade Posterior Mount
Patient Weight	Up to 330 lbs. (150 kg)	Up to 180 lbs. (82 kg)	Up to 330 lbs. (150 kg)	Up to 275 lbs. (125 kg)
Foot Weight	19.4 oz. (550 g)	19.4 oz. (550 g)	28 oz. (800 g)	35 oz. (1000 g)
Clearance	10.75 in. (27.3 cm)	2.0 in. (5.1 cm)	2.0 in. (5.1 cm)	3 in. (7.6 cm)
Mounting Height	10.75 in. (27.3 cm)	10.5 in. – 18.5 in. (26.7 – 47.0 cm)	15.5 in. – 21.5 in. (39.4 – 54.6 cm)	15.5 in. – 26 in. (39.4 – 66.0 cm)
Width at Minimum Build Height		1.75 in. (44 mm)	2.25 in. (57 mm)	2.25 in. (57 mm)
Warranty	12 months	12 months	12 months	6 months
Patient Trial	Yes	No	No	No

## Ordering Information

Patient Weight	Running Blade Direct Mount	Pediatric Running Blade Posterior Mount	Running Blade Posterior Mount	Sprint Blade Posterior Mount*
50 – 74 lbs. (23 – 33 kg)		<b>175-120-3A</b>		
75 – 99 lbs. (34 – 44 kg)	<b>175-130-3A</b>	<b>175-120-2A</b>		
100 – 124 lbs. (45 – 56 kg)	<b>175-130-2A</b>	<b>175-120-A</b>		<b>177-1X0-X**</b>
125 – 149 lbs. (57 – 68 kg)	<b>175-130-A</b>	<b>175-120-B</b>	<b>175-110-A</b>	<b>177-1X0-X**</b>
150 – 179 lbs. (69 – 81 kg)	<b>175-130-B</b>	<b>175-120-C</b>	<b>175-110-B</b>	<b>177-1X0-X**</b>
180 – 209 lbs. (82 – 95 kg)	<b>175-130-C</b>		<b>175-110-C</b>	<b>177-1X0-X**</b>
210 – 239 lbs. (96 – 109 kg)	<b>175-130-D</b>		<b>175-110-D</b>	<b>177-1X0-X**</b>
240 – 275 lbs. (110 – 125 kg)	<b>175-130-E</b>		<b>175-110-E</b>	<b>177-1X0-X**</b>
276 – 330 lbs. (125 – 150 kg)	<b>175-130-F</b>		<b>175-110-F</b>	

\*Obsidian Sprint does not include track sole

\*\*Obsidian Sprint requires order form available at [www.fillauer.com](http://www.fillauer.com)

## Accessories

Description	Running Blade Direct Mount	Pediatric Running Blade Posterior Mount	Running Blade Posterior Mount
Adjustable Alignment Bracket		<b>180-10-5000</b>	<b>180-10-4000</b>
Posterior Mounting Bracket		<b>180-10-3000</b>	<b>180-10-2010</b>
Replacement Running Tread	<b>100-40-4000</b>	<b>100-40-4000</b>	<b>100-40-4000</b>

*Fillauer*®



[www.fillauer.com](http://www.fillauer.com)

**Fillauer LLC**

2710 Amnicola Highway  
Chattanooga, TN 37406  
423.624.0946  
Fax 423.629.7936

© 2019 Fillauer LLC  
AD374/07-15-19/10-01-20