

MC Waterproof Collar for Quick Disconnect Wrist



Fillauer[®]
Motion Control

MC Waterproof Collar for Quick Disconnect Wrist



Product Description

The Motion Control Waterproof Collar for Quick Disconnect Wrist creates a waterproof seal at the wrist of a myoelectric prosthesis. When used in conjunction with the Motion Control End Cap this allows the user of a Motion Control ETD or ETD2 to now submerge their transradial prosthesis beyond the wrist.

Wearers can confidently submerge their hand into a sink of sudsy dishes or grab that trophy fish trying to get away without the fear of water ingress through the wrist unit. Yet, they can still interchange their waterproof device with a non-waterproof hand for an evening out.*

*The Waterproof Collar for Quick Disconnect Wrist is compatible with the ETD, ETD2 and the Motion Control Waterproof Endcap. Contact Motion Control for details.

Product Features

- Waterproof to IPX7 (tested at 1 meter for 30 minutes)
- Any manufacturer's non-waterproof devices may be interchanged with the Quick Disconnect wrist
- Convenient Maintenance Schedule
- Available with ETD and ETD2

Warranty

1 year limited warranty

Satisfaction Guarantee

30 days from date of shipment

Ordering Information

www.UtahArm.com/Catalog

Fillauer
Motion Control®



www.fillauer.com

Motion Control, Inc.

115 N. Wright Brothers Drive
Salt Lake City, UT 84116
Phone 801.326.3434
Fax 801.978.0848
Toll Free 888.696.2767

© 2020 Motion Control, Inc.
WQD Sell Sheet 07-29-2020