



Absolute Suspension Sleeve

Fillauer®

Absolute Suspension Sleeve

The Absolute Suspension Sleeve optimizes comfort, durability, and range of motion. The sleeve is made of thermoplastic elastomer and covered with a reinforced Lycra® type material. The pre-flexed design allows maximum range of motion and alleviates posterior bunching when the knee is flexed. The Absolute Suspension Sleeve is available in multiple sizes and your choice of beige or black.

Product Features

- Maximum comfort
- Slides against clothing
- Available in black or beige fabric
- Made of thermoplastic elastomer and covered with a reinforced Lycra® type material.
- The inside of the pre-flexed sleeve is reinforced with.

Warranty

90 day warranty

Absolute Pre-flexed Suspension Sleeve, 3 mm - Extra Soft

Size	Distal Circumference	Proximal Circumference	Beige	Black
Small	21 - 25 cm (8.3 - 9.8 in.)	26 - 31 cm (10.2 - 12.2 in.)	803471314	803461314
Medium	26 - 32 cm (10.2 - 12.6 in.)	31 - 38 cm (12.2 - 15.0 in.)	803471318	803461318
Medium+	30 - 36 cm (11.8 - 14.2 in.)	38 - 44 cm (15.0 - 17.3 in.)	803471322	803461322
Large	32 - 40 cm (12.6 - 15.8 in.)	44 - 50 cm (17.3 - 19.7 in.)	803471326	803461326
Large+	34 - 42 cm (13.4 - 16.5 in.)	50 - 60 cm (19.7 - 23.6 in.)	803471332	803461332
X-Large	40 - 49 cm (15.7 - 19.3 in.)	60 - 70 cm (23.6 - 27.5 in.)	803471340	803461340



Absolute Straight Suspension Sleeve 3 mm - Extra Soft

Size	Distal Circumference	Proximal Circumference	Beige	Black
X-Small	10 - 17 cm (3.9 - 6.7 in.)	16 - 24 cm (6.3 - 9.4 in.)	803371301	803361301
Small	15 - 22 cm (5.9 - 8.7 in.)	23 - 34 cm (9.1 - 13.4 in.)	803371302	803361302
Medium	21 - 30 cm (8.3 - 11.8 in.)	29 - 43 cm (11.4 - 16.9 in.)	803371303	803361303
Medium+	24 - 35 cm (9.4 - 13.8 in.)	33 - 47 cm (13 - 18.5 in.)	803371304	803361304
Large	28 - 42 cm (11 - 16.5 in.)	36 - 51 cm (14.2 - 20.1 in.)	803371305	803361305
Large+	31 - 46 cm (12.2 - 18.1 in.)	39 - 55 cm (15.4 - 21.7 in.)	803371306	803361306
X-Large	34 - 50 cm (13.4 - 19.7 in.)	42 - 60 cm (16.5 - 23.6 in.)	803371307	803361307
XX-Large	36 - 55 cm (14.2 - 21.7 in.)	44 - 70 cm (17.3 - 27.6 in.)	803371308	803361308



Fillauer



www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936

© 2018 Fillauer LLC
AD344/06-06-16/01-02-18