



ProCover

Patent pending.

Fillauer®

ProCover

The Fillauer ProCover is made of durable and resilient polyurethane foam that enhances performance and protects the foot module. The split design makes the cover easy to install and remove for ease of maintenance and cleaning. The ProCover is intended for use with the AllPro, AllPro PM, or Wave Comfort 3.

Product Features

- Durable design and material is ideal for high activity patients
- Split design allows for easy maintenance, cleaning, or switching with other cosmetic covers
- Design-to-fit shape securely fits the foot plate for a more responsive feel
- Polyurethane foam maintains superior energy return over time
- Unique shape enhances foot rollover performance

Specifications

Compatible feet: AllPro, AllPro PM, and Wave Comfort 3

Sizes: 22 – 30 cm

Weight (26 cm): 9.1 oz. (259 g)

Material: Polyurethane foam

Activity level: K3 – K4

Warranty

6 month warranty

Ordering Information

ProCover purchase includes two hardware kits and one soling sheet.

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
ProCover	182-10-1022	182-10-1023	182-10-1024	182-10-1025	182-10-1026	182-10-1027	182-10-1028	182-10-1029	182-10-1030

Replacement Parts

	Quantity	Size	Part Number
Hardware Replacement Kit	4 Connector Pins, 2 Velcro® Straps	22 – 25 cm	182-10-2300
Hardware Replacement Kit	4 Connector Pins, 2 Velcro® Straps	26 – 30 cm	182-10-2350

Patent pending.



Fillauer®

www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936

© 2018 Fillauer LLC
AD389/12-11-18