

Custom Pediatric Formula Order Form



2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
info@fillauer.com

Company

Account No.

Address Line 1

Address Line 2

City, State, Zip

Practitioner

Phone

Email

Purchase Order

Patient Information

Amputation Level

Patient ID

Patient Age

Patient Height

Patient Weight

Ordering Information

Custom Pylon Part Number

Foot Plate Part Number


Pediatric Posterior Mounting Bracket
(Included: 180-10-3000)

Custom Pylon Part Number

Patient Weight		Pylon
30 – 45 lbs.	14 – 20 kg	178-120-5A-C
46 – 60 lbs.	21 – 27 kg	178-120-4A-C
61 – 75 lbs.	28 – 34 kg	178-120-3A-C
76 – 99 lbs.	35 – 45 kg	178-120-2A-C
100 – 125 lbs.	45 – 57 kg	178-120-A-C

Foot Plate Part Number

Foot Size	Left	Right
16 cm	178-120-16L	178-120-16R
17 cm	178-120-17L	178-120-17R
18 cm	178-120-18L	178-120-18R
19 cm	178-120-19L	178-120-19R
20 cm	178-120-20L	178-120-20R
21 cm	178-120-21L	178-120-21R
22 cm	178-120-22L	178-120-22R
23 cm	178-120-23L	178-120-23R
24 cm	178-120-24L	178-120-24R



Custom Pylon Length

Measure the maximum height to avoid popliteal impingement in full flexion. (Normally this will be the posterior trim line.)

Must be 9 – 16 in. (23 – 41 cm)

Additional Instructions
