

ProControl 2



Fillauer[®]
Motion Control

ProControl 2 from Motion Control



Product Description

The ProControl 2 provides the ability to myoelectrically control both the terminal device(s) and the wrist rotator from one modular controller. Requires Motion Control Standard Hand or ETD.

Using microprocessor technology, you can tailor the inputs and the outputs to the patient's needs. Switching between the hand and the wrist is easy – choose from three options. A variety of available harnesses and/or wire adapters allows for almost any battery or input combination.

Satisfaction Guarantee

30 days from date of shipment

Warranty

2 year limited warranty

Product Features

- Control Terminal Device and Wrist with one controller
- Use a variety of inputs
- Dual or Single Channel EMG
- AutoCal - allows the wearer to the adjust sensitivity
- Ideal for controlling (optional) Hosmer E-TWO Electric Elbow
- 2 year limited warranty

Optional Features

Input/Switch Options

- EMG
- Linear Potentiometer
- Touch Pad
- Hand/Wrist switching, or direct switch control of wrist

Terminal Devices

- MC Motor Direct Electric Terminal Device (ETD)
- MC Motor Direct Hand

Battery Options

- Internal single cell or dual (split) cell
- Other companies' batteries

Specifications

Operating Temperature: 0° to 44° C (32° to 110° F)

Transport & Storage Temperature: -18° to 60° C (0° to 140° F)

Operating Voltage: 7.2 Vdc

More specifications at www.UtahArm.com

Ordering Information

www.UtahArm.com/Catalog

Fillauer
Motion Control

www.fillauer.com

Motion Control, Inc.
115 N Wright Brothers Drive
Salt Lake City, UT 84116, USA
Phone 801.326.3434
Fax 801.978.0848
Toll Free 888.696.2767

ProControl 2 Brochure 07-28-2020
© 2020 Motion Control, Inc.