



Dual Purpose Expansion Arbors

	<p>Small Arbor</p> <p>225330 1/2-13 Threads 225335 5/8-11 Threads 225330MM M16 Threads 225345 Slim Blue Coneiak Replacement</p>	<p>Small Sanding Sleeves</p> <p>D24 Small 24 Grit D80 Small 80 Grit D120 Small 120 Grit</p>
	<p>Medium Arbor</p> <p>225320 1/2-13 Threads 225325 5/8-11 Threads 225320MM M16 Threads 225340 Tapered Red Coneiak Replacement</p>	<p>Medium Sanding Sleeves</p> <p>C24 Medium 24 Grit C80 Medium 80 Grit C120 Medium 120 Grit C180 Medium 180 Grit C240 Medium 240 Grit C320 Medium 320 Grit</p>

Patent pending

Maintenance

- Arbor should be cleaned/wiped off with alcohol weekly to keep the surface of the arbor tacky
- Visually inspect arbor weekly to make sure it has not been damaged

Warranty

This product has a 6 month warranty against manufacturing defects.

At no time should the expansion arbor be used without the proper sanding sleeve listed above. Doing so could cause the arbor to fail and will void the warranty.

Polishing arbors wear with use. Wear is dependent upon the material being polished, the pressure applied by the user, and the length of use. Products returned with extensive wear will not be accepted for warranty credit. Polishing arbors are recommended for thermoplastic materials only. Warranty will be voided if used on thermoset materials (i.e. acrylic resins, laminations, etc.).

Please Note: Air bubbles may become visible after use but will not affect product performance.

Fillauer®

www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
 Chattanooga, TN 37406
 423.624.0946
 Fax 423.629.7936

© 2020 Fillauer LLC
 PM300/04-21-15/12-03-20