



**Even more
flexible.**

AllPro DS Foot

**Foot
Weight**
25.6 oz.
720 g

**Maximum
Patient Weight**
250 lbs.
113 kg

**Build
Height**
8.75 / 10.75 in.
22.2 / 27.3 cm

**Activity
Level**
K3 / K4

Fillauer®



60° of Axial Rotation
The integrated DuraShock pylon provides rotation and shock absorption.

Shorter & Lighter
Integrating the DuraShock is 10% lighter and 0.75 in. (1.9 cm) shorter than adding a DuraShock separately.



More Alignment Options
The RPI (10.75 in only) allows two different placements of the alignment in reference to the foot length. Either 30% or 35% of foot length.





Perfect Heel Stiffness

The replaceable heel bumper let's you fine-tune heel stiffness for an even better experience.

Industry-Leading Flexibility

The split-toe, carbon composite design provides multiaxial movement and dynamic response.


AllPro DS

Specifications

Weight limit: 250 lbs. (113 kg)
 Heel height: $\frac{3}{8}$ in. (10 mm)
 Product weight: 25.6 oz. (720 g) Size 27 / 10.75 in.
 Warranty: 24 months for AllPro foot unit
 12 months for DuraShock pylon

Build Heights

8.75 in. (22.2 cm)
 10.75 in. (27.3 cm)



Ordering Information

Foot includes Spectra Sock. Foot Shell sold separately. DuraShock unit will be red, blue, or black based on weight category.

139-106-XX-XXX AllPro DS foot, 8.75 in. pylon, foot sizes 22 – 30 cm

139-108-XX-XXX AllPro DS foot, 10.75 in. pylon, foot sizes 22 – 30 cm

Patient Weight		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100 – 119 lbs.	45 – 54 kg	22-A2R	23-A2R	24-A2R	25-B4R	26-B4R				
120 – 139 lbs.	54 – 63 kg	22-A3R	23-A3R	24-A3R	25-B4R	26-B4R	27-C6B	28-C6B		
140 – 159 lbs.	64 – 72 kg	22-B4R	23-B4R	24-B4R	25-B4B	26-B4B	27-C6B	28-C6B	29-D7K	30-D7K
160 – 179 lbs.	73 – 81 kg	22-B4R	23-B4R	24-B4R	25-C5B	26-C5B	27-C6K	28-C6K	29-D7K	30-D7K
180 – 199 lbs.	82 – 90 kg	22-C5B	23-C5B	24-C5B	25-C5K	26-C5K	27-D7K	28-D7K	29-D7K	30-D7K
200 – 224 lbs.	91 – 102 kg	22-C5B	23-C5B	24-C5B	25-D6K	26-D6K	27-D7K	28-D7K	29-E8K	30-E8K
225 – 250 lbs.	102 – 113 kg				25-D6K	26-D6K	27-E8K	28-E8K	29-E8K	30-E8K

Foot Shell

Ordering Information

XX XX 13 CC 3 Foot Shell



CC 03

CC 09

CC 13

CC 99

To order, select the side (left or right) and foot length (22 – 30 cm) from the chart below. Then, choose the color (CC 03, 09, 13, or 99).

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30

Fillauer[®]

www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936

© 2021 Fillauer LLC
AD391/02-23-21