



Totally
waterproof.
Truly
revolutionary!

Introducing the TASKA HandGen2

Pushing the boundaries of waterproof even further, the TASKA HandGen2 is available now.

For more information, contact Fillauer or visit taskaprosthetics.com

Fillauer Motion Control
Exclusive distributor
for U.S. & Canada
801.326.3434

do more

 **TASKA™**

Fillauer®

Key features

Resettable

Knuckle Release

The releasable knuckles protect the fingers from damage and are easily reset.

Flexible fingers
Lateral spread allows the hand to effectively grip a wide variety of objects firmly.

Shock absorbing Knuckle Block
Gives the hand durability and allows the fingers to splay, offering a more natural grip.

Motorized thumb rotation

Motorized rotation mimics real hand movement and improves grip speed.

Grip Cycle Buttons
Easy to use buttons provide simple access to grip changes and other features.

Flexible wrist

Built-in flexible wrist with three lockable positions or free flexion mode.

New to the TASKA HandGen2

Release Button

Prevents accidental disconnection and provides unlimited wrist rotation.

Waterproof Quick Disconnect Wrist

Compatible with industry-standard QD socket connections, including wrist rotators.



Fillauer Motion Control
Exclusive distributor for U.S. & Canada
801.326.3434

do more

 **TASKA**[™]

Fillauer[®]