

Orthotics





In 1914, George Fillauer Sr. founded the Red Star Pharmacy in Chattanooga, Tennessee. In 1934, his son and fellow CPO, Carlton Fillauer expanded his father's business to establish Fillauer Orthopedic as an innovator in the orthotics and prosthetics industry. Today, the Fillauer Companies are headed by Karl Fillauer as Chairman and his son Michael as the CEO, since their inception, the Fillauer Companies have always been headed by CPO's, and have grown to be a globally renowned leader in the development, manufacturing, distribution and application of orthotic and prosthetic products. The expansion has mainly taken place by acquisitions of smaller manufacturing companies with a unique profile. Over the years we have become nine companies under the Fillauer umbrella and during the 100 year celebration in 2014 a decision was taken on how the company would be steered into the next century. Therefore a consolidation took place in 2015 where all the companies became Fillauer.

Fillauer Europe was founded as Centri in 1948. Located just outside Stockholm, Sweden, Centri joined the Fillauer Companies 2002. While Centri maintained the close, cooperative relationships it traditionally had with prosthetists in Scandinavia, as Fillauer Europe it has now become the European headquarters of Fillauer. With its extensive portfolio of products, Fillauer Europe is committed to working with prosthetists both in Europe, but also in other regions of the world in order to create only the best solutions for patients with all manner of prosthetic and orthotic needs. In addition, Fillauer Europe plays a lead Research & Development role within Fillauer. Fillauer Europe's products include upper and lower extremity prosthetics and the Dynamic Walk which has catapulted us to the forefront of the drop foot orthotics market. Our R&D department also works in partnership with many of Sweden's leading hospitals to develop innovative and high-tech machinery.



Index

Fillauer Product Warranty	6
Returns, Adjustments, and Credits	7
Catalog symbol legend	9

Dynamic Walk

Dynamic Walk - Size Charts	11
Dynamic Walk Standard, Double Double	11
Dynamic Walk - Standard	13
Dynamic Walk - Wide Ankle	14
Dynamic Walk - Double-Double	15
Dynamic Walk - Wide Ankle Double-Double	16
Dynamic Walk - Single Side Medial	17
Dynamic Walk - Single Side Lateral	18
Dynamic Walk - Single Side Anterior	19
Dynamic Walk - Shoe Mount	20
Dynamic Walk - Modular	21
Dynamic Walk - Mix & Match	22
Dynamic Walk - Mix & Match	23
Custom Calf Height or Foot Length	23
Order form - Dynamic Walk Custom Calf Height or Foot Length	24
Dynamic Walk - Custom Made	25
Order form - Dynamic Walk Custom Made Adult and Pediatric	26

Orthotic Joints

Modular Integrated Ankle Joint, Double Action	29
Modular Integrated Ankle Joint, Adjustable Range of Motion	29
Integrated Ankle Joint, Double Action	30
Integrated Ankle Joint, Adjustable range of motion	30
Integrated Ankle Joint, Dorsiflexion Assist	31
Split Tail Stirrups	32
Solid Stirrups - Wide Tongue	32
Short Leg Upright - Double Action	33
Short Leg Upright - limited/ free motion	33
Short Leg Upright - Toe Pick Up	34
Ankle Joint - Limited/ Free Motion	35
Ankle Joint - Toe Pick Up	35
Ankle Joint - Double Action	36
Stirrups - Standard	37
Stirrups - Split Tongue	38
Stirrups	39
Long Tongue & Wide Tongue	39
Rectangle Caliper Plates	40
Insert Stirrup For Plastic Orthotics	40
Polycar-C™ Ankle Inserts	41
Dorsiflexion Assist Ankle Joint Kits With Polycar-C	42

PDA Ankle Joint Kits With Polycar-C	43
Steplock®	44
Interlock™	45
2400 Joints	46
Bar Kit	47
2400 Bar Kit	47
Close Contour Pediatric Knee Joint	48
Release Kit SBK	48
Release Kit RBK	48
Release Kit - UBK	49
Lever Release Kit - LRK	49
Close Contour	50
Trigger Release for Camlocks	51
Zippy Quick Release	51
Drop Lock Aluminum Hip Joint Assemblies	52
Drop Lock Aluminum Knee Joint Assemblies	52
Drop Lock Low Profile Aluminum Knee Joint Assemblies	53
Spring Loaded Lever Lift	53

RGO & IRGO

RGO/IRGO Overview	55
Small Hooped Cable RGO	56
Large Hooped Cable RGO without Abduction	57
Large Hooped Cable RGO with Abduction	58
Hooped Cable RGO Cables and Housings	59
Latch Knob for Small Hooped Cable RGO	60
Push Button for Hooped Cable RGO	60
RGO II System for Hooped Cable RGO with Abduction	61
Small Horizontal Cable RGO	62
Medium Horizontal Cable RGO without Abduction	63
Medium Horizontal Cable RGO with Abduction	64
Large Horizontal Cable RGO without Abduction	65
Large Horizontal Cable RGO with Abduction	66
Horizontal RGO Cables	67
Latch Knob for Small Horizontal Cable and Rocker Bar RGOs	68
Push Button for Horizontal Cable and Rocker Bar RGOs	69
RGO II System for Horizontal Cable and Rocker Bar RGOs	69
Specialty Connectors for ISOCENTRIC® RGO	70
Cable Track Kit	71
Aluminum Pelvic Bands	71
Lower Abduction Bar Assembly	72
Lower Bar Assembly without Abduction	72
Knee Joints	73
Pelvic Bands	73
External Ankle Foot Orthosis (AFOs)	74
Standard In-Shoe Ankle Foot Orthosis (AFO)	74
Gauges	74



Spinal Braces

Providence Nighttime Brace	76
Scoliosis Full Time Brace	77
LSO - Single Opening	77
LSO - Bi-valve	78
TLSO - Single Opening (A.O.O.)	78
TLSO - Bi-valve	79
Flex Foam™ I	79
Flex Foam II	80
Flex Foam III	80
Pectus Brace	81
Order Form - Scoliosis/ Kyphosis Brace	83
Order Form - LSO/ TLSO, Single or Bivalve Opening	84
Order Form - Flex Foam	85
Order Form - Pectus Brace	86
Scoliosis Module	87
A.O.O. Spinal Brace	88
Basic Strap	89
Anterior Stay	89
Welding Rod	89
Spring Steel Strap	90
Chafes with Triangular Body Clasp	90

Hip Joints - Scottish Rite	106
Pelvic Bands and Thigh Bars Scottish Rite	107
Telescoping Spreader Bars Scottish Rite	107

Insoles

TFC - Thermoformable Composite Insole	92
---------------------------------------	----

PEVA

PEVA Shoe - Open Toe Boots - Firm Heel Cap	94
Clip-On Bar	96
Denis Browne Splint	96
Dobbs Bar	97
Dobbs Bar Spring Assisted Quick Release	97

Abduction Orthosis

COD SU*PER 3D Hip Joint	99
A.A.O. - Anti-Adduction Orthosis	99
Scottish Rite Kit	100
Thrust Bearing Hip Joint Assemblies	101
Thrust Bearing Hip Joint Components Scottish Rite	102
Motion Control Plate and Drop Lock Scottish Rite	103
Hip Joint Assemblies, Stainless Steel Scottish Rite	103
Component Kits with Telescoping Spreader Bar - Scottish Rite	104
Fixed 25° and 35° Abduction Thrust Bearing Hip Joints - Scottish Rite	105
Variable Abduction Thrust Bearing	



Fillauer Product Warranty

Fillauer items under warranty will be replaced or repaired (at Fillauer's discretion) at no charge, if the material or workmanship is found to be defective. Normal wear is not warranted. The warranty will be void if the item has been: a. fabricated or installed outside Fillauer's recommendations b. modified c. exposed to a corrosive environment d. if used outside of the products intended use. Below is a summary of warranties.

ORTHOTIC, PROSTHETIC COMPONENTS AND EQUIPMENT

are warranted for 12 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor with the following exceptions:

SUSPENSION SLEEVES, SILICONE GLOVES, SPINAL BRACES & CUSTOM LINERS/SLEEVES

are warranted for 3 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor.

Custom Fabrication Items, Prosthetic Liners, Dream Skin, Foot Shells, Cosmetic Gloves, Passive Hands and Trautman Tools & Accessories are warranted for 6 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor.

HOSMER HOOKS

are warranted for 18 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor.

ELBOWS

are warranted for 18 – 24 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor. See specific product warranty for more information.

FILLAUER FEET

are warranted for 12, 24, or 36 months

12 months: Running Blade

24 months: AllPro, K2, Motion Foot® MX, Formula and Raize

36 months: Aeris Activity, Aeris Performance 2, Aeris Performance LP, Element, Element DS, Ibex XD, and Wave Comfort 3

COAPT & TASKA

Are warranted for 24 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor.

All warranties are voided by alteration or misuse of the product.

Returns, Adjustments, and Credits

Inspect all goods immediately upon receipt for damage or errors in shipping. Call Customer Service within 3 business days in the event of errors or damage. Credit will not be issued for custom items made to specification or non-stock items older than 30 days, unless Fillauer shipped in error or there is defect in materials or workmanship.

Prior approval and a return authorization number must be obtained from Customer Service for all goods returned. Credit for returned merchandise is subject to inspection and may be subject to a 15% restocking fee for merchandise in resalable condition. Goods returned without prior approval may be subject to a 25% restocking fee. No credit shall be issued for goods that are more than 12 months old. Returns must be returned in unopened sellable packaging accompanied by a copy of original invoice, order number, or PO number. The return authorization number, provided by Customer Service, must appear on the shipping label.



Dynamic Walk

Orthotic Joints

RGO & IRGO

Spinal Braces

Insoles

Peva

Abduction Orthosis

Tiny Titans

Catalog symbol legend

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

Color symbols

 CC03/Tan

 CC99/Black

 Grey

 Green

 Pink

 Orange



Dynamic Walk

Specification symbols

SE

Manufactured in Sweden

US

Manufactured in USA

ES

Manufactured in Spain

Dynamic Walk - Size Charts

Dynamic Walk Standard, Double Double

SIZE	LENGTH	WIDTH	HEIGHT	CALF CIRCUMFERENCE
				
Small	9.1 in. (230 mm)	2.6 in. (67 mm)*	12.6 in. (320 mm)	9.4 - 13.0 in. (240 - 330 mm)
Medium	9.8 in. (250 mm)	2.8 in. (72 mm)*	14.0 in. (355 mm)	13.4 - 16.5 in. (340 - 420 mm)
Large	10.6 in. (270 mm)	3.1 in. (78 mm)*	14.6 in. (370 mm)	16.9 - 20.5 in. (430 - 520 mm)

*Width measured approximately one centimeter above the foot plate

Dynamic Walk- Wide Ankle, Wide Ankle Double - Double

SIZE	LENGTH	WIDTH*	HEIGHT	CALF CIRCUMFERENCE
				
Small	9.1 in. (230 mm)	2.8 in. (72 mm)*	13.2 in. (320mm)	9.4 - 13.0 in. (240 - 330 mm)
Medium	9.8 in. (250 mm)	3.1 in. (80 mm)*	14.0 in. (355 mm)	13.4 - 16.5 in. (340 - 420 mm)
Large	10.6 in. (270 mm)	3.4 in. (87 mm)*	15.0 in. (370 mm)	16.9 - 20.5 in. (430 - 520 mm)

*Width measured approximately one centimeter above the foot plate



Single Side Medial/Lateral

SIZE	LENGTH	HEIGHT	CALF CIRCUMFERENCE
			
Small	9.1 in. (230 mm)	13.4 in. (340 mm)	9.4 - 13.0 in. (240 - 330 mm)
Medium	9.8 in. (250 mm)	14.6 in. (370 mm)	13.4 - 16.5 in. (340 - 420 mm)
Large	10.6 in. (270 mm)	15.2 in. (385 mm)	16.9 - 20.5 in. (430 - 520 mm)

Single Side Anterior

SIZE	LENGTH	HEIGHT	CALF CIRCUMFERENCE
			
Small	9.2 in. (235 mm)	14.1 in. (360 mm)	9.4 - 13.0 in. (240 - 330 mm)
Medium	10.1 in. (257 mm)	15.1 in. (383 mm)	13.4 - 16.5 in. (340 - 420 mm)
Large	10.8 in. (276 mm)	16 in. (407 mm)	16.9 - 20.5 in. (430 - 520 mm)

Shoe Mount

SIZE	LENGTH	HEIGHT	CALF CIRCUMFERENCE
			
Small	6.4 in. (163mm)	14.2 in. (360 mm)	9.4 - 13.0 in. (240 - 330 mm)
Medium	6.4 in. (163 mm)	14.9 in. (380 mm)	13.4 - 16.5 in. (340 - 420 mm)
Large	6.4 in. (163 mm)	15.7 in. (400 mm)	16.9 - 20.5 in. (430 - 520 mm)

Modular

SIZE	LENGTH	HEIGHT
		
Small	9.1 in. (230 mm)	12.6 in. (320 mm)
Medium	9.8 in. (250 mm)	13.4 in. (340 mm)
Large	10.6 in. (270 mm)	14.6 in. (370 mm)

Dynamic Walk - Standard

SE

Product Features

- Lightweight
- Medial/ Laterally stable
- Flexible
- High strength
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700111121K	700111122K
Medium	700111141K	700111142K
Large	700111161K	700111162K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700111121	700111122
Medium	700111141	700111142
Large	700111161	700111162

PADDING

DESCRIPTION	PRODUCT NUMBER
Small	700100220
Medium	700100240
Large	700100260





Dynamic Walk - Wide Ankle

SE

Product Features

- Increased width between the PEEK-rods
- Lightweight
- Flexible
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700111321K	700111322K
Medium	700111341K	700111342K
Large	700111361K	700111362K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700111321	700111322
Medium	700111341	700111342
Large	700111361	700111362

PADDING

DESCRIPTION	PRODUCT NUMBER
Small	700100220
Medium	700100240
Large	700100260



Dynamic Walk - Double-Double

SE

Product Features

- Four PEEK-rods for stronger lift
- Allows for everyday activities
- Flexible
- Lightweight and high strength
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700117121K	700117122K
Medium	700117141K	700117142K
Large	700117161K	700117162K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700117121	700117122
Medium	700117141	700117142
Large	700117161	700117162

PADDING

DESCRIPTION	PRODUCT NUMBER
Small	700100220
Medium	700100240
Large	700100260





Dynamic Walk - Wide Ankle Double-Double

SE



Product Features

- Four PEEK-rods for stronger lift
- Allows for everyday activities
- Flexible
- Lightweight and high strength
- Warranty: 12 months
- Size descriptions on page 11-12

Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700117321K	700117322K
Medium	700117341K	700117342K
Large	700117361K	700117362K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700117321	700117322
Medium	700117341	700117342
Large	700117361	700117362

PADDING

DESCRIPTION	PRODUCT NUMBER
Small	700100220
Medium	700100240
Large	700100260



Dynamic Walk - Single Side Medial

SE

Product Features

- Lightweight
- Easy donning and doffing
- Available with PEEK-rods on either medial or lateral side
- Fits easily into normal shoes
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700113121K	700113122K
Medium	700113141K	700113142K
Large	700113161K	700113162K

Dynamic Walk **without** Padding

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700113121	700113122
Medium	700113141	700113142
Large	700113161	700113162

PADDING MEDIAL

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700100322	700100321
Medium	700100342	700100341
Large	700100362	700100361





Dynamic Walk - Single Side Lateral

SE

Product Features

- Lightweight
- Easy donning and doffing
- PEEK-rods on lateral side
- Fits easily into normal shoes
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700114121K	700114122K
Medium	700114141K	700114142K
Large	700114161K	700114162K

Dynamic Walk **without** Padding

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700114121	700114122
Medium	700114141	700114142
Large	700114161	700114162

PADDING LATERAL

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700100321	700100322
Medium	700100341	700100342
Large	700100361	700100362



Dynamic Walk - Single Side Anterior

SE

SpecificationsProduct Features

- Lightweight
- Position in neutral
- High strength
- Fits easily into normal shoes
- Warranty:12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700115121K	700115122K
Medium	700115141K	700115142K
Large	700115161K	700115162K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700115121	700115122
Medium	700115141	700115142
Large	700115161	700115162

PADDING

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700100421	700100422
Medium	700100441	700100442
Large	700100461	700100462



Dynamic Walk - Shoe Mount

SE US

Product Features

- Available with custom height
- Lightweight and high strength
- Assemble yourself or let us do it for you
- NOTE: single side lateral only
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

Dynamic Walk **with** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700110721K	700110722K
Medium	700110741K	700110742K
Large	700110761K	700110762K

Dynamic Walk **without** Padding

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700110721	700110722
Medium	700110741	700110742
Large	700110761	700110762

PADDING

SIZE	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	700100321	700100322
Medium	700100341	700100342
Large	700100361	700100362



Dynamic Walk - Modular

SE US

The Dynamic Walk Modular is an extension device that adds all the features and benefits of the Dynamic Walk to any knee orthosis. The Dynamic Walk Modular prevents distal migration of the knee orthosis while also providing some “toe pick-up” to allow for better ground clearance.

Product Features

- Select from three sizes
- Available in dorsi-assistor straight version
- Easy attachment to orthosis
- Attachment section:
Inner dimensions: 3.7 x 0.2 in. (95 x 6 mm)
Outer dimensions: 4.7 x 0.6 in. (120 x 16 mm)
- Warranty: 12 months
- Size descriptions on page 11-12



Ordering information

PRODUCT NUMBER DORSI-ASSIST

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	023889	023888
Medium	023891	023890
Large	023893	023892

PRODUCT NUMBER STRAIGHT

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Small	023895	023894
Medium	023897	023896
Large	023899	023898



Dynamic Walk - Mix & Match

SE

Product Features

- Select from model and size of foot plate
- Select from model and size of calf band
- Height is determined by model and size of calf section. See corresponding model on previous pages
- Lightweight and high strength
- Warranty: 6 months
- Size descriptions on page 11-12



Ordering information

When ordering please include the following information:

PRODUCT NUMBER

DESCRIPTION	MODEL	CALFSIZE	FOOTSIZE	SIDE
Dynamic Walk - Mix & Match - 70012	X	X	X	X
Standard	1			
Wide Ankle	2			
Single Side Medial	3			
Single Side Lateral	4			
Double-Double	5			
Small Calf		2		
Medium Calf		4		
Large Calf		6		
Small Foot			2	
Medium Foot			4	
Large Foot			6	
Left				1
Right				2

Example: Single side lateral, large calf, medium foot, right side: 700 124 642

If you need to change size on the footpart or change the height please use Custom Calf height or Foot length instead on page 21-22

FOOT SECTION	SMALL	MEDIUM	LARGE
All types except Shoe Mount	230 mm	250 mm	270 mm

CALF SECTION	SMALL	MEDIUM	LARGE
Standard/Double Double	320 mm	355 mm	370 mm
Wide Ankle Standard/Wide Ankle Double Double	335 mm	355 mm	380 mm
Single Side Lateral/Single Side Medial/ Shoe Mount	340 mm	370 mm	385 mm

Dynamic Walk - Mix & Match Custom Calf Height or Foot Length

SE

Product Features

- Select from model and size of foot plate
- Select from model and size of calf band
- Height is determined by model and size of calf section. (See charts below)
- Lightweight and high strength
- Warranty: 6 months



Custom Calf Height or Foot length

FOOT SECTION	SMALL	MEDIUM	LARGE
All types except Shoe Mount	230 mm +/- 13 mm	250 mm +/- 13 mm	270 mm +/- 13 mm
CALF SECTION	SMALL	MEDIUM	LARGE
Standard/Double Double	320 mm +/- 25 mm	355 mm +/- 30 mm	370 mm +/- 35 mm
Wide Ankle Standard/Wide Ankle Double Double	335 mm +/- 25 mm	355 mm +/- 30 mm	380 mm +/- 35 mm
Single Side Lateral/Single Side Medial/ Shoe Mount	340 mm +/- 25 mm	370 mm +/- 30 mm	385 mm +/- 35 mm

Ordering information

Use order form on page 22

DESCRIPTION	PRODUCT NUMBER
Dynamic Walk - Mix & Match Custom Calf Height or Foot Length	700100105

NOTE: Padding is sold separately.



Fillauer®

Order form - Dynamic Walk Custom Calf Height or Foot Length

Customer: Date:

Patient reference/ order number: Ordered by:

Phone: Email:

Patient information:

Age: _____ Weight (kg): _____

Activity level:

☐ Low ☐ Medium ☐ High

Side:

☐ Left ☐ Right ☐ BilateralType (choose **one**): 700 100 105

- ☐ Standard
- ☐ Standard Wide Ankle
- ☐ Double - Double
- ☐ Double - Double Wide Ankle
- ☐ Single Side Lateral
- ☐ Single Side Medial
- ☐ Shoe Mount - Custom Height

Foot Plate

- ☐ Small Increase/Decrease Length (mm)
- ☐ Medium + _____
- ☐ Large - _____

Calf Band

Height is determined by model and size of calf section. (See chart below)

- ☐ Small Increase/Decrease Height (mm)
- ☐ Medium + _____
- ☐ Large - _____


Material Options:

Foot Plate:

- ☐ Standard
- ☐ Rigid foot plate (for diabetics)
- ☐ Flexible foot plate (high activity users)

Calf band:

- ☐ Standard
- ☐ Rigid calf band (more stability)



Measurements (mm):

A _____ B _____

Explanation:

A - Height of orthosis
B - Foot length

FOOT LENGTH	SMALL	MEDIUM	LARGE
All types except Shoe Mount	230 mm +/- 13 mm	250 mm +/- 13 mm	270 mm +/- 13 mm
CALF HEIGHT	SMALL	MEDIUM	LARGE
Standard	320 mm +/- 25 mm	355 mm +/- 30 mm	370 mm +/- 35 mm
Single Side Lateral /Single Side Medial, Double Double	340 mm +/- 25 mm	370 mm +/- 30 mm	385 mm +/- 35 mm
Shoe Mount	360 mm +/- 25 mm	380 mm +/- 30 mm	400 mm +/- 35 mm

NOTE: Padding is ordered separately.

Dynamic Walk - Custom Made

SE

US

The Dynamic Walk Custom Made is recommended for patients in need of something utterly special or when the Standard, Single Side or Mix & Match features will not fit correctly. Whether the need is for a customisation of the lifting force or size, Fillauer's Dynamic Walk can be Custom Made to fit perfectly.

Warranty: 6 months

Please send us:

- Positive Cast
- Negative Cast
- Positive Foam

Product Features

- Customisable size and lift
- Perfect fit from patient cast
- Lightweight and high strength
- Order form on page 24

Ordering information

PRODUCT NUMBER

DESCRIPTION	PRODUCT NUMBER
Dynamic Walk - Custom Made	700100100
Negative Cast	700100146*
Positive Foam	700100147*

*add this product number when sending negative cast or foam

Casting/Mold Instructions

- Send a negative cast or a positive scanned mold
- The cast/mold should be taken seated with 90 degrees in the ankle and with as neutral (medio/lateral) position in the ankle joint is possible
- Mark pressure points and sensitive areas
- Please make a smooth cast with 4 layers of plaster, the plaster must be 5 cm longer than the desired height of the finished orthosis
- Reinforce the cast with 2 layers of plaster
- If you send a positive cast be sure that the rebar is so far in as some cm from the proximal end of the foot/toes
- Fill in the whole measurement form and send it together with the cast/mold
- The cast will be discarded 6 months after delivery

NOTE: Padding is sold separately.



Pediatric size Dynamic Walk



Order form - Dynamic Walk Custom Made Adult and Pediatric

Customer: Date & Ordered by:

Patient reference/ order number: (mark cast/model, trace of insert, with this reference)

Phone: Email:

Casting/Mold Instructions

Dynamic walk is not a corrective orthosis. The cast/mold should be possible to align with 90° in the ankle and close to neutral in the frontal and transverse plane, to achieve good function. If this is not achievable you need to contact Fillauer prior to order: order@fillauer.com

- If the patient uses a foot orthotic, cast/mold with the foot orthotic on.
- Send a negative/positive plaster cast or a positive foam model
- Mark pressure points and sensitive areas and if trimlines should be different from standard
- Please make a smooth cast with **4 layers** of plaster, the cast/mold must be **5 cm longer** than the desired height of the finished orthosis
- Reinforce cast if necessary
- If you send a positive cast be sure that the rebar reaches the forefoot and is sticking out >10 cm from the proximal end
- If possible, attach a trace of the shoe insert, the brace is going into
- Fill in the whole measurement form and send it together with the cast/mold
- The cast will be discarded 6 months after delivery

Patient information:

 Age: ____ Weight (kg): ____ Side: ☒ Left ☐ Right

Activity level:

☐ Low ☐ Medium ☐ High

☐ Cast/mold with foot orthotic on

Base of Production:

☐ Negative plaster cast 700 100 146 ☐ Positive plaster cast 700 100 145 ☐ Positive Foam Model 700 100 147

 Type (choose **one**): 700 100 100

- ☐
- Standard
-
- ☐
- Single Side Lateral
-
- ☐
- Single Side Medial
-
- ☐
- Double Double
-
- ☐
- Shoe Mount

Material Options:

Foot Plate:	Calf band:	
☐ Standard	☐ Standard	
☐ Rigid foot plate (for diabetics)	☐ Rigid calf band (more stability)	
☐ Flexible foot plate (high activity users)		

Peek thickness:


☐ 6 mm (Adult Standard) ☐ 5 mm (standard pediatric/less assist) ☐ 4 mm (only small children)

For help on making product selections please contact Fillauer to discuss

Individual padding is not available for custom orthotics, but material can be ordered with the orthotic to be adjusted by the clinician.

☐ Sheet of padding: 700100149

Additional information:


 Measurements: ☐ mm ☐ inches
 Desired measurements of finished brace

A ____ B ____ C ____

D ____ E ____

Explanation:

- A - Height of orthosis (usually just above the highest point of the calf)
 B - Heel height
 C - Foot plate length
 D - Foot plate width at MTP
 E - Foot plate width at the middle of the heel





Orthotic Joints

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

Modular Integrated Ankle Joint, Double Action

US

Product Features

- Sold in pairs
- Square bar slot for easy fabrication
- Nickel plated for corrosion resistance
- Available in steel or aluminum alloy
- Offers spring-loaded dorsi-assist combined with a stepless adjustable dorsi-flexion stop
- Warranty: 12 months



Ordering information

MATERIAL	PROD. NO. 12 MM BAR SLOT	PROD. NO. 16 MM BAR SLOT	PROD. NO. 19 MM BAR SLOT
Steel	IA12-EU	IA16-EU	IA19-EU
Aluminum alloy	IA12-AL-EU	IA16-AL-EU	IA19-AL-EU

Modular Integrated Ankle Joint, Adjustable Range of Motion

US

Product Features

- Sold in pairs
- Square bar slot for easy fabrication
- Available in steel or aluminum alloy
- Features anterior and posterior ROM adjustment
- Accepts double action stirrup head or solid stirrup (not included)
- Possible to initially lock ROM and then gradually increase it as patient progresses
- Warranty: 12 months



Ordering information

MATERIAL	PROD. NO. 12 MM BAR SLOT	PROD. NO. 16 MM BAR SLOT	PROD. NO. 19 MM BAR SLOT
Steel	IA12ARM-EU	IA16ARM-EU	IA19ARM-EU
Aluminum alloy	IA12ARM-AL-EU	IA16ARM-AL-EU	IA19ARM-AL-EU

Integrated Ankle Joint, Double Action

US

Product Features

- Sold in pairs
- Dorsi-flexion and/ or plantar flexion assist
- Nickel plated aluminum for corrosion resistance
- Includes all mounting hardware, stirrups and alignment fixture
- Warranty: 12 months



Ordering information

DESCRIPTION	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
Integrated Ankle Joint - Double Action	IA-CH-DA-AL	IA-YTH-DA-AL	IA-AD-DA-AL

Integrated Ankle Joint Adjustable range of motion

US

Product Features

- Sold in pairs
- Uses a double action stirrup
- Stepless adjustability
- Includes all mounting hardware stirrups and alignment fixture
- Aluminum
- Warranty: 12 months



Ordering information

DESCRIPTION	PROD. NO. SMALL	PROD. NO. LARGE
Integrated Ankle Joint - Adjustable Range of Motion	IA-SM-ARM-AL	IA-LG-ARM-AL

Integrated Ankle Joint Dorsiflexion Assist


US

Product Features

- Sold in pairs
- Low profile, lightweight alternative to modular system
- Includes all mounting hardware, stirrups and alignment fixture
- Aluminum
- Warranty: 12 months

Ordering information

DESCRIPTION	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
Integrated Ankle Joint- Dorsiflexion Assist	IA-CH-DF-AL	IA-YTH-DF-AL	IA-AD-DF-AL



Split Tail Stirrups

US

Product Features

- Unique “Y” shape of the split-tail is perfectly suited for use in plastic NYUCB shoe insert configurations
- Special order drilled or drilled and tapped
- Sold in pairs
- Warranty: 12 months



Ordering information

DESCRIPTION	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
Split Tail Stirrups	IAR-CH-DA-ST	IAR-YTH-DA-ST	IAR-AD-DA-ST

Solid Stirrups - Wide Tongue

US

Product Features

- For use with double action ankle joints and adjustable range of motion ankle joints
- Warranty: 12 months



Ordering information

SIZE	WIDTH	PRODUCT NUMBER
Adult	9.5 in. (241 mm)	ST-WTDA9.5
Adult	10.25 in. (260 mm)	ST-WTD10.25
Adult	11 in. (279 mm)	ST-WTDA11

Short Leg Upright - Double Action

US

Product Features

- Aluminum
- For lamination
- Sold individually
- Length: 16 in. (40.6 cm)
- Product weight: 4.5 oz. (128 g)
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
Large	020518
Large with Positioning Pins	020519
SPARE PARTS LARGE	SPARE PARTS LARGE WITH POSITION PINS
010406 - Short Leg Double Action Large Bar	010406 - Short Leg Double Action Large Bar
880332 - Set Screw, 9/32"-24, Socket	880334 - Positioning Pin Threaded
880730 - Joint Nut, 1/4"-28	880730 - Joint Nut, 1/4"-28
880731 - Joint Screw, 1/4"-28	880731 - Joint Screw, 1/4"-28
881301 - Ball, 1/4 in. Diameter Chrome	
881408 - Compression Spring, .2250.0 × 1.3 Long	

Short Leg Upright - limited/ free motion

US

Product Features

- Aluminum
- Sold individually
- Length small: 8.5 in. (21.6 cm)
- Length medium: 12.5 in. (31.8 cm)
- Length large: 17.0 in. (43.0 cm)
- Warranty: 12 months



Ordering information

DESCRIPTION	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
Limited motion	-	020514	020515
Free motion	020509	020510	020511
SPARE PARTS LIMITED MOTION	SPARE PARTS FREE MOTION		
010416 Small Short Leg Limited Motion Brace	010417 Small Short Leg Free Motion Brace		
010410 Medium Short Leg Limited Motion Brace	0104101 Medium Short Leg Free Motion Brace		
010404 Large Short Leg Limited Motion Brace	010405 Large Short Leg Free Motion Brace		
880731 Joint Screw, 1/4"-28 × 0.33 Long	020804 Screw, 10-32 × 1/4 in., Truss Head Slotted		
880730 Joint Nut, 1/4"-28 × 0.280	020805 Joint Nut, 10-32 × 0.217		

Short Leg Upright - Toe Pick Up

US

Product Features

- Aluminum
- Sold individually
- Length large: 17.0 in. (43.0 cm)
- Warranty: 12 months



Ordering information

SIZE	PROD. NO. LEFT	PROD. NO. RIGHT
Small	020501	020500
Large	020506	020504
SPARE PARTS SMALL		SPARE PARTS LARGE
010414	Right Short Leg Toe Lift Brace	010402 Right Short Leg Toe Lift Brace
010415	Left Short Leg Toe Lift Brace	010403 Left Short Leg Toe Lift Brace
880330	Socket Set Screw, 1/4-28	880332 Socket Set Screw, 9/32-24
881405	Spring, 0.208 × 1.17 in. Long	881404 Spring, 0.230 × 1.8 in. Long
881300	Ball, 7/32 in. Diameter Chrome	881301 Ball, 1/4 in. Diameter Chrome
020805	Joint Nut, 10-32 x 0.217 in.	880730 Joint Nut, 1/4"-28 x 0.280
020804	Joint Screw, 10-32 x 1/4 in. Truss Head	880731 Joint Screw, 1/4"-28 x 0.33

Ankle Joint - Limited/ Free Motion

US

Product Features

- Sold each
- Stainless steel
- Warranty: 12 months

Ordering information

LIMITED MOTION

DESCRIPTION	PRODUCT NUMBER
Limited Motion, 5/8 in. (16 mm) bar slot	020862
Limited Motion, 3/4 in. (19 mm) bar slot	020866

SPARE PARTS

010438	3/4 in. (19 mm) Ankle, Joint Limited Motion	880731	Joint Screw, 1/4"-28 × 0.33 in.
010432	5/8 in. (16 mm) Ankle Joint Limited Motion	880730	Joint Nut, 1/4"-28 × 0.28 in.

FREE MOTION

DESCRIPTION	PRODUCT NUMBER
Free Motion, 5/8 in. (16 mm) bar slot	020872
Free Motion, 3/4 in. (19 mm) bar slot	020876

SPARE PARTS

010433	5/8 in. (16 mm) Ankle Joint Free Motion	880731	Joint Screw, 1/4"-28 × 0.33 in.
010437	3/4 in. (19 mm) Ankle Joint, Free Motion	880730	Joint Nut, 1/4"-28 × 0.28 in.



Ankle Joint - Toe Pick Up

US

Product Features

- Sold each
- Openings for a 5/8 in. (16 mm) wide bar
- Stainless steel
- Warranty: 12 months

Ordering information

DESCRIPTION	PRODUCT NUMBER LEFT	PRODUCT NUMBER RIGHT
Toe Pick Up, 5/8 in. (16 mm) bar slot	020855	020854

SPARE PARTS

010430	Right Ankle Joint Toe Pick Up	880731	Joint Screw, 1/4"-28 × 0.33 in.
010431	Left Ankle Joint Toe Pick Up	880730	Joint Nut, 1/4"-28 × 0.28 in.
881301	Ball, 1/4 in. Diameter, Chrome	881406	Spring, 0.23 × 1.4 in. Long
880332	Socket Set Screw, 9/32-24 Flat Point		



Ankle Joint - Double Action

US

Product Features

- Sold each
- Stainless steel
- Warranty: 12 months

Ordering information



DOUBLE ACTION 5/8 IN. AND 3/4 IN.

DESCRIPTION	PRODUCT NUMBER
Double Action, 5/8 in. (16 mm) Channel	020882
Double Action, 3/4 in. (19 mm) Channel	020886

SPARE PARTS

010434	5/8 in. (16 mm) Model Ankle Joint, Mach Double Action	880730	Joint Nut, 1/4"-28 × 0.28
010439	3/4 in. (19 mm) Model Ankle Joint, Mach Double Action	881301	Ball, 1/4 in. Diameter Chrome
880332	Socket Set Screw, 9/32-24 Flat Point	881408	Spring , 0.225 × 1.36 in. Long
880731	Joint Screw, 1/4"-28 × 0.33 in. Long		

DOUBLE ACTION THREADED PINS AND POSITIONING PINS

DESCRIPTION	PRODUCT NUMBER
Double Action with Threaded Pins, 5/8 in. (16 mm) Channel	020884
Double Action with Positioning Pins 3/4 in. (19 mm) Channel	020888

SPARE PARTS

010434	5/8 in. (16 mm) Model Ankle Joint, Mach Double Action	880731	Joint Screw, 1/4"-28 × 0.33 in.
010439	3/4 in. (19 mm) Model Ankle Joint, Mach Double Action	880730	Joint Nut, 1/4"-28 × 0.28 in.
880334	Positioning Pin, with Thread		

Stirrups - Standard

US

Product Features

- Sold each
- Has an opening for a 5/8 in. (16 mm) wide bar
- Warranty: 12 months

LM/FM



WIDTH	PROD. NO
8-3/4 in. (222 mm)	010816
9-1/2 in. (241 mm)	010817
10-1/4 in. (260 mm)	010818

Double Action



WIDTH	PROD. NO
8-3/4 in. (222 mm)	010830
9-1/2 in. (241 mm)	010831
10-1/4 in. (260 mm)	010832
11 in. (279 mm)	010878

Toe Pick Up



WIDTH	PROD. NO
8-3/4 in. (222 mm)	010802
9-1/2 in. (241 mm)	010803
10-1/4 in. (260 mm)	010804
11 in. (279 mm)	010877

Stirrups - Split Tongue

US

Product Features

- Sold each
- Warranty: 12 months

LM/FM



WIDTH	PROD. NO
1/8 x 5/8 in. (3.2 x 16 mm)	010863
1/8 x 3/4 in. (3.2 x 19 mm)	010864
1/8 x 7/8 in. (3.2 x 22 mm)	010865

Double Action



WIDTH	PROD. NO
1/8 x 5/8 in. (3.2 x 16 mm)	010871
1/8 x 3/4 in. (3.2 x 19 mm)	010872
1/8 x 7/8 in. (3.2 x 22 mm)	010873

Toe Pick Up



WIDTH	PROD. NO
1/8 x 5/8 in. (3.2 x 16 mm)	010855
1/8 x 3/4 in. (3.2 x 19 mm)	010856
1/8 x 7/8 in. (3.2 x 22 mm)	010857

Stirrups

Long Tongue & Wide Tongue

US

Product Features

- Sold each
- Has an opening for a 16 mm wide bar
- Warranty: 12 months

Long Tongue



LM/ FM LARGE

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010842

DOUBLE ACTION LARGE

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010845

TOE PICK UP

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010839

Wide Tongue



LM/ FM LARGE

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010851

DOUBLE ACTION LARGE

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010854

TOE PICK UP

WIDTH	PROD. NO
10-1/4 in. (260 mm)	010848

Rectangle Caliper Plates

US

Product Features

- Zinc plated
- Sold each



Ordering information

SIZE	LENGTH	PROD. NO. LEFT	PROD. NO. RIGHT
1/8 x 5/8 in. Uprights	3.5 in. (88.9 mm)	020248	020230
1/8 x 3/4 in. Uprights	3-3/4 in. (95.3 mm)	020263	020255

Insert Stirrup For Plastic Orthotics

US

Product Features

- Stainless steel
- Sold in pair, medial and lateral
- For plastic laminations
- Warranty: 12 months



without offset

Ordering information

INSERT STIRRUP WITHOUT OFFSET

DESCRIPTION	2-3/8 IN. (60 MM)	2-3/4 IN. (70MM)
Stirrup, Insert, Dubble Action	020096	020094

Polycar-C™ Ankle Inserts

US

Product Features

- High performance reinforced composite
- Heating temp. 350°-400°F (175°C - 200°C)
- Easy to mold
- High carbon content
- Thermobonds to polypropylene, copolymer and metal

Ordering information

DESCRIPTION	SIZE	PRODUCT NUMBER
PolyCar-C™ insert, 0.06 in. (1.5 mm)	Small	076630
PolyCar-C™ insert, 0.09 in. (2.3 mm)	Medium	076632
PolyCar-C™ insert, 0.09 in. (2.3 mm)	Large	076634
PolyCar-C™ insert, 0.13 in. (3.2 mm)	X-Large	076636



Dorsiflexion Assist Ankle Joint Kits With Polycar-C

US

Specifications:

Weight limit: 181 lbs. (82 kg)
Warranty: 12 months

Product Features

- Indicated for conditions requiring toe lift
- Aluminum head, stainless steel, stirrup, right and left, pair

Ordering information

DESCRIPTION		PRODUCT NUMBER	
Dorsiflexion Assist Kit		020836	
SPARE PARTS			
A	020830 - Dorsiflexion ankle joint body, right	F	020802 - Pre-cut PolyCar-C™, .060 in. composite, 2 each
	020834 - Dorsiflexion ankle joint body, left	G	020804 - Screw, 10-32 x 1/4 in., truss HD, SS, 2 each
B	020760 - AFO alignment tool, 1 each	H	020832 - Dorsiflexion stirrup insert, right, 1 each, LG
C	880330 - Socket set screw 1/4"-28 in.		020828 - Dorsiflexion stirrup insert, left, 1 each, LG
D	020807 - Spring, 0.21 in. O.D., 2 each	I	020805 - Joint nut, 2 each 10-32 x 0.217 in.
E	881300 - Ball, chrome, 7/32 in., 2 each	J	885345 - Allen wrench, 1/8 in. short



PDA Ankle Joint Kits With Polycar-C

US

Specifications:

Weight limit: 181 lbs. (82 kg)

Warranty: 12 months

Product Features

- Plantar/ dorsiflexion assist
- Indicated for conditions requiring rigid ankle support and positioning
- Aluminum head, stainless steel, stirrup, right and left, pair

Ordering information

DESCRIPTION	PRODUCT NUMBER
PDA Ankle Joint Kit, large	020894
SPARE PARTS	
A 020890 - PDA ankle joint body, 2 ea.	G 020807 - Spring, .21 in. O.D., 2 each
B 020760 - AFO alignment tool, 1 ea,	H 809763 - Dowel pin, 1/8 x 7/8 in., 4 each
C 020802 - Pre-cut polycar-C™, .060 in. composite, 2 each	I 880330 - Spring retainer screw, 2 each
D 020804 - Screw, 10-32 x 1/4 in., truss HD, SS, 2 each	J 020826 - Double action stirrup insert, 2 each, Large
E 020805 - Joint nut, 2 each 10-32 x 0.217 in.	K 885345 - Allen wrench, 1/8 in.
F 881300 - Ball, chrome, 7/32 in., 2 each	





Steplock®

US

Product Features

- Steplock® orthotic joint offers the patient security even though full extension can not be reached
- Perfect joint for treating a flexion contracture
- The joint has 9 different positions, locks in flexion and allows movement in extension
- Release kit available, see page 41 - 42
- Nickel plated alloy
- Straight - medial and lateral straight
- Left medial offset (pairs) - left medial contoured, lateral straight
- Right medial offset (pairs) - right medial contoured, lateral straight
- Double contour - medial and lateral contoured
- Warranty: 12 months

Ordering information

JOINT PAIRS WITH BARS

DESCRIPTION	3/16 X 1/2 IN. (5 X 12 MM)	3/16 X 5/8 IN. (5 X 16 MM)	1/4 X 5/8 IN. (6 X 16 MM)	3/16 X 3/4 IN. (5 X 19 MM)	1/4 X 3/4 IN. (6 X 19 MM)
Straight	SL-512KS	SL-516KS	SL-616KS	SL-519KS	SL-619KS
Left medial contour	SL-512KL	SL-516KL	SL-616KL	SL-519KL	SL-619KL
Right medial contour	SL-512KR	SL-516KR	SL-616KR	SL-519KR	SL-619KR
Double contour	SL-512KB	SL-516KB	SL-616KB	SL-519KB	SL-619KB
Hip joints	SL-512H	SL-516H	SL-616H	SL-519H	SL-619H

JOINT PAIRS WITHOUT BARS

DESCRIPTION	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT	3/4 IN. (19 MM) BAR SLOT
Straight	SL-12JS	SL-58JS	SL-34JS
Left medial contour	SL-12JL	SL-58JL	SL-34JL
Right medial contour	SL-12JR	SL-58JR	SL-34JR
Double contour	SL-12JB	SL-58JB	SL-34JB

SINGLE JOINT WITH BARS

DESCRIPTION	3/16 X 1/2 IN. (5 X 12 MM)	3/16 X 5/8 IN. (5 X 16 MM)	1/4 X 5/8 IN. (6 X 16 MM)	3/16 X 3/4 IN. (5 X 19 MM)	1/4 X 3/4 IN. (6 X 19 MM)
Left medial contour	SL-512KOL	SL-516KOL	SL-616KOL	SL-519KOL	SL-619KOL
Right medial contour	SL-512KOR	SL-516KOR	SL-616KOR	SL-519KOR	SL-619KOR
Straight left lateral	SL-512KSL	SL-516KSL	SL-616KSL	SL-519KSL	SL-619KSL
Straight right lateral	SL-512KSR	SL-516KSR	SL-616KSR	SL-519KSR	SL-619KSR

SINGLE JOINT WITHOUT BARS

DESCRIPTION	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT	3/4 IN. (19 MM) BAR SLOT
Left medial contour	SL-12JOL	SL-58JOL	SL-34JOL
Right contour	SL-12JOR	SL-58JOR	SL-34JOR
Straight left lateral	SL-12JSL	SL-58JSL	SL-34JSL
Straight right lateral	SL-12JSR	SL-58JSR	SL-34JSR



Interlock™

US

Product Features

- Interlock is a drop lock alternative that is modular
- Release kit available, see page 41 - 42
- Nickel plated alloy
- Straight - medial and lateral straight
- Left medial offset (pairs) - left medial contoured, lateral straight
- Right medial offset (pairs) - right medial contoured, lateral straight
- Double contour - medial and lateral contoured
- Warranty: 12 months

Ordering information

JOINT PAIRS WITH BARS

DESCRIPTION	3/16 X 1/2 IN. (5 X 12 MM)	3/16 X 5/8 IN. (5 X 16 MM)	1/4 X 5/8 IN. (6 X 16 MM)	3/16 X 3/4 IN. (5 X 19 MM)	1/4 X 3/4 IN. (6 X 19 MM)
Straight	IL-512KS	IL-516KS	IL-616KS	IL-519KS	IL-619KS
Left Medial Contour	IL-512KL	IL-516KL	IL-616KL	IL-519KL	IL-619KL
Right Medial Contour	IL-512KR	IL-516KR	IL-616KR	IL-519KR	IL-619KR
Double Contour	IL-512KB	IL-516KB	IL-616KB	IL-519KB	IL-619KB
Hip Joints	IL-512H	IL-516H	IL-616H	-	-

JOINT PAIRS WITHOUT BARS

DESCRIPTION	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT	3/4 IN. (19 MM) BAR SLOT
Straight	IL-12JS	IL-58JS	IL-34JS
Left Medial Contour	IL-12JL	IL-58JL	IL-34JL
Right Medial Contour	IL-12JR	IL-58JR	IL-34JR
Double Contour	IL-12JB	IL-58JB	IL-34JB

SINGLE JOINT WITH BARS

DESCRIPTION	3/16 X 1/2 IN. (5 X 12 MM)	3/16 X 5/8 IN. (5 X 16 MM)	1/4 X 5/8 IN. (6 X 16 MM)	3/15 X 3/4 IN. (5 X 19 MM)	1/4 X 3/4 IN. (6 X 19 MM)
Left Medial Contour	IL-512KOL	IL-516KOL	IL-616KOL	IL-519KOL	IL-619KOL
Right Medial Contour	IL-512KOR	IL-516KOR	IL-616KOR	IL-519KOR	IL-619KOR
Straight Left Lateral	IL-512KSL	IL-516KSL	IL-616KSL	IL-519KSL	IL-619KSL
Straight Right Lateral	IL-512KSR	IL-516KSR	IL-616KSR	IL-519KSR	IL-619KSR

SINGLE JOINT WITHOUT BARS

DESCRIPTION	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT	3/4 IN. (19 MM) BAR SLOT
Left Medial Contour	IL-12JOL	IL-58JOL	IL-34JOL
Right Medial Contour	IL-12JOR	IL-58JOR	IL-34JOR
Straight Left Lateral	IL-12JSL	IL-58JSL	IL-34JSL
Right Lateral	IL-12JSR	IL-58JSR	IL-34JSR






2400 Joints

US

Product Features

- Low profile lines give the 2400 series less weight and a sleeker appearance
- Sold without uprights
- Available in single position full extension Unilock, ratcheting Steplock® or adjustable stop free-motion
- Notched neck on distal for easy to bend custom medial/ lateral contours, and fewer parts to stock
- Built-in lock retaining feature in both minimal and posterior offset
- Release kit is available, see page 42
- Warranty: 12 months



Ordering information

2400 JOINT PAIRS	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT	3/4 IN. (19 MM) BAR SLOT
UniLock (UL) Minimal Offset (MO)	2400-ULMO5K	2400-ULMO6K	2400-ULMO7K
UniLock (UL) Posterior Offset (PO)	2400-ULPO5K	2400-ULPO6K	2400-ULPO7K
StepLock® (SL) Minimal Offset (MO)	2400-SLMO5K	2400-SLMO6K	2400-SLMO7K
StepLock® (SL) Posterior Offset (PO)	2400-SLPO5K	2400-SLPO6K	2400-SLPO7K
Free Motion Knee Minimal Offset (MO)	2400-FM-MO-5K	2400-FM-MO-6K	2400-FM-MO-7K
Free Motion Knee Posterior Offset (PO)	2400-FM-PO-5K	2400-FM-PO-6K	2400-FM-PO-7K
Reverse Steplock®	2400-R-SLMO5-PR	2400-R-SLMO6-PR	-

REVERSE STEPLOCK® SINGLE JOINT

SIZE	1/2 IN. (12 MM) BAR SLOT	5/8 IN. (16 MM) BAR SLOT
Left	2400-R-SLMO5L-KT	2400-R-SLMO6L-KT
Right	2400-R-SLMO5R-KT	2400-R-SLMO6R-KT



Reverse
Steplock®

Bar Kit

US

Product Features

- Comes in a four piece kit
- Material: Aluminum



Ordering information


DESCRIPTION	3/16 X 1/2 IN. (5 X 12 MM)	3/16 X 5/8 IN. (5 X 16 MM)	3/16 X 3/4 IN. (5 X 19 MM)	1/4 X 5/8 (6 X 16 MM)	1/4 X 3/4 (6 X 19 MM)
Bar Kit	BARKIT512	BARKIT516	BARKIT519	BARKIT616	BARKIT619

2400 Bar Kit

US

Product Features

- 3/16 in. (5 mm) thick bar
- Comes in a four piece kit
- Material: Aluminum



Ordering information

DESCRIPTION	PRODUCT NUMBER
3/16 x 1/2 in. (5 x 12 mm)	2400-BAR-KIT-5
1/4 x 5/8 in. (6 x 16 mm)	2400-BAR-KIT-6
3/16 x 3/4 in. (7 x 19 mm)	2400-BAR-KIT-7

Close Contour Pediatric Knee Joint

US

Product Features

- Miniature pediatric Steplock® knee joint
- 120° range of motion with 9 intermediate locking positions
- Rugged nickel plated alloy construction
- Can be contoured within 1/2 in. (12 mm) of knee center
- Innovative lock retainer allows free motion
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
Close Contour Joint Pairs - Steplock, 4.5	CC4-5SL

Release Kit SBK

US

Product Features

- For use on original series Steplock® only
- Features a soldered stainless steel bail (joints, uprights, plastics not included)
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
Bail Kit for Regular Step Lock, IL	SBK

Release Kit RBK

US

Product Features

- For use on original series Interlock and Steplock
- Features a stainless steel rod bail (joints, uprights, plastics not included)
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
Bail Kit Rod for Steplock + IL	RBK

Release Kit - UBK

US

Product Features

- For use on 2400/ 2755 series Unilock
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
.IMAC Bail Kit 2400/2755.Unilock Knee.Joint.....	UBK

Lever Release Kit - LRK

US

Product Features

- The lever release kit provides a remote-operated unlocking mechanism that connects the medial and lateral joints, making it easier to lock and unlock any spring-loaded joint (joints, uprights and plastics not included)
- Easy to fabricate
- Strong
- Warranty: 12 months



Ordering information

DESCRIPTION	PRODUCT NUMBER
. Fits the Interlock and Steplock size 3/16 x 1/2 in. (5 x 12 mm)	LRK512SM
. Fits the Interlock and Steplock all size except size 3/16 x 1/2 in. (5 x 12 mm)	LRK1619

Close Contour


US

Product Features

- Elbow - forward extension ratcheting and reverse flexion ratcheting for positioning
- Wrist ratcheting from 90° plantar flexion to 30° dorsi flexion
- Includes free motion bar joint for opposite side
- Steplock features a built-in lock retainer
- Unparalleled strength to weight ratio with OTS alloy
- Nickel plated finish for durability
- Warranty: 12 months

Ordering information

DESCRIPTION	PROD. NO. LEFT LATERAL	PROD. NO. RIGHT LATERAL
Extension Elbow	CC-4ESL-LT	CC-4ESL-RT
Reverse Elbow (flexion)	CC-4RSL-LT	CC-4RSL-RT
Wrist	CC-4WSL-LT	CC-4WSL-RT





Trigger Release for Camlocks

US


Product Features

- Proximal release for sitting
- Easy use from standing position
- Great for RGO's
- For medium and large camlock joints




Ordering information


DESCRIPTION	PRODUCT NUMBER
Trigger Release, with Dual Cable	023750
Trigger Release, with Bail Rod	023751
Trigger Release, with Single Cable	023753



023750



023751



023753

Reference pages 47-48 Cam Locks

Zippy Quick Release

US

Product Features

- Allows easy removal and reattachment from AFO upper assemblies



Ordering information

DESCRIPTION	PRODUCT NUMBER
Zippy Quick Release, with Cleat	020700
Zippy Quick Release, with Threaded Mounting Holes in Cleat	020710

Drop Lock Aluminum Hip Joint Assemblies

US

Product Features

- Sold in pairs
- Warranty: 12 months



Ordering information

DESCRIPTION	BAR	PRODUCT NUMBER
Small	3/16 x 1/2 in. (5 x 13 mm)	050401
Medium	1/4 x 5/8 in. (6 x 16 mm)	050402
Large	1/4 x 3/4 in. (6 x 19 mm)	050403

Drop Lock Aluminum Knee Joint Assemblies

US

Product Features

- Sold in pairs
- Warranty: 12 months



Ordering information

DESCRIPTION	BAR	PRODUCT NUMBER
Small	3/16 x 1/2 in. (5 x 13 mm)	023800
Medium	1/4 x 5/8 in. (6 x 16 mm)	023801
Large	1/4 x 3/4 in. (6 x 19 mm)	023802



Drop Lock Low Profile Aluminum Knee Joint Assemblies

US

Product Features

- Sold in pairs
- Stainless steel joint
- Aluminum bars
- Warranty: 12 months

Ordering information

MEDIUM JOINT 3/16 X 5/8 (5 X 16 MM) BAR

DESCRIPTION	PROD. NO. LEFT	PROD. NO. RIGHT	PROD. NO. LEFT & RIGHT
Straight	-	-	023134
Medial And Lateral Contoured	-	-	023136
Medial Contoured	023138	023139	-

LARGE JOINT 3/16 X 3/4 (5 X 19 MM) BAR

DESCRIPTION	PROD. NO. LEFT	PROD. NO. RIGHT	PROD. NO. LEFT & RIGHT
Straight	-	-	023140
Medial And Lateral Contoured	-	-	023142
Medial Contoured	023144	023145	023143



Spring Loaded Lever Lift

US

Product Features

- For use with Drop Lock. Increases lever arm for easy unlocking

Ordering information



DESCRIPTION	PRODUCT NUMBER
Complete Assembly	029033
SPARE PARTS	
Rod- for spring loaded lever	029009
Spring- for spring loaded lever	029017
Guide Post- for spring loaded lever	029025

RGO & IRGO

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain



RGO/IRGO Overview

An orthotic solution that provides ambulation for the paraplegic injured who is unable to stand and/or walk due to total or partial loss of function of the lower extremities. All designs can be used with either internal or external KAFO's

Applications: Child through adult



Hoop Cable Design

- With its familiarity and simplicity, the Hooped Cable Design is the simplest and easiest to use and adjust.
- The Hooped Cable Design offers greater friction in use than the other designs.
- Popular choice for spinal lesions.
- Page 51 - 58

Horizontal Cable Design

- This design offers cosmetically most discretion in a low maintenance style. The Horizontal Cable design can be used with either a thermoplastic or a rigid pelvic band of standard or butterfly style.
- The design is often combined with a molded plastic LSO (lumbo-sacral orthosis) or TLSO (thoraco-lumbo-sacral orthosis) jacket, depending on the level of the lesion or control desired.
- Page 59 - 66

IRGO® Design

- The Isocentri RGO (IRGO) is the heavy duty alternative. It offers a very strong mechanism with a thicker Rocker Bar reciprocator that is welded to the sidebars.
- Maintenance free with very low friction makes it a very good choice for the ones with no control or very low control of the lower extremities.
- Page 67 - 69

Small Hooped Cable RGO

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Small XLong includes 3/16 x 5/8 x 9 in. (4.8 x 15.9 x 228.6 mm) Lower Bars	028040

SPARE PARTS


- A 027946 - Joint Screw, 1/4"-28 x 11/16 in., 2 Each
- B 027979 - Thrust Bearing
- C 027987 - Thrust Race
- D **Lower Bar Assembly, Standard or XLong Bars, 1 Pair Included**
- Standard, Right - 026880
- Standard, Left - 026881
- XLong, Right - 026884
- XLong, Left - 026885
- E **Coupling Plate, 1 Pair Included**
- 026427 - Coupling Plate, Right
- 026435 - Coupling Plate, Left
- F 026708 - Spacer
- G 027961 - Set Screw, 8-32 x 3/8 in., 2 Each
- H 027904 - Upper Hip Joint Bar
- I 027953 - Joint Nut
- J 027725 - Swivel Nut
- K 027727 - Swivel Coupling
- L 027860 - Screw, 5-40 x .203 in., 4 Each
- M 028888 - Cable, with Threaded Terminals, Anterior, 18 in. (457 mm)
- N 028914 - Cable, Housing, Anterior, 14 in. (356 mm)
- O 028886 - Cable, with Threaded Terminals, Posterior, 16 in. (406 mm)
- P 028912 - Cable Housing, Posterior, 12 in. (305 mm)
- Q 885140 - Hex Nut, 5-40
- R. **028050KT - Conduit Bracket Kit**
026666 - Posterior Conduit Bracket Assembly, 2 Each
- R. 028605 - Posterior Conduit Bracket
- S. 026478 - Conduit Bracket Clamp
- U. 886870 - Rivets, Stainless Steel 1/8 x 5/8 in. Truss Head, 8 Each
- V. 880574 - Conduit Bracket Screws, 6-32 x 3/8 in., 4 Each
026674 - Anterior Conduit Bracket Assembly, 2 Each
- W. 028613 - Anterior Conduit Bracket
- X. 026478 - Conduit Bracket Clamp
- Y. 880048 - Flat Head Screw, 8-32 x 1/2 in., 2 Each
- V. 880574 - Conduit Bracket Screws, 6-32 x 3/8 in., 4 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml) Tube
- 027420 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)



R. Conduit Bracket Kit





Large Hooped Cable RGO without Abduction

US

Ordering information

DESCRIPTION

Large XLong, includes 1/2 x 7/8 x 12 in. (12.7 x 22.2 x 304.8 mm) Lower Bars

PRODUCT NUMBER

028046


SPARE PARTS

- A. 027771 - Joint Screw, 7/16 x .843 in., 2 Each
- B. 027797 - Thrust Bearing
- C. 027805 - Thrust Race
- D. **Lower Bar Assembly, Standard or XLong Bars, 1 Pair**
 - 026900 - XLong, Right
 - 026901 - XLong, Left
- E. **Coupling Plate, 1 Pair Included**
 - 026601 - 2 Position Coupling Plate, Right
 - 026603 - 2 Position Coupling Plate, Left
- F. 028647 - Spacer
- G. 027896 - Set Screw, 10-32 x 3/8 in., 2 Each
- H. 027730 - Upper Hip Joint Bar
- I. 028399 - Joint Nut
- J. 027726 - Swivel Nut
- K. 027728 - Swivel Coupling
- L. 027862 - Screw, 6-32 x .350 in., 4 Each
- M. 885100 - Nut
- N. 028830 - Cable, with Threaded Terminals, Anterior, 30 in. (762 mm)
- O. 028866 - Cable Housing, Anterior, 26 in. (660 mm)
- P. 028828 - Cable, with Threaded Terminals, Posterior, 28 in. (711 mm)
- Q. 028864 - Cable Housing, Posterior, 24 in. (610 mm)
- R. **028056KT - Conduit Bracket Kit**
 - 026864 - Posterior Conduit Bracket Assembly, 2 Each**
 - S. 026628 - Posterior Conduit Bracket
 - T. 027877 - Conduit Bracket Clamp
 - U. 028936 - Support Bar
 - V. 886872 - Rivets, Stainless Steel 5/32 x 5/8 in. Truss Head, 8 Each
 - W. 880809 - Screws, 8-32 x 1/2 in., 4 Each
 - X. 026660 - Conduit Bracket Plate (Not Pictured)
 - 026674 - Anterior Conduit Bracket Assembly, 2 Each**
 - W. 028613 - Anterior Conduit Bracket
 - X. 026478 - Conduit Bracket Clamp
 - Y. 880048 - Flat Head Screw, 8-32 x 1/2 in., 2 Each
 - V. 880574 - Screws, 6-32 x 3/8 in., 4 Each


Also Included But Not Pictured:

889002 - Loctite #242, .02 oz. (0.6 ml) Tube

027420 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)



R. Conduit Bracket Kit



Large Hooped Cable RGO with Abduction

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Large, includes 1/2 x 3/4 x 11-3/8 in. (12.7 x 19.1 x 288.9 mm) Lower Bars	026162

SPARE PARTS

- A. 027771 - Joint Screw, 7/16 x .843 in., 2 Each
- B. 027797 - Thrust Bearing
- C. 027805 - Thrust Race
- D. **Lower Bar Assembly, RGO II System with Abduction, 1 Pair**
 - 026196 - RGO II System, Right
 - 026198 - RGO II System, Left
- E. **Coupling Plate, 1 Pair Included**
 - 026601 - 2 Position Coupling Plate, Right
 - 026603 - 2 Position Coupling Plate, Left
- F. 028647 - Spacer
- G. 027896 - Set Screw, 10-32 x 3/8 in., 2 Each
- H. 027730 - Upper Hip Joint Bar
- I. 028399 - Joint Nut
- J. 027726 - Swivel Nut
- K. 027728 - Swivel Coupling
- L. 027862 - Screw, 6-32 x .350 in., 4 Each
- M. 885100 - Nut
- N. 028830 - Cable, with Threaded Terminals, Anterior, 30 in. (762 mm)
- O. 028866 - Cable Housing, Anterior, 26 in. (660 mm)
- P. 028828 - Cable, with Threaded Terminals, Posterior, 28 in. (711 mm)
- Q. 028864 - Cable Housing, Posterior, 24 in. (610 mm)
- R. **028056KT - Conduit Bracket Kit**
 - 026864 - Posterior Conduit Bracket Assembly, 2 Each**
 - S. 026628 - Posterior Conduit Bracket
 - T. 027877 - Conduit Bracket Clamp
 - U. 028936 - Support Bar
 - V. 886872 - Rivets, Stainless Steel 5/32 x 5/8 in., Truss Head, 8 Each
 - W. 880809 - Screws, 8-32 x 1/2 in., 4 Each
 - X. 026660 - Conduit Bracket Plate (Not Pictured)
 - 026674 - Anterior Conduit Bracket Assembly, 2 Each**
 - W. 028613 Anterior Conduit Bracket
 - X. 026478 Conduit Bracket Clamp
 - Y. 880048 Flat Head Screw, 8-32 x 1/2 in., 2 Each
 - V. 880574 Screws, 6-32 x 3/8 in., 4 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml) Tube
- 027420 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)



R. Conduit Bracket Kit



Hooped Cable RGO Cables and Housings

US

Ordering information



SMALL HIP JOINTS, 3/32 IN. (2.4 MM) CABLE WITH TERMINALS

CABLE	HOUSING	PROD. NO. CABLE	PROD. NO. HOUSING
13 in. (330mm)	9 in. (229 mm)	028883	028909
14 in. (356 mm)	10 in. (254 mm)	028884	028910
15 in. (381 mm)	11 in. (279 mm)	028885	028911
16 in. (406 mm)	12 in. (305 mm)	028886	028912
17 in. (432 mm)	13 in. (330 mm)	028887	028913
18 in. (457 mm)	14 in. (356 mm)	028888	028914
19 in. (483 mm)	15 in. (381 mm)	028889	028915
20 in. (508 mm)	16 in. (406 mm)	028890	028916
21 in. (533 mm)	17 in. (432 mm)	028891	028917
22 in. (559 mm)	18 in. (457 mm)	028892	028918
23 in. (584 mm)	19 in. (483 mm)	028893	028919
24 in. (610 mm)	20 in. (508 mm)	028894	028920
25 in. (635 mm)	21 in. (533 mm)	028895	028921
26 in. (660 mm)	22 in. (559 mm)	028896	028922

MEDIUM AND LARGE HIP JOINTS, 1/8 IN. (3.2 MM) CABLE WITH TERMINALS

CABLE	HOUSING	PROD. NO. CABLE	PROD. NO. HOUSING
12 in. (305 mm)	8 in. (203 mm)		028908
15 in. (381 mm)	11 in. (279 mm)	028815	
16 in. (406 mm)	12 in. (305 mm)	028816	028852
17 in. (432 mm)	13 in. (330 mm)	028817	028853
18 in. (457 mm)	14 in. (356 mm)	028818	028854
19 in. (483 mm)	15 in. (381 mm)	028819	028855
20 in. (508 mm)	16 in. (406 mm)	028820	028856
21 in. (533 mm)	17 in. (432 mm)	028821	028857
22 in. (559 mm)	18 in. (457 mm)	028822	028858
23 in. (584 mm)	19 in. (483 mm)	028823	028859
24 in. (610 mm)	20 in. (508 mm)	028824	028860
25 in. (635 mm)	21 in. (533 mm)	028825	028861
26 in. (660 mm)	22 in. (559 mm)	028826	028862
27 in. (686 mm)	23 in. (584 mm)	028827	028863
28 in. (711 mm)	24 in. (610 mm)	028828	028864
29 in. (737 mm)	25 in. (635 mm)	028829	028865
30 in. (762 mm)	26 in. (660 mm)	028830	028866
31 in. (787 mm)	27 in. (686 mm)	028831	028867
32 in. (813 mm)	28 in. (711 mm)	028832	028868
33 in. (838 mm)	29 in. (737 mm)	028833	028869
34 in. (864 mm)	30 in. (762 mm)	028834	028870
35 in. (889 mm)	31 in. (787 mm)	028835	028871
36 in. (914 mm)	32 in. (813 mm)	028836	028872
37 in. (940 mm)	33 in. (838 mm)	028837	028873

Latch Knob for Small Hooped Cable RGO

US

Product Features

- For small Hooped Cable RGOs only
- Warranty: 12 months



Ordering information

DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
Standard	3/16 x 5/8 x 5 in. (4.8 x 15.9 x 127 mm)	028062	028060
XLong	3/16 x 5/8 x 9 in. (4.8 x 15.9 x 228.6 mm)	028066	028064

Push Button for Hooped Cable RGO

US

Product Features

- For medium and large Hooped Cable RGOs
- Push button flexion lock release
- Two step coupling plate to assist standing
- Warranty: 12 months



Ordering information


DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
Medium Standard	1/4 x 3/4 x 6 in. (6.4 x 19.1 x 152.4 mm)	028082	028080
Medium XLong	1/4 x 3/4 x 10 in. (6.4 x 19.1 x 254 mm)	028086	028084
Large Standard	5/16 x 7/8 x 8 in. (7.9 x 22.2 x 203.2 mm)	028092	028090
Large XLong	5/16 x 7/8 x 12 in. (7.9 x 22.2 x 304.8 mm)	028096	028094

RGO II System for Hooped Cable RGO with Abduction

US

Product Features

- Available for medium and large RGO sizes only
- Abduction joint with ring lock release
- Push button flexion lock release
- Two step coupling plate to assist standing
- Automatic relocking with internal spring
- Long, heavy duty lower hip joint bar, compatible with existing upper bars
- Warranty: 12 months



Ordering information

DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
Medium	7/16 x 3/4 x 7 in. (11.1 x 19.1 x 177.8 mm)	026750	026752
Large	1/2 x 3/4 x 10 in. (12.7 x 19.1 x 254 mm)	026756	026754

.....

Small Horizontal Cable RGO

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Small Standard, includes 3/16 x 5/8 x 5 in. (4.8 x 15.9 x 127 mm) Lower Bars	028728
Small XLong includes 3/16 x 5/8 x 9 in. (4.8 x 15.9 x 228.6 mm) Lower Bars	028728XL

SPARE PARTS

- A 027946 - Joint Screw, 1/4-28 x 11/16 in., 2 Each
- B 027979 - Thrust Bearing
- C 027987 - Thrust Race
- D **Lower Bar Assembly, Standard or XLong Bars, 1 Pair Included**
 - Standard, Right - 026880
 - Standard, Left - 026881
 - XLong, Right - 026884
 - XLong, Left - 026885
- E **Coupling Plate, 1 Pair Included**
 - 026669 - Coupling Plate, Right
 - 026668 - Coupling Plate, Left
- F 026708 - Spacer
- G 027961 - Set Screw, 8-32 x 3/8 in., 2 Each
- H 027904 - Upper Hip Joint Bar
- I 027953 - Joint Nut
- J 026550 - Rocker Arm Base, with Screw
- K **Rocker Arm, 1 Pair Included in Assembly**
 - 028580 - Rocker Arm, Right
 - 028581 - Rocker Arm, Left
- L 027045 - Thrust Bearing
- M 026450 - Joint Screw
- N 026447 - Swivel Nut
- O 026452 - Swivel Coupling
- P 880917 - Screw, 8-32 x 3/8 in., 4 Each
- Q 885105 - Jam Nut, 8-32, Small Pattern - 1/4 in. Hex
- R **026732 - Rod End Assembly, 2 Each**
 - S. 880738 - Rod End, Female
 - T. 880784 - Rod End, Male
 - U. 43005 - Washer
 - V. 880212 - Screw, 10-32 x 5/8 in., 4 Each
 - W. 885133 - Jam Nut, 10-32 - 5/16 in. Hex
 - X. 803300 - Extension Stop, 2 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml) Tube
- 028575 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)
- 886887 - Rivets



R. Turnbuckle Assembly



Medium Horizontal Cable RGO without Abduction

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Medium XLong, includes 7/16 x 3/4 x 10 in. (11.1 x 19.1 x 254 mm), Lower Bars	028722XL

SPARE PARTS


- A. 026534 - Joint Screw, 3/8 x 13/16 in., 2 Each
- B. 027045 - Thrust Bearing
- C. 027052 - Thrust Race
- D. **Lower Bar Assembly, Standard or XLong Bars, 1 Pair Included**
 - XLong, Right - 026892
 - XLong, Left - 026893
- E. **Coupling Plate, 1 Pair Included**
 - 026184 - 2 Position Coupling Plate, Right
 - 026186 - 2 Position Coupling Plate, Left
- F. 028621 - Spacer
- G. 027961 - Set Screw, 8-32 x 3/8 in., 2 Each
- H. 026955 - Upper Hip Joint Bar
- I. 026542 - Joint Nut
- J. 026550 - Rocker Arm Base, with Screw
- K. **Rocker Arm, 1 Pair Included in Assembly**
 - 028578 - Rocker Arm, Right
 - 028579 - Rocker Arm, Left
- L. 027045 - Thrust Bearing
- M. 026450 - Joint Screw
- N. 026447 - Swivel Nut
- O. 026452 - Swivel Coupling
- P. 880917 - Screw, 8-32 x 3/8 in., 4 Each
- Q. 885105 - Jam Nut, 8-32, Small Pattern - 1/4 in. Hex
- R. **026732 - Rod End Assembly, 2 Each**
 - S. 880738 - Rod End, Female
 - T. 880784 - Rod End, Male
 - U. 43005 - Washer
 - V. 880212 - Screw, 10-32 x 5/8 in., 4 Each
 - W. 885133 - Nut, 10-32 - 5/16 in. Hex
 - X. 803300 - Extension Stop, 2 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml.) Tube
- 028575 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)
- 886887 - Rivets



R. Turnbuckle Assembly



Medium Horizontal Cable RGO with Abduction

US

Ordering information

DESCRIPTION	PRODUCT NUMBER
Medium, includes 7/16 x 3/4 x 8-1/4 (11.1 x 19.1 x 209.6 mm) RGO II System Lower	028752

SPARE PARTS


- A. 026534 - Joint Screw, 3/8 x 13/16 in., 2 Each
- B. 027045 - Thrust Bearing
- C. 027052 - Thrust Race
- D. **Lower Bar Assembly, RGO II System with Abduction, 1 Pair Included**
 - 026156 - RGO II System, Right
 - 026154 - RGO II System, Left
- E. **Coupling Plate, 1 Pair Included**
 - 026184 - 2 Position Coupling Plate, Right
 - 026186 - 2 Position Coupling Plate, Left
- F. 028621 - Spacer
- G. 027961 - Set Screw, 8-32 x 3/8 in., 2 Each
- H. 026955 - Upper Hip Joint Bar
- I. 026542 - Joint Nut
- J. 026550 - Rocker Arm Base, with Screw
- K. **Rocker Arm, 1 Pair Included in Assembly**
 - 028578 - Rocker Arm, Right
 - 028579 - Rocker Arm, Left
- L. 027045 - Thrust Bearing
- M. 026450 - Joint Screw
- N. 026447 - Swivel Nut
- O. 026452 - Swivel Coupling
- P. 880917 - Screw, 8-32 x 3/8 in., 4 Each
- Q. 885105 - Nut, 8-32, Small Pattern - 1/4 in. Hex
- R. **026732 - Rod End Assembly, 2 Each**
 - S. 880738 - Rod End, Female
 - T. 880784 - Rod End, Male
 - U. 43005 - Washer
 - V. 880212 - Screw, 10-32 x 5/8 in., 4 Each
 - W. 885133 - Nut, 10-32 - 5/16 in. Hex
 - X. 803300 - Extension Stop, 2 Each

Also Included But Not Pictured:

889002 - Loctite #242, .02 oz. (0.6 ml.) Tube

028575 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)

886887 - Rivets



R. Turnbuckle Assembly



Large Horizontal Cable RGO without Abduction

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Large XLong, includes s 1/2 x 7/8 x 12 in. (12.7 x 22.2 x 304.8 mm) Lower Bars	028726XL

SPARE PARTS


- A. 027771 - Joint Screw, 7/16 x .843 in., 2 Each
- B. 027797 - Thrust Bearing
- C. 027805 - Thrust Race
- D. **Lower Bar Assembly, Standard or XLong Bars, 1 Pair Included**
 - 026900- XLong, Right
 - 026901- XLong, Left
- E. **Coupling Plate, 1 Pair Included**
 - 026706 - 2 Position Coupling Plate, Right
 - 026710 - 2 Position Coupling Plate, Left
- F. 028647 - Spacer
- G. 027896 - Screw, 10-32 x 3/8 in., 2 Each
- H. 027730 - Upper Hip Joint Bar
- I. 028399 - Joint Nut
- J. 026550 - Rocker Arm Base, with Screw
- K. **Rocker Arm, 1 Pair Included in Assembly**
 - 028578 - Rocker Arm, Right
 - 028579 - Rocker Arm, Left
- L. 027045 - Thrust Bearing
- M. 026450 - Joint Screw
- N. 026447 - Swivel Nut
- O. 026452 - Swivel Coupling
- P. 880917 - Screw, 8-32 x 3/8 in., 4 Each
- Q. 885105 - Nut, 8-32, Small Pattern - 1/4 in. Hex
- R. **026732 - Rod End Assembly, 2 Each**
 - S. 880738 - Rod End, Female
 - T. 880784 - Rod End, Male
 - U. 43005 - Washer
 - V. 880212 - Screw, 10-32 x 5/8 in., 4 Each
 - W. 885133 - Nut, 10-32 - 5/16 in. Hex
 - X. 803300 - Extension Stop, 2 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml.) Tube
- 028575 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)
- 886887 - Rivets



R. Turnbuckle Assembly



Large Horizontal Cable RGO with Abduction

US

Ordering information


DESCRIPTION	PRODUCT NUMBER
Large 1/2 x 3/4 x 11-3/8 in. (12.7 x 19.1 x 289 mm) RGO II System Lower Bars	028754

SPARE PARTS

- A. 027771 - Joint Screw, 7/16 x .843 in., 2 Each
- B. 027797 - Thrust Bearing
- C. 027805 - Thrust Race
- D. **Lower Bar Assembly, RGO II System with Abduction, 1 Pair Included**
 - 026196 - RGO II System, Right
 - 026198 - RGO II System, Left
- E. **Coupling Plate, 1 Pair Included**
 - 026706 - 2 Position Coupling Plate, Right
 - 026710 - 2 Position Coupling Plate, Left
- F. 028647 - Spacer
- G. 027896 - Screw, 10-32 x 3/8 in., 2 Each
- H. 027730 - Upper Hip Joint Bar
- I. 028399 - Joint Nut
- J. 026550 - Rocker Arm Base, with Screw
- K. **Rocker Arm, 1 Pair Included in Assembly**
 - 028578 - Rocker Arm, Right
 - 028579 - Rocker Arm, Left
- L. 027045 - Thrust Bearing
- M. 026450 - Joint Screw
- N. 026447 - Swivel Nut
- O. 026452 - Swivel Coupling
- P. 880917 - Screw, 8-32 x 3/8 in., 4 Each
- Q. 885105 - Nut, 8-32, Small Pattern - 1/4 in. Hex
- R. **026732 - Rod End Assembly, 2 Each**
 - S. 880738 - Rod End, Female
 - T. 880784 - Rod End, Male
 - U. 43005 - Washer
 - V. 880212 - Screw, 10-32 x 5/8 in., 4 Each
 - W. 885133 - Nut, 10-32 - 5/16 in. Hex
 - X. 803300 - Extension Stop, 2 Each

Also Included But Not Pictured:

- 889002 - Loctite #242, .02 oz. (0.6 ml.) Tube
- 028575 - Outriggers, 2 Each (For Plastic Girdle Fabrication, Refer to RGO Manual)
- 886887 - Rivets




R. Turnbuckle Assembly

Horizontal RGO Cables

US

Product Features

When fabricating the Horizontal Cable System, it is impossible for Orthotists to predict the necessary length of cables. Fabricators must complete the basic structure then measure to determine which length of cables are needed. The cable starter kit provides orthotists with a selection of cables on hand helping to avoid lengthy assembly processes.



Ordering information

TYPE "A" CABLE (NOT ADJUSTABLE)

SIZE	PRODUCT NUMBER
3-3/4 - 4-1/2 in. (9.5 - 11.4 cm)	027200
4-1/4 - 5 in. (10.8 - 12.7 cm)	027202
4-3/4 - 5-1/2 in. (12.1 - 14 cm)	027204
5-1/4 - 6 in. (13.3 - 15.2 cm)	027206
5-3/4 - 6-1/2 in. (14.6 - 16.5 cm)	027208
6-1/4 - 7 in. (15.9 - 17.8 cm)	027210
6-3/4 - 7-1/2 in. (17.1 - 19.1 cm)	027212
7-1/4 - 8 in. (18.4 - 20.3 cm)	027214

TYPE "B" CABLE (ADJUSTABLE)

SIZE	PRODUCT NUMBER
8 - 9 in. (20.3 - 22.9 cm)	027220
9 - 10 in. (22.9 - 25.4 cm)	027222
10 - 11 in. (25.4 - 28 cm)	027224
11 - 12 in. (28 - 30.5 cm)	027226
12 - 13 in. (30.5 - 33 cm)	027228
13 - 14 in. (33 - 35.6 cm)	027230
14 - 15 in. (35.6 - 38 cm)	027232
15 - 16 in. (38 - 40.6 cm)	027234
16 - 17 in. (40.6 - 43.2 cm)	027236
Set Screw, 10-32 x 1-1/2 in. For Type "B" Cables	880209

CABLE SET

SIZE	PRODUCT NUMBER
Cable Starter Set, Includes 2 Each of Type "B" Cables Listed Below	027240

Latch Knob for Small Horizontal Cable and Rocker Bar RGOs

US

Product Features

- For small RGOs, Horizontal Cable or Rocker Bar designs only
- Extra long lower bar assembly by special order only
- Warranty: 12 months



Ordering information

DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
XLong	3/16 x 5/8 x 9 in. (4.8 x 15.9 x 228.6 mm)	028706XL	028702XL



Push Button for Horizontal Cable and Rocker Bar RGOs

US

Product Features

- For medium and large Horizontal Cable and Rocker Bar RGOs
- Push button flexion lock release
- Two step coupling plate to assist standing
- Extra long lower bar assembly by special order only
- Warranty: 12 months



Ordering information

DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
Medium XLong	1/4 x 3/4 x 10 in. (6.4 x 19.1 x 254 mm)	028710XL	028708XL
Large XLong	5/16 x 7/8 x 12 in. (7.9 x 22.2 x 304.8 mm)	028714XL	028712XL

RGO II System for Horizontal Cable and Rocker Bar RGOs

US

Product Features

- Available for medium and large Horizontal Cable and Rocker Bar RGOs
- Abduction joint with ring lock release
- Push button flexion lock release
- Two step coupling plate to assist standing
- Automatic relocking with internal spring
- Long, heavy duty lower hip joint bar, compatible with existing upper bars
- Warranty: 12 months



Ordering information


DESCRIPTION	LOWER BAR	PROD. NO. LEFT	PROD. NO. RIGHT
Medium	7/16 x 3/4 x 7 in. (11.1 x 19.1 x 177.8 mm)	026150	026152
Large	1/2 x 3/4 x 10 in. (12.7 x 19.1 x 254 mm)	026192	026190

Specialty Connectors for ISOCENTRIC® RGO

US

Product Features

- Available in stem size 3/16 in. (4.8 mm) or 1/4 in. (6.4 mm)
- Durable design and easy to use
- **A)** Standard connector (104 and 104X), dialed 4 - 6 cm
- **B)** Connector with extensions for large hip flexion angle (106C), dialed 6 - 8 cm
- **C)** "Easy adjust" connector, allows for fuss-free adjustment of torso angle at hips (106)



Ordering information

DESCRIPTION	PRODUCT NUMBER
Standard Connector (A)	104
Standard Connector, Heavy Duty (A)	104X
Standard Connector, Female	62
Standard Connector, Male	63
Extended Connector, Top	63T
Extended Connector, Bottom	63B
Easy Adjust Connector Assembly (C)	106
Easy Adjust Connector Assembly, Custom (B)	106C
Easy Adjust Tie Rod End, Male	EA1
Easy Adjust, Turnbuckle	EA2
Easy Adjust, Set Screw	EA3
Easy Adjust Tie Rod End, Female	EA4
Standard Connector Extender	EA5

Cable Track Kit

US

Product Features

Use when fabricating RGO metal pelvic bands. Cut the cable track to size, per fabricating instructions in RGO manual. When permanently attached to RGO pelvic band the cable track supports the cables and prevents them from rubbing the patient.



Ordering information

SIZE	PRODUCT NUMBER
Cable Track Kit, For Metal Pelvic Bands	028730
SPARE PARTS	
026650 - Cable Track Molding, Copolymer (1 Each)	
885135 - Screw, Socket-Flat Head, 8-32 x 3/8 in. (4 Each)	

Aluminum Pelvic Bands

US

Product Features

- Full size patterns for cutting bands and contouring bands available on request
- Manual, measurement sheets and information on central fabrication services available on request



Ordering information

PELVIC CIRCUMFERENCE	BAND SIZE	PRODUCT NUMBER
3/16 x 1-3/4 x 14-1/2 in. (4.8 x 44.5 x 368.3 mm)	19 – 21 in. (483 - 533 mm)	027498
3/16 x 1-3/4 x 15-1/2 in. (4.8 x 44.5 x 393.7 mm)	21 – 23 in. (533 - 584 mm)	027500
3/16 x 1-3/4 x 16-3/4 in. (4.8 x 44.5 x 425.5 mm)	23 – 25-1/2 in. (584 - 648 mm)	027502
3/16 x 1-3/4 x 18 in. (4.8 x 44.5 x 457.2 mm)	25-1/2 – 28 in. (648 - 711 mm)	027504
1/4 x 1-3/4 x 19-1/4 in. (6.4 x 44.5 x 489.0 mm)	28 – 30-1/2 in. (711 - 775 mm)	027506
1/4 x 1-3/4 x 21 in. (6.4 x 44.5 x 533.4 mm)	31 – 34 in. (787 - 864 mm)	027508
1/4 x 1-3/4 x 22-1/2 in. (6.4 x 44.5 x 571.5 mm)	34 – 37 in. (864 - 940 mm)	027510
1/4 x 1-3/4 x 24 in. (6.4 x 44.5 x 609.6 mm)	37 – 40 in. (940 - 1016 mm)	027512
1/4 x 1-3/4 x 26 in. (6.4 x 44.5 x 660.4 mm)	40 – 43 in. (1016 - 1092 mm)	027514

Lower Abduction Bar Assembly

US



Ordering information

SIZE	PRODUCT NUMBER RIGHT	PRODUCT NUMBER LEFT
Medium	026156	026154
Large	026196	026198

Lower Bar Assembly without Abduction

US



Latch Knob Style For Small Sizes



Push Button Style For Medium and Large Sizes

Ordering information

SIZE	PRODUCT NUMBER RIGHT	PRODUCT NUMBER LEFT
Small, Latch Knob	026880	026881
Small, Latch Knob X-long	026884	026885
Medium, Push Button X-long	026892	026893
Large, Push Button X-long	026900	026901

Knee Joints

US

Product Features

- Constructed to be very stiff to lateral forces
- Heavy duty joints can be ordered with custom angle of pre-flexion (extra charge applies)
- Patients with hip and knee contractures can special order joints

Ordering information

DESCRIPTION	PROD. NO. PAIR	PROD. NO. LEFT	PROD. NO. RIGHT
Standard 1/4 x 3/4 in. (6.4 x 19.1 mm)	023802	208SL	208SR
Extra Duty, Pair 1/4 x 1 in. (6.4 x 25.4 mm)	208X	208XL	208XR
Extra Extra Duty, Pair 3/8 x 1 in. (9.5 x 25.4 mm)	208XX	208XXL	208XXR




Pelvic Bands

US

Product Features

- Specialty pelvic bands for other than RGO usage; great for HKAFOs, Legg-Perthes, C.P. and other heavy uses
- Custom made to cast and measurements
- Heavy duty welded high strength aluminum construction
- Can be ordered “Standard” or “Light Duty”
- Can be ordered with free or lockable hip joints, as well as adjustable hip flexion angle
- Pelvic band can be made regular or light duty



Ordering information

DESCRIPTION	PRODUCT NUMBER
Pelvic Band (Custom)	E4

External Ankle Foot Orthosis (AFOs)

US

Product Features

- As part of a RGO, these AFOs improve balance greatly, allow for easy don/doff, and eliminate footwear problems because they allow a patient to wear street shoes
- Work well as floor reaction AFOs
- Custom made from casts or measurements
- Available in three colors of polypropylene (blue, black, and natural)

Non-articulated



Ordering information

DESCRIPTION	PRODUCT NUMBER
Non-Articulated	211N

Standard In-Shoe Ankle Foot Orthosis (AFO)

US

Product Features

- Custom made from casts or measurements
- Available in three colors of polypropylene (blue, black, and natural)

Ordering information


DESCRIPTION	PRODUCT NUMBER
Standard In-Shoe	211S

Gauges

US

Product Features

- Anatomical measuring tool
- Measures in inches and metric
- Adjustable Pelvic Casting Frame is available by request



Ordering information

DESCRIPTION	PRODUCT NUMBER
A. 17 in. (432 mm) Universal A-P & M-L Gauge	220336
B. Original AP-ML Gauge	220335
C. Casting Fixture Gauge	220355



Spinal Braces

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

Providence Nighttime Brace

US

Product Features

- Over corrective
- Controls opposite forces & both directly lateral and rotational forces
- Polypropylene material
- Special measurement board
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Providence Brace, finished	688500000
Measurement Board incl. Blocks (adult)	688000101
Extra Blocks (adult)	688000102
Pediatrics Blocks	688000104

For order form see page 88



Scoliosis Full Time Brace

US

Product Features

- Full time bracing system for the treatment of idiopathic scoliosis
- This orthosis relies on a series of pressures and void areas to exert corrective forces on the spinal column
- Indications: lumbar curves, thoracic curves, double major curves and excessive lordosis
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Scoliosis Module, Unfinished	688100000
For order form see page 89	

LSO - Single Opening

US

Product Features

- Optimal circumferential control
- High degree of adjustability
- Easy donning and doffing
- Available with posterior, anterior and lateral opening
- Indications: spondylolisthesis, lumbar fusion, degenerative disc disease, laminectomy, sports injuries, compression fractures and spondylitis
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Single Opening, Finished	688600010
Single Opening, Unfinished	688600011

For order form see page 90

LSO - Bi-valve

US

Product Features

- This particular design provides the maximum degree of control with comprehensive adjustability due to volume changes
- Easy donning and doffing
- Often used in post-operative situations
- Indications: post-operative lumbar surgery, herniated disc, low back pain, compression fracture, spondylolisthesis, degenerative disc disease
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Bi-valve, Finished	688600020
Bi-valve, Unfinished	688600021

For order form see page 90

TLSO - Single Opening (A.O.O.)

US

Product Features

- Optimal circumferential control
- High profile controls spinal pathologies in the lumbar and thoracic regions
- Available with posterior/- anterior/- or lateral opening
- Indications: spinal stabilization, osteoporosis, compression fracture, thoracic immobilization, degenerative disc disease and spondylitis
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Single Opening, Finished	688600012
Single Opening, Unfinished	688600013

For order form see page 90



TLSO - Bi-valve

US

Product Feature

- Maximum adjustability while maintaining an intimate and comfortable fit
- It provides uniform rotational control while stabilizing the pelvis
- Indications: post-operative fusion, kyphosis, compression fracture, osteoporosis, burst Fx, spinal stenosis, spinal muscular atrophy and laminectomy
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Bi-valve, Finished	688600022
Bi-valve, Unfinished	688600023

For order form see page 90

Flex Foam™ I

US

Product Features

- The Flex Foam™ I offers a softer inner layer of foam combined with a dense outer foam for the relief of bony prominences
- Rigid external frame. The frame is detachable and can be strategically adjusted as indicated by patient needs
- Available with single- and bivalve openings
- Available in LSO and TLSO
- Indications: cerebral palsy, osteoporosis, myelomeningocele, wheelchair seating support, scoliotic deformities, muscular dystrophy, cancer patients, geriatrics, lumbar fusion, degenerative disc disease, spondylolisthesis, spinal stenosis
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Single Opening, Finished	688600030
Single Opening, Unfinished	688600031
Bivalve Opening, Finished	688600032
Bivalve Opening, Unfinished	688600033

For order form see page 91

Flex Foam II

US

Product Features

- The Flex Foam™ II offers a rigid interior frame sandwiched between foam layers, allowing the wearer a greater degree of comfort while providing the appropriate rigidity for support
- Available with single- or bivalve openings
- Available in LSO and TLSO
- Indications: osteoporosis, kyphotic control, spinal stenosis and post-operative spinal surgery, herniated disc, compression fracture, spondylolysis
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Single Opening, Finished	688600040
Single Opening, Unfinished	688600041
Bivalve Opening, Finished	688600042
Bivalve Opening, Unfinished	688600043

For order form see page 91

Flex Foam III

US

Product Features

- The Flex Foam™ III is a soft but firm construction made of "laminated foam" for less rigid support
- Available in LSO and TLSO
- Indications: cerebral palsy, osteoporosis, myelomeningocele, wheelchair seating support, scoliotic deformities, muscular dystrophy, cancer patients, geriatrics
- For available transfer paper please see: <https://spinaltech.com/transfer-patterns>



Ordering information

DESCRIPTION	PRODUCT NUMBER
Single Opening, Finished	688600050
Single Opening, Unfinished	688600051
Bivalve Opening, Finished	688600052
Bivalve Opening, Unfinished	688600053

For order form see page 91

Pectus Brace

US

Product Features

- Provides anterior/posterior corrective compression over affected site
- Custom to measurements
- Compression adjusted with ratchet mechanism
- Closed cell foam pads for ease of cleaning
- Indications: pectus carinatum



Ordering information

DESCRIPTION	PRODUCT NUMBER
Pectus Brace	688600060
For order form see page 91	

Fillauer

E-mail: order@fillauer.com
Phone: +46 8 505 332 00

Order form - Providence

Customer: Fillauer's PO nr:

Ordered by: Patient number:

Email: Date:

Shipping instructions from Spinal to Fillauer: UPS express saver on Fillauer's account


Anterior opening, with tongue:		Yes ☐	No ☐
Style:	Finished ☒	Straps ☒	
Patient is a current or previous brace wearer:	Patient number stated with last Fillauer order: _____		
Male: ☐	Female: ☐	Age: _____	Hgt (cm): _____ Wgt (kg): _____
X-rays attached by email:		Yes ☒	
Transfer paper choice 1: _____			
Transfer paper choice 2: _____			
Additional information: _____			

Measurements cm:		Curve Analysis:	
Cir		Thoracic Curve: L or R ☐	
Axilla	☐	Thoracic Apex: T _____	
Xyphoid	☐	Lumbar Curve: L or R ☐	
Waist	☐	Lumbar Apex: L _____	
ASIS	☐	Thoracolumbar Curve: L or R ☐	
Troch	☐	Thoracolumbar Apex: L1 or T12 ☐	
Xyphoid Waist/Row 0 Pubis Length and A/P measurements taken with patient on board and pressure pads in place.			

Pad Placement with Board

	Column (down)	Row (across)
Stabilizing pads:		
Axilla	_____	_____
Trochanter Pad	_____	_____
	Column (down)	Row (across)
Pressure pads:		Distance
Lumbar Pad	_____	_____ mm
Thoracic Pad	_____	_____ mm

☐ Pediatric blocks were used.



Spinal Use Only

M	P	F	S	Inspected by	Mold #
---	---	---	---	--------------	--------



Order form Scoliosis/ Kyphosis Brace

Fillauer

Kung Hans väg 2
192 68 Sollentuna
Sweden

Phone: +46 8 505 332 00
E-mail: order@fillauer.com

Customer:

Fillauer's PO nr.:

Patient number:

Date & Ordered by:

Phone/email:

Shipping instruction from Spinal to Fillauer:
UPS Express Saver on Fillauer's account

Opening: Anterior ☐ Posterior ☐		Milwaukee orthosis: ☐	Neck: _____
Style: Unfinished ☒ Straps (for brace) ☒		Milwaukee anterior upright (with posterior reinforcement and straps)	ML: _____ AP: _____
Pads: Soft ☐ Firm ☐		Low Profile TLSO orthosis (Kyphosis): (with adjustable anterior sternal shield, posterior reinforcements and straps) ☐	
Lumbar _____ Troch _____ Thoracic _____ Axilla _____		Instructions:	
Male/ Female ☐ ☐ Age: _____ Height(cm): _____ Weight(kg): _____			
Lordosis: ☐ 5° ☐ 15° ☐ Other: _____			
Diagnosis:			

Shape/ Reliefs:
(e; bony prom, etc)

Office use only

M	P
F	S
Inspected by:	

Measurements taken: Standing ☐ Supine ☐

Abdominal reliefs:
(if required - choose Pend or Full followed by size)

Pend ☐
Full ☐
Small ☐
Med ☐
Lg ☐
X-Lg ☐

Topo ☐

Labels on diagram: Cir., ML, Breast Cir., Cup Size, Chest, Xyphoid, Waist, ASIS, Topo, 1" Above Trochanter, Sternal Notch, Nipple Line, Xyphoid II, Apex, Axilla, Occipital, Superior Scapula, Inferior Scapula, Coccyx, Troc, G. Fold or to Firm Seat, Symphysis Pubis, ASIS, Xyphoid, Pend, Full, Small, Med, Lg, X-Lg, Topo.

Fillauer

Order form LSO/TLSO, Single or Bivalve Opening

Fillauer

Kung Hans väg 2
192 68 Sollentuna
Sweden

Phone: +46 8 505 332 00
E-mail: order@fillauer.com

Customer:

Fillauer Europe's PO#:

Patient number:

Date & Ordered by:

Phone/email:

Shipping instruction from Spinal to Fillauer Europe:
UPS Express Saver on Fillauer's account

1.	☐ Single opening	☐ Anterior overlap	Anterior w/ tongue	☐ Tongue "Polyethen" ☐ Tongue "Firm foam"	☐ Posterior opening	Lateral opening ☐ L ☐ R
	☐ Bi-valve opening	Smooth overlap	☐ Ant. over Post. ☐ Ant. into Post.	☐ Step overlap	☐ Open with tongue	
2.	☐ High Profile (TLSO)	☐ Finished	☐ Lined	☐ Straps		
	☐ Low Profile (LSO)	☐ Unfinished	☐ Unlined			
3.	☐ Male	Age:	Height(cm):	Weight(kg):		
	☐ Female					

Diagnosis:

Shape/ Reliefs:
(boney, prom. etc.)

Instructions:

4.	Lordosis: ☐ 15° ☐ Other:	Office use only		
		<table border="1"> <tr> <td>M</td> <td>P</td> </tr> <tr> <td>F</td> <td>S</td> </tr> </table>	M	P
M	P			
F	S			

Measurements taken (in cm): Supine ☐ Standing ☐

Abdominal reliefs:
(if required - choose Pend or Full followed by size)

Pend ☐
Full ☐
Small ☐
Med ☐
Lg ☐
X-Lg ☐

Topo ☐
Abdominal Binder ☐

Office Use Only				
M	P	F	S	Inspected by
				Mold #

Fillauer



Order form Flex Foam

Fillauer.

Kung Hans Väg 2
19268 Sollentuna
Sweden

Phone: +46 8 505 332 00
E-mail: order@fillauer.com

Fillauer Europe's PO#:

Patient number:

Date & Ordered by:

Phone/email:

Shipping instruction from Spinal to Fillauer Europe :
UPS Express Saver on Fillauer Europe's account

1.	☐ External frame	☐ Internal frame "Frame between foam layers"	☐ No frame "Laminated foam"			
2.	☐ Single opening	☐ Anterior w/ tongue	☐ Tongue "Polyethen" ☐ Tongue "Firm foam"	☐ Posterior opening	☐ Anterior overlap	Lateral opening ☐ L ☐ R
	☐ Bi-valve opening	☐ Smooth Overlap	☐ Anterior over Post. Anterior into Post.	☐ Step overlap	☐ Open with tongues	
3.	☐ High Profile (TLSO) ☐ Low Profile (LSO)		☐ Finished ☐ Unfinished	Plastic thickness Standard "3mm" Other: _____		Liner thickness Standard "6mm" Other: _____
	Male/ Female ☐ ☐	Age: _____	Height(cm): _____	Weight(kg): _____		
Lordosis: ☐ 15° ☐ Other: _____ Measurements taken (in cm): ☐ Standing ☐ Supine						
Diagnosis: _____						
4.	Instructions:			Shape/ Reliefs: (bone, prom. etc.)		

Abdominal reliefs:
(if required - choose Pend or Full followed by size)

Pend ☐
Full ☐
Small ☐
Med ☐
Lg ☐
X-Lg ☐

Abdominal Binder ☐

Office Use Only

M	P	F	S	Inspected by	Mold #
---	---	---	---	--------------	--------

Fillauer.

Order form Pectus Brace

Kung Hans väg 2
192 68 Sollentuna
Sweden

Phone: +46 8 505 332 00

E-mail: fillauer@fillauer.com

Fillauer Europe's PO#:

Customer:

Patient name:

Ordered by:

Phone:

E-mail:


Shipping instruction:

UPS Express Saver on Fillauer Europe's account

Straps:

Ratchet

Special instructions:



A - M/L at Level of Deformity

B - Clearly Mark Deformity

C - Circumference at Level of Deformity

D - Distance from Apex of Deformity to Sternal Notch

Fillauer.

Scoliosis Module

US

Product Features

- The Scoliosis Module System is a full time bracing system for the treatment of idiopathic scoliosis. It uses a system of pressure and void areas to exert corrective forces upon the spinal column.
- 22 different sizes
- Available in 5° or 15° lordosis



Ordering information

XYPHOID	WAIST	HIP	ASIS	PROD. NO.
48 cm	41 cm	53 cm	13 cm	688215001
51 cm	41 cm	58 cm	15 cm	688215002
53 cm	43 cm	61 cm	15 cm	688215003
56 cm	46 cm	66 cm	18 cm	688215004
56 cm	46 cm	71 cm	20 cm	688215005
58 cm	46 cm	71 cm	20 cm	688215006
61 cm	46 cm	76 cm	20 cm	688215007
64 cm	51 cm	79 cm	23 cm	688215008
66 cm	56 cm	79 cm	23 cm	688215009
66 cm	51 cm	81 cm	23 cm	688215010
66 cm	58 cm	81 cm	23 cm	688215011
71 cm	56 cm	84 cm	23 cm	688215012
69 cm	53 cm	86 cm	25 cm	688215013
69 cm	58 cm	86 cm	23 cm	688215014
74 cm	58 cm	89 cm	25 cm	688215015
74 cm	64 cm	89 cm	23 cm	688215016
71 cm	61 cm	91 cm	25 cm	688215017
76 cm	69 cm	91 cm	25 cm	688215018
76 cm	61 cm	94 cm	28 cm	688215019
84 cm	71 cm	94 cm	25 cm	688215020
79 cm	69 cm	97 cm	28 cm	688215021
86 cm	76 cm	97 cm	28 cm	688215022

A.O.O. Spinal Brace

US

Product Features

- The Anterior Overlap Opening (AOO) Module is a supportive bracing system which comes in a variety of common sizes. As it does not require custom sizing and manufacturing, it can be purchased on demand and without the cost or lead time of custom pieces
- Available in 16 different sizes
- Lined and unlined
- 15° lordosis
- Warranty: 12 months



Ordering information

FEMALE SIZES/CM

	XYPHOID	WAIST	HIP	PROD. NO. LINED	PROD. NO. UNLINED
F33	74	66	84	688015033	688115033
F35	79	71	89	688015035	688115035
F37	84	76	94	688015037	688115037
F39	89	81	99	688015039	688115039
F41	94	86	104	688015041	688115041
F43	99	91	109	688015043	688115043
F45	104	97	114	688015045	688115045
F47	109	102	119	688015047	688115047

MALE SIZES/CM

	XYPHOID	WAIST	HIP	PROD. NO. LINED	PROD. NO. UNLINED
M33	81	71	84	688015133	688115133
M35	86	76	89	688015135	688115135
M37	91	81	94	688015137	688115137
M39	97	86	99	688015139	688115139
M41	102	91	104	688015141	688115141
M43	107	97	109	688015143	688115143
M45	112	102	114	688015145	688115145
M47	117	107	119	688015147	688115147

SPRING STEEL STRAP KIT

DESCRIPTION	PRODUCT NUMBER
15 cm	688000106
17,5 cm	688000107
20 cm	688000108
22,5 cm	688000109





Basic Strap

US



Product Features

- Basic strap with velcro
- Sold individually
- Spare part Providence

Ordering information

SIZE	PRODUCT NUMBER
25 x 254 mm	688000137
38 x 305 mm	688000130
50 x 305 mm	688000131
38 x 355 mm	688000132
50 x 355 mm	688000014

Anterior Stay

US



Product Features

- Sold individually

Ordering information

SIZE	PRODUCT NUMBER
16 cm	688000016
20 cm	688000020
26 cm	688000026

Welding Rod

US



Ordering information

SIZE	PRODUCT NUMBER
Diameter 3 mm, Length 100 cm	688000103
Diameter 5 mm, Length 100 cm	688000105

Spring Steel Strap

US

Product Features

- Sold in pairs



Ordering information

SIZE	PRODUCT NUMBER
2.54 x 15 cm	688000006
2.54 x 20 cm	688000008
2.54 x 25.4 cm	688000010
2.54 x 30.5 cm	688000012
2.54 x 35.5 cm	688000014

Chafes with Triangular Body Clasp

US

Product Features

- Sold individually
- Spare parts for Providence Braces
- Works with Spring Steel Strap
- Works with Basic Strap



Ordering information

SIZE	PRODUCT NUMBER
25 mm	688000136
38 mm	688000134
50 mm	688000135



Insoles

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

TFC - Thermoformable Composite Insole

SE

Product Features

- Flexible insoles are 0.06 in. (1.5 mm) thick and are better for shaping double curved, for example around the heel
- Rigid insoles are 0.1 in. (2.5 mm) thick and made to be used as a rolling insole for diabetic patients. Made for shaping single curves and is harder to shape around the heel
- Sold in pairs



Ordering information

TFC INSOLE FLEXIBLE, COMFIL®

SIZE	EUROPE	US WOMEN	US MEN	PRODUCT NUMBER
35 - 37		5 - 6.5	3.5 - 5	700030121
38 - 40		7.5 - 9.5	6 - 7.5	700030122
41 - 43		9.5 - 10.5	8 - 9	700030123
44 - 46		12 - 13.5	10.5 - 12	700030124

TFC INSOLE RIGID, COMFIL®

SIZE	EUROPE	US WOMEN	US MEN	PRODUCT NUMBER
35 - 37		5 - 6.5	3.5 - 5	700030141
38 - 40		7.5 - 9.5	6 - 7.5	700030142
41 - 43		9.5 - 10.	8 - 9	700030143
44 - 46		12 - 13.5	10.5 - 12	700030144



PEVA

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

PEVA Shoe - Open Toe Boots - Firm Heel Cap

ES

Specifications

Warranty: 12 month

Product Features

- Full leather lining
- Padded and free standing tongue
- Extra strengthening of the heel cap
- Heel inspection window



Ordering information

STRAIGHT MEDIUM LAST

SIZE	LAST LENGTH	PROD. NO. PAIR	PROD. NO. LEFT	PROD. NO. RIGHT
12	2.9 in. (73 mm)	697800012	697800112	697800212
13	3.1 in. (80 mm)	697800013	697800113	697800213
14	3.4 in. (86 mm)	697800014	697800114	697800214
15	3.7 in. (93 mm)	697800015	697800115	697800215
16	3.9 in. (100 mm)	697800016	697800116	697800216
17	4.2 in. (106 mm)	697800017	697800117	697800217
18	4.4 in. (113 mm)	697800018	697800018	697800218
19	4.7 in. (120 mm)	697800019	697800119	697800219
20	5.0 in. (126 mm)	697800020	697800120	697800220
21	5.2 in. (133 mm)	697800021	697800121	697800221
22	5.5 in. (140 mm)	697800022	697800122	697800222
23	5.8 in. (147 mm)	697800023	697800123	697800223
24	6.0 in. (153 mm)	697800024	697800124	697800224
25	6.2 in. (158 mm)	697800025	697800125	697800225
26	6.5 in. (165 mm)	697800026	697800126	697800226
27	6.8 in. (173 mm)	697800027	697800127	697800227

**STRAIGHT NARROW LAST**

SIZE	LAST LENGTH	PROD. NO. PAIR	PROD. NO. LEFT	PROD. NO. RIGHT
12	2.9 in. (73 mm)	697805012	697805112	697805212
13	3.1 in. (80 mm)	697805013	697805113	697805213
14	3.4 in. (86 mm)	697805014	697805114	697805214
15	3.7 in. (93 mm)	697805015	697805115	697805215
16	3.9 in. (100 mm)	697805016	697805116	697805216
17	4.2 in. (106 mm)	697805017	697805117	697805217
18	4.4 in. (113 mm)	697805018	697805118	697805218
19	4.7 in. (120 mm)	697805019	697805119	697805219
20	5.0 in. (126 mm)	697805020	697805120	697805220
21	5.2 in. (133 mm)	697805021	697805121	697805221
22	5.5 in. (140 mm)	697805022	697805122	697805222
23	5.8 in. (147 mm)	697805023	697805123	697805223
24	6.0 in. (153 mm)	697805024	697805124	697805224
25	6.2 in. (158 mm)	697805025	697805125	697805225
26	6.5 in. (165 mm)	697805026	697805126	697805226
27	6.8 in. (173 mm)	697805027	697805127	697805227

ABDUCTION LAST

SIZE	LAST LENGTH	PROD. NO. PAIR	PROD. NO. LEFT	PROD. NO. RIGHT
12	2.9 in. (73 mm)	697803012	697803112	697803212
13	3.1 in. (80 mm)	697803013	697803113	697803213
14	3.4 in. (86 mm)	697803014	697803114	697803214
15	3.7 in. (93 mm)	697803015	697803115	697803215
16	3.9 in. (100 mm)	697803016	697803116	697803216
17	4.2 in. (106 mm)	697803017	697803117	697803217
18	4.4 in. (113 mm)	697803018	697803118	697803218
19	4.7 in. (120 mm)	697803019	697803119	697803219
20	5.0 in. (126 mm)	697803020	697803120	697803220
21	5.2 in. (133 mm)	697803021	697803121	697803221
22	5.5 in. (140 mm)	697803022	697803122	697803222
23	5.8 in. (147 mm)	697803022	697803123	697803223
24	6.0 in. (153 mm)	697803024	697803124	697803224
25	6.2 in. (158 mm)	697803025	697803125	697803225
26	6.5 in. (165 mm)	697803026	697803126	697803226
27	6.8 in. (173 mm)	697803027	697803127	697803227

Clip-On Bar

US

Specifications

Warranty: 12 month

Product Features

- Adjustable bar
- Plates that clip on to bar
- Available with 0° or 10° dorsal flexion



Ordering information

DESCRIPTION	LENGTH	PRODUCT NUMBER
Clip-On Bar 0° dorsal flexion	6 - 8 in. (152 - 203 mm)	697000020
Clip-On Bar 0° dorsal flexion	8 - 12 in. (203 - 305 mm)	697000021
Clip-On Bar 10° dorsal flexion	6 - 8 in. (152 - 203 mm)	697000022
Clip-On Bar 10° dorsal flexion	8 - 12 in. (203 - 305 mm)	697000023

Denis Browne Splint

US

Specifications

Warranty: 12 month

Product Features

- Adjustable bar
- Holes for splint adapter
- Includes protractor wrench, hex key and screws



Ordering information

SIZE	LENGTH	PRODUCT NUMBER
Small Footplates 3.0 in. (75 mm)	7-11 in. (177 - 279 mm)	012195
Small Footplates 3.0 in. (75 mm)	9-15 in. (228 - 381 mm)	012203
Large Footplates 3.9 in. (100 mm)	7-11 in. (177 - 279 mm)	012196
Large Footplates 3.9 in. (100 mm)	9-15 in. (228 - 381 mm)	012204

Dobbs Bar

US

Product Features

- Adjustable bar
- Less blisters and red spots on the back of the heels, because the brace moves with the child
- Children can sleep with one leg straight and the other bent which reduces stress on the child’s knees
- Warranty: 12 months



Ordering information

SIZE	LENGTH	PRODUCT NUMBER
Plates 3.9 in. (100 mm)	8 - 12 in. (203 - 304 mm)	697000010

Dobbs Bar Spring Assisted Quick Release

US

Product Features

- For Mitchell Sandal
- Degree dial with a simple turn of screw
- Forefoot extension enable standing in brace
- Adjustable width
- Warranty: 12 months



Ordering information

SIZE	LENGTH	PRODUCT NUMBER
Small	6 - 9 in. (152 - 229 mm)	697000008
Regular	9 - 11 in. (229 - 279 mm)	697000009



Abduction Orthosis

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

ES Manufactured in Spain

COD SU*PER 3D Hip Joint

US

Product Features

- Limits ab-/adduction and internal/external rotation
- To be used with thigh cuffs or KAFO's
- Allows full hip flexion and extension
- Warranty: 12 months

Ordering information

DESCRIPTION	PRODUCT NUMBER
Quick Disconnect	S3DQD



A.A.O. - Anti-Adduction Orthosis

US

Product Features

- Aluminum joint with stainless steel adjustment sliders
- Cuffs are 0.16 in. (4 mm) copolymer
- Silicone medial condylar wedges for comfort included
- Prevents scissoring during gait
- Warranty: 12 months

Ordering information

DESCRIPTION	PRODUCT NUMBER
Small Copolymer Cuffs, 6.5 - 10 in. (165 - 254 mm)	052110
Medium Copolymer Cuffs, 10 in. and above (254 mm and above)	052112
Custom Fabricated A.A.O.	511514
A.A.O. joint only	052100

SPARE PARTS - CONDYLAR WEDGES WITH VELCRO

- A** 127553 - Small
B 127556 - Medium
C 127558 - Large





Scottish Rite Kit

US

Product Features

- Hip abduction spreader bar kit (thigh cuffs not included)
- Warranty: 12 months



Ordering information

FOR THIGH CUFF ATTACHMENT


DESCRIPTION	PRODUCT NUMBER
Small, 7 - 11 in. (178 - 279 mm) Adjustable with Medial Thigh Bar, small	050322
Large, 9 - 15 in. (229 - 381 mm) Adjustable with Medial Thigh Bar, large	050326

FOR BILATERAL KAFO

DESCRIPTION	PRODUCT NUMBER
Small, 7 - 11 in. (178 - 279 mm) Adjustable, small	050320
Large, 9 - 15 in. (229 - 381 mm) Adjustable, large	050324

SUBBORTHOLEN THIGH CUFFS (PAIRS)

DESCRIPTION	PRODUCT NUMBER
Circumference up to 12-1/2 in. (32 cm) (including 2 velcro)	050021
Circumference up to 12-5/8 in. (34.5 cm) (including 2 velcro)	050032
Circumference up to 13-5/8 in. (37 cm) (including 3 velcro)	050039
Circumference up to 14-5/8 in. (39.6 cm) (including 3 velcro)	050047



Thrust Bearing Hip Joint Assemblies



US

Specifications

SIZE	A	B	C	D	E	F	G
Small	13/16	6	31/32	5	5/8	3/16	1-7/8
Medium	15/16	7-1/2	1-1/8	5-1/2	13/16	1/4	1-7/8
Large	1-1/8	8	1-7/16	6	7/8	5/16	2-1/4

Product Features

- Joint size selected on basis of patient weight and height
- Abduction angles at 25° and 35° correspond to bilateral hip abduction angle of 70° and 90° respectively
- Thrust bearings for low friction
- Preset abduction angle
- Various motion control options available




Ordering information

SIZE	HEAD DIAMETER	PROD. NO. 0° ABDUCTION*	PROD. NO. 10° ABDUCTION*	PROD. NO. 25° ABDUCTION*	PROD. NO. 35° ABDUCTION*
Small	31/32 in. (24.6 mm)	058642	058644	058646	058648
Medium	1-1/8 in. (28.6 mm)	058545	058552	058560	058578
Large	1-7/16 in. (36.5 mm)	058586	058594	058602	058610

* Lateral surface of thigh is not parallel to long axis of femur. Abduction angle of femur equals joint angle plus 10°

NOTE: During final assembly of joint, apply #242 Loctite to threads of joint screw and set screw


SIZE	PROD. NO. VARIABLE ABDUCTION	PROD. NO. VARIABLE ABDUCTION WITH XLONG 12 IN. LOWER HIP JOINT BAR
Small	050728	—
Medium	050732	050758
Large	050757	050759





Thrust Bearing Hip Joint Components Scottish Rite

US



Ordering information

	DESCRIPTION	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
A	Joint Screw	027946	026534	027771
B	Thrust Bearing (2 each)	027979	027045	027797
C	Thrust Race (2 each)	027987	027052	027805
D	Spacer (Plastic)	883078	883076	883074
E	Spacer (Inner)	026708	028621	028647
F	Upper Hip Joint Bar	027904	026955	027730
G	Joint Nut	027953	026542	028399
H	Set Screw	027961	027961	027896
I	Lower Hip Joint Bar, 0°	050890	050716	050708
I	Lower Hip Joint Bar, 10°	050893	050724	050542
I	Lower Hip Joint Bar, 25°	050895	050617	050559
I	Lower Hip Joint Bar, 35°	050897	050625	050567
J	Motion Control Plate	808528	808529	808537
K	Stop Screw	880670	880670	880672
L	Nut	881003	881003	881011
M	Variable Abduction Lower Hip Joint Bar	050772	050774	050776
M	Variable Abduction Lower Hip Joint Bar, XLong	-	050778	050780
N	Drop Lock	050869	050870	050874
N	Abd, Screw, Set, 1/4"-28 x 7/16 in., (Not Shown)	-	880303	880305
N	Abd. Screw, Adj., 1/4"-28 x 3/4 in., (Not Shown)	-	880306	880306
N	Abd, Screw, Set, 10-32 x 3/8 in., (Not Shown)	880025	880311	880312

Motion Control Plate and Drop Lock Scottish Rite

US

Ordering information

CIRCUMFERENCE RANGE	PROD. NO. SMALL	PROD. NO. MEDIUM	PROD. NO. LARGE
Motion Control Plate	058528	058529	058537
Drop Lock Only	050869	050870	050874
Drop Lock and Control Plate	050876	050878	050882



Hip Joint Assemblies, Stainless Steel Scottish Rite

US

Ordering information


DESCRIPTION	MEASUREMENTS	PRODUCT NUMBER
Small	1/8 x 5/8 x 5-3/8 in. (3.2 x 15.9 x 136.5 mm)	050336
Large	3/16 x 3/4 x 8-1/2 in. (4.8 x 19.1 x 215.9 mm)	050070
Screw, Fillister Head (Sold Separately)	10-32 x 5/16 in.	880906
Screw, Pan Head (Sold Separately)	1/4"-28 x 7/16 in.	029181

SPARE PARTS - SMALL

- A** 050468 - Upper
- B** 029223 - Joint Nut, 3/8 x 23/64 in.
- C** 050351 - Lower, 1/8 x 5/8 x 5-3/8 in.

SPARE PARTS - LARGE

- A** 050484 - Upper, 3/8 in. Offset
- B** 029223 - Joint Nut, 3/8 x 23/64 in.
- C** 050096 - Lower, 3/16 x 3/4 x 8-1/2 in.





Component Kits with Telescoping Spreader Bar - Scottish Rite

US

Product Features


- Refer to Central Fabrication Price List for custom orthosis
- Measurement sheet available
- Reprints available of Dr. Wood W. Lovell's paper, The Scottish Rite Hospital Orthosis for Legg-Perthes Disease

Ordering information

SIZE	PATIENT HEIGHT	PRODUCT NUMBER
Small, 1-2 years	< 36 in. (< 914 mm)	050435
Small Partial Kit, 1-2 years*	< 36 in. (< 914 mm)	050393
Large	> 36 in. (> 914 mm)	050427

SPARE PARTS - SMALL KIT & PARTIAL SMALL KIT (NOTE: PARTIAL KIT ONLY INCLUDES PARTS A & B)	SPARE PARTS - LARGE KIT
A 050336 - Hip Joint Assembly, 2 Each	A 050070 - Hip Joint Assembly, 2 Each
B 050344 - Telescoping Spreader Bar Assembly	B 050310 - Telescoping Spreader Bar Assembly
C 027708 - Pelvic Band, 6 x 1-1/2 x 20 in.	C 027710 - Pelvic Band, 3/16 x 1-3/4 x 20 in.
D 050385 - Medial Thigh Bar, 2 Each	D 050161 - Medial Thigh Bar, 2 Each
E 050375 - Lateral Thigh Bar, 2 Each	E 050139 - Lateral Thigh Bar, 2 Each
F 885906 - Button Head Socket Cap Screw, 10 Each	F 885906 - Screws, 10-32 x 1/4 in., Button Head Socket Cap Screw, 10 Each
	885904- Screws, 10-32 x 5/16 in., Button Head Socket Cap Screw, 6 Each

* Only includes part A and B




Fixed 25° and 35° Abduction Thrust Bearing Hip Joints - Scottish Rite

US

Product Features

- Contact Central Fabrication for custom orthosis
- Measurement sheet available
- Reprints available of Dr. Wood W. Lovell's paper, The Scottish Rite Hospital Orthosis for Legg-Perthes Disease



Ordering information

25° ABDUCTION COMPONENT KIT

SIZE	PRODUCT NUMBER
Small	050454
Medium	050810
Large	050826

SPARE PARTS - SMALL

A 027708 - Pelvic Band	C 050375 - Lateral Thigh Bars, 2 Each
B 058646 - Thrust Bearing Hip Joint 25°, 2 Each	D 885904 - Screws, 10-32 x 5/16 in., Button Head Socket Cap Screw, 6 Each

SPARE PARTS - MEDIUM

A 027708 - Pelvic Band	C 050139 - Lateral Thigh Bars, 2 Each
B 058560 - Thrust Bearing Hip Joint 25°, 2 Each	D 885902 - Screws, 10-32 x 3/8 in., Button Head Socket Cap Screw, 6 Each

SPARE PARTS - LARGE

A 027710 - Pelvic Band, 3/16 x 1-3/4 x 20 in.	C 050139 - Lateral Thigh Bars, 2 Each
B 058602 - Thrust Bearing Hip Joint 25°, 2 Each	D 885900 - Screws, 10-32 x 1/2 in., Button Head Socket Cap Screw, 6 Each

35° ABDUCTION COMPONENT KIT

SIZE	PRODUCT NUMBER
Small	050456
Medium	050818
Large	050834

SPARE PARTS - SMALL

A 027708 - Pelvic Band, 3/16 x 1-1/2 x 20 in.	C 050375 - Lateral Thigh Bars, 2 Each
B 058648 - Thrust Bearing Hip Joint 35°, 2 Each	D 885904 - Screws, 10-32 x 5/16 in., Button Head Socket Cap Screw, 6 Each

SPARE PARTS - MEDIUM

A 027708 - Pelvic Band, 3/16 x 1-1/2 x 20 in.	C 050139 - Lateral Thigh Bars, 2 Each
B 058578 - Thrust Bearing Hip Joint 35°, 2 Each	D 885902 - Screws, 10-32 x 3/8 in., Button Head Socket Cap Screw, 6 Each

SPARE PARTS - LARGE


A 027710 - Pelvic Band, 3/16 x 1-3/4 x 20 in.	C 050139 - Lateral Thigh Bars, 2 Each
B 058610 - Thrust Bearing Hip Joint 35°, 2 Each	D 885900 - Screws, 10-32 x 1/2 in., Button Head Socket Cap Screw, 6 Each

Variable Abduction Thrust Bearing Hip Joints - Scottish Rite

US

Product Features

- Contact Central Fabrication for custom orthosis
- Measurement sheet available
- Reprints available of Dr. Wood W. Lovell's paper, The Scottish Rite Hospital Orthosis for Legg-Perthes Disease




Ordering information

SIZE	PRODUCT NUMBER
Small	050458
Medium	050842
Large	050850
SPARE PARTS - SMALL	
A 027708 - Pelvic Band, 3/16 x 1-1/2 x 20 in.	C 050375 - Lateral Thigh Bars, 2 Each
B 050728 - Thrust Bearing Hip Joint, 2 Each	D 885902 - Screws, 10-32 x 3/8 in., Button Head Socket Cap Screw, 6 Eac
SPARE PARTS - MEDIUM	
A 027708 - Pelvic Band, 3/16 x 1-1/2 x 20 in.	C 050139 - Lateral Thigh Bars, 2 Each
B 058560 - Thrust Bearing Hip Joint, 2 Each	D 885902 - Screws, 10-32 x 3/8 in., Button Head Socket Cap Screw, 6 Eac
SPARE PARTS - LARGE	
A 027710 - Pelvic Band, 3/16 x 1-3/4 x 20 in.	C 050139 - Lateral Thigh Bars, 2 Each
B 050757 - Thrust Bearing Hip Joint, 2 Each	D 885900 - Screws, 10-32 x 1/2 in., Button Head Socket Cap Screw, 6 Each

Pelvic Bands and Thigh Bars Scottish Rite

US



Ordering information

DESCRIPTION	MEASUREMENTS	PRODUCT NUMBER
A Pelvic Band, Aluminum	3/16 x 1-1/2 x 20-1/2 in. (4.8 x 38.1 x 520.7 mm)	027708
A Pelvic Band, Aluminum	3/16 x 1-3/4 x 24 in. (4.8 x 44.5 x 609.6 mm)	027710
B Lateral Thigh Bar with Loop	1/8 x 5/8 x 3-3/8 in. (3.2 x 15.9 x 85.7 mm)	050375
B Lateral Thigh Bar with Loop	1/8 x 3/4 x 3-1/2 in. (3.2 x 19.1 x 88.9 mm)	050137
B Lateral Thigh Bar with Loop	1/8 x 3/4 x 5-1/2 in. (3.2 x 19.1 x 139.7 mm)	050139
C Medial Thigh Bar	1/8 x 5/8 x 2-1/2 in. (3.2 x 15.9 x 63.5 mm)	050385
C Medial Thigh Bar	1/8 x 3/4 x 4-1/2 in. (3.2 x 19.1 x 114.3 mm)	050161

Telescoping Spreader Bars Scottish Rite

US



Ordering information

DESCRIPTION	ROD ENDS	PRODUCT NUMBER
12 in. (304.8 mm) Small	- 5/16 in. with 3/8 in., Stainless Steel Tubing - 1/2 in. Stainless Steel Tubing with Two Medial Attachment	050344
16-1/2 in. (419.1 mm) Long	- 3/8 in. Rod End with 7/16 in. Stainless Steel Tubing - 3/8 in. Rod End with 9/16 in. Stainless Steel Tubing with Two Medial Attachment Bars	050310



Index

2400 Bar Kit	47	Bearing Hip Joints - Scottish Rite	105
2400 Joints	46	Flex Foam II	80
A.A.O. - Anti-Adduction Orthosis	99	Flex Foam III	80
Aluminum Pelvic Bands	71	Flex Foam™ I	79
Ankle Joint - Double Action	36	Gauges	74
Ankle Joint - Limited/ Free Motion	35	Hip Joint Assemblies, Stainless Steel	
Ankle Joint - Toe Pick Up	35	Scottish Rite	103
Anterior Stay	89	Hoop Cable RGO Cables and Housings	59
A.O.O. Spinal Brace	88	Horizontal RGO Cables	67
Bar Kit	47	Insert Stirrup For Plastic Orthotics	40
Basic Strap	89	Integrated Ankle Joint	
Cable Track Kit	71	Adjustable range of motion	30
Catalog symbol legend	9	Integrated Ankle Joint	
Chafes with Triangular Body Clasp	90	Dorsiflexion Assist	31
Clip-On Bar	96	Integrated Ankle Joint, Double Action	30
Close Contour	50	Interlock™	45
Close Contour Pediatric Knee Joint	48	Knee Joints	73
COD SU*PER 3D Hip Joint	99	Large Hooped Cable RGO with Abduction	58
Component Kits with Telescoping Spreader Bar - Scottish Rite	104	Large Hooped Cable RGO without Abduction	57
Denis Browne Splint	96	Large Horizontal Cable RGO with Abduction	66
Dobbs Bar	97	Large Horizontal Cable RGO without Abduction	65
Dobbs Bar Spring Assisted Quick Release	97	Latch Knob for Small Hooped Cable RGO	60
Dorsiflexion Assist Ankle Joint Kits With Polycar-C	42	Latch Knob for Small Horizontal Cable and Rocker Bar RGOs	68
Drop Lock Aluminum Hip Joint		Lever Release Kit - LRK	49
Assemblies	52	Lower Abduction Bar Assembly	72
Drop Lock Aluminum Knee Joint		Lower Bar Assembly without	
Assemblies	52	Abduction	72
Drop Lock Low Profile Aluminum Knee Joint Assemblies	53	LSO - Bi-valve	78
Dynamic Walk - Custom Made	25	LSO - Single Opening	77
Dynamic Walk - Double-Double	15	Medium Horizontal Cable RGO with Abduction	64
Dynamic Walk - Mix & Match	22	Medium Horizontal Cable RGO without Abduction	63
Dynamic Walk - Mix & Match		Modular Integrated Ankle Joint,	
Custom Calf Height or Foot Length	23	Adjustable Range of Motion	29
Dynamic Walk - Modular	21	Modular Integrated Ankle Joint,	
Dynamic Walk - Shoe Mount	20	Double Action	29
Dynamic Walk - Single Side Anterior	19	Motion Control Plate and Drop Lock	
Dynamic Walk - Single Side Lateral	18	Scottish Rite	103
Dynamic Walk - Single Side Medial	17	Order form - Dynamic Walk Custom Calf Height or Foot Length	24
Dynamic Walk - Size Charts	11	Order form - Dynamic Walk Custom Made Adult and Pediatric	26
Dynamic Walk - Standard	13	Order Form - Flex Foam	85
Dynamic Walk Standard, Double Double	11	Order Form - LSO/ TLSO, Single or Bivalve Opening	84
Dynamic Walk - Wide Ankle	14	Order Form - Pectus Brace	86
Dynamic Walk - Wide Ankle Double-Double	16	Order Form - Scoliosis/ Kyphosis Brace	83
External Ankle Foot Orthosis (AFOs)	74	PDA Ankle Joint Kits With Polycar-C	43
Fillauer Product Warranty	6	Pectus Brace	81
Fixed 25° and 35° Abduction Thrust		Pelvic Bands	73
		Pelvic Bands and Thigh Bars	



Scottish Rite	107	Zippy Quick Release	51
PEVA Shoe - Open Toe Boots - Firm Heel Cap	94		
Polycar-C™ Ankle Inserts	41		
Providence Nighttime Brace	76		
Push Button for Hooped Cable RGO	60		
Push Button for Horizontal Cable and Rocker Bar RGOs	69		
Rectangle Caliper Plates	40		
Release Kit RBK	48		
Release Kit SBK	48		
Release Kit - UBK	49		
Returns, Adjustments, and Credits	7		
RGO II System for Hooped Cable RGO with Abduction	61		
RGO II System for Horizontal Cable and Rocker Bar RGOs	69		
RGO/IRGO Overview	55		
Scoliosis Full Time Brace	77		
Scoliosis Module	87		
Scottish Rite Kit	100		
Short Leg Upright - Double Action	33		
Short Leg Upright - limited/ free motion	33		
Short Leg Upright - Toe Pick Up	34		
Small Hooped Cable RGO	56		
Small Horizontal Cable RGO	62		
Solid Stirrups - Wide Tongue	32		
Specialty Connectors for ISOCENTRIC® RGO	70		
Split Tail Stirrups	32		
Spring Loaded Lever Lift	53		
Spring Steel Strap	90		
Standard In-Shoe Ankle Foot Orthosis (AFO)	74		
Steplock®	44		
Stirrups			
Long Tongue & Wide Tongue	39		
Stirrups - Split Tongue	38		
Stirrups - Standard	37		
Telescoping Spreader Bars			
Scottish Rite	107		
TFC - Thermoformable Composite Insole	92		
Thrust Bearing Hip Joint Assemblies	101		
Thrust Bearing Hip Joint Components			
Scottish Rite	102		
TLSO - Bi-valve	79		
TLSO - Single Opening (A.O.O.)	78		
Trigger Release for Camlocks	51		
Variable Abduction Thrust Bearing			
Hip Joints - Scottish Rite	106		
Welding Rod	89		

