

Product Order Form



FACILITY INFORMATION	
Clinician name:	
Company name:	
Shipping address:	
Phone:	
Email:	
Purchase order #:	

PATIENT INFORMATION:	
Trial reference (if applicable):	
Insurance provider (if applicable):	
History:	<input type="checkbox"/> New prosthetic user <input type="checkbox"/> Prosthetic renewal
Previous device (if applicable):	

PLEASE SPECIFY THE MODEL OF OTHER COMPONENTRY USED IN THE PROSTHESIS (if applicable):

Control system:		Wrist connector/rotator:	
Elbow:		Other:	



DIMENSION			
	7¼" hand	8¼" hand	
A	Middle finger tip to hand base (fully open)	179mm	180mm
B1	Hand base to wrist base (QDW)	17mm	17mm
B2	Hand base to wrist base (LPW inc. Lamination Collar)	12mm	12mm
C	Diameter at wrist	54mm	54mm
D	Thumb tip to hand base	123mm	123mm
E	Palm width	82mm	88mm
F	Palm circumference	7¼ in (197mm)	8¼ in (210mm)

HAND MODEL

(see table to right for full list of items included with the TASKA Hand)

Size:	<input type="checkbox"/> 7¼ (medium) <input type="checkbox"/> 8¼ (large)
Side:	<input type="checkbox"/> Left <input type="checkbox"/> Right
Wrist type:	<input type="checkbox"/> Quick Disconnect <input type="checkbox"/> Low Profile
Colour:	<input type="checkbox"/> Black <input type="checkbox"/> White <input type="checkbox"/> Sand

INCLUDED IN HAND PACKAGE*

(Quick Disconnect Wrist QDW & Low Profile Wrist LPW)

	QDW	LPW
2000 mAh Battery, Chargers, Power Switch	•	•
2 year standard warranty	•	•
Annual service with each year of warranty (inc. ext.)	•	•
2 x Bluetooth adaptors	•	•
Replacement Grip Tips and Faceplates	•	•
Low Profile Lamination Ring, Switch Block	n/a	•

WARRANTY EXTENSION (sold separately)

1 year 2 year 3 year

QUICK DISCONNECT WRIST COMPONENTS (sold separately to TASKA Hand)

<input type="checkbox"/> QDW Lamination Collar Kit (4-band)	Included in kit: (hand-side) 2 x TASKA Seal Ring (socket-side) TASKA Lamination Collar, Coaxial plug (4-band)
<input type="checkbox"/> QDW Lamination Collar Kit (6-band)	Included in kit: (hand-side) 2 x TASKA Seal Ring (socket-side) TASKA Lamination Collar, Coaxial plug (6-band)
<input type="checkbox"/> MC Seal Ring Kit	For waterproof connection to compatible Motion Control electric wrist rotators