

## Order form - Dynamic Walk Custom Calf Height or Foot Length

Customer: ..... Date: .....  
 Patient reference/ order number: ..... Ordered by: .....  
 Phone: ..... Email: .....

**Patient information:**

Age: \_\_\_\_\_ Weight (kg): \_\_\_\_\_

**Activity level:**

Low      Medium      High

**Side:**

Left      Right      Bilateral

Type (choose **one**): 700 100 105

- Standard
- Standard Wide Ankle
- Double - Double
- Double - Double Wide Ankle
- Single Side Lateral
- Single Side Medial
- Shoe Mount - Custom Height

**Foot Plate**

Small	Increase/Decrease Length (mm)
Medium	+ _____
Large	- _____


**Calf Band**

Height is determined by model and size of calf section. (See chart below)

Small	Increase/Decrease Height (mm)
Medium	+ _____
Large	- _____

**Material Options:**

<b>Foot Plate:</b>	<b>Calf band:</b>
Standard	Standard
Rigid foot plate (for diabetics)	Rigid calf band (more stability)
Flexible foot plate (high activity users)	



Measurements (mm):

A \_\_\_\_\_ B \_\_\_\_\_

Explanation:  
 A - Height of orthosis  
 B - Foot length

FOOT LENGTH	SMALL	MEDIUM	LARGE
All types except Shoe Mount	230 mm +/- 13 mm	250 mm +/- 13 mm	270 mm +/- 13 mm
CALF CUFF SELECTION	SMALL	MEDIUM	LARGE
Calf Circumference	240 - 330 mm	340 - 420 mm	430 - 520 mm
CALF HEIGHT	SMALL	MEDIUM	LARGE
Standard / Standard Wide Ankle	330 mm +/- 25 mm	355 mm +/- 30 mm	380 mm +/- 35 mm
Double Double / Double Double Wide Ankle	345 mm +/- 25 mm	380 mm +/- 30 mm	400 mm +/- 35 mm
Single Side Lateral / Medial	345 mm +/- 25 mm	370 mm +/- 30 mm	395 mm +/- 35 mm
Shoe Mount	360 mm +/- 25 mm	385 mm +/- 30 mm	410 mm +/- 35 mm

**NOTE: Padding is ordered separately.**