

DuraTwist

Features & Benefits

- Reduction of shear forces on the residual limb during rotary motions
- Large range of transverse plane motion ($\pm 20^\circ$ rotation at 20 N·m torque, $\pm 30^\circ$ max)
- 1 – 2 mm vertical compression
- Lightweight
- Small size
- Waterproof (1 meter for 30 minutes)

Specifications

Product Number	DT002	DT001
Color	Blue	Black
Patient Weight	< 180 lbs. (82 kg)	< 250 lbs. (113 kg)
Activity Level	K2 – K3	K2 – K3
Build Height	2.625 in. (6.6 cm)	2.625 in. (6.6 cm)
Product Weight	8.8 oz. (250 g)	8.8 oz. (250 g)
Torque Setting	15 N·m	15 N·m
Range of Motion	20° Internal / External Rotation	
Materials	Titanium and Urethane Rubber	
Durable	Meets ISO-10328 Standard	
Waterproof	This device is rated for submersion to 1 meter for 30 minutes	
Warranty	12 months	



DT002



DT001

Fillauer®

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946