

Reimagining Possibilities

The **Dynamic Walk AFO** is designed to allow people with dorsiflexor paresis or paralysis to confidently regain their stride. Whether their goal is to walk, hike, jog, or play, the Dynamic Walk's flexible dorsi-assist joints allow them to move freely. Featuring a heat-moldable, trimmable calf band and foot plate, the Dynamic Walk ensures a comfortable fit every time.

For more information, please visit **www.fillauer.se** or contact us at

fillauer@fillauer.com



Use the QR Code to learn more about the Dynamic Walk Series





Dynamic Walk Series



Standard



Standard Wide Ankle



Double-Double



Wide Ankle Double-Double



Single Side Lateral



Single Side Medial



Single Side Anterior



Shoe Mount



Modular



Mix & Match



Custom Calf Height or Foot Length



Custom

How to choose the right Dynamic Walk

To achieve a successful Dynamic Walk fitting, a thorough patient analysis is key and should include the following: other disabilities, muscle strength, stability in ankle and knee joint, walking pattern and usage. Once this information has been gathered and analysed you can choose the most suitable Dynamic Walk.

	ML STABILITY	WIDE ANKLE	STRONGER LIFT	KNEE HYREREX- TENSION	WEAK QUADRICEPS	LESS HAND/ARM FUNCTION	OEDEMA PROBLEMS
Standard	Х			Х			
Wide Ankle	Х	Х		Х			
Double-Double	Х		Х	Х			
Wide Ankle Double- Double	Х	Х	Х	Х			
Single Side Medial				Х		Х	Х
Single Side Laterial				Х		Х	Х
Shoe Mount				Х		х	Х
Single Side Anterior					Х	Х	Х

